


Please note: The receiver block mounting hardware is seating system specific. Visually inspect the side rails for proper selection or contact Technical Service for more information.



Item	Part No.	Name	Qty
A1	FRMASMB8586	ARMREST,ASSY,RECEIVER BLOCK,QUANTUM SINGLE POST HEIGHT ADJUSTABLE,D-1354-019	1
B1	HDW142473	SCREW,SHCS,BLACK ZINC,M6-1.0 X 18MM,ISO 4762 CLASS 12.9,W/NYLON PATCH	2
C1	MEC131114	SQUARE NUT,T-SLOT, M6X1.0	2
D1	SCRSHCS1005	SCREW, SHCS,BLACK ZINC,1/4-20 X 0.75"	2
E1	DWR1360A002	T-NUT, 1/4-20 X 0.5625",A-1360-002	2