

Model Number	JAZZYAIR2
Weight Capacity	300 lbs.
Turning Radius	24.25"
Overall Width ³	23.25"
Overall Length ³	33.21"/ 42" (with footplate)
Maximum Speed ²	Up to 4 mph
Ground Clearance ³	1" at motor bracket
Front Wheels	6" solid casters
Drive Wheels	10" solid
Rear Wheels	6" solid casters
Drivetrain	Two motor, Mid-Wheel 6® drive
Braking System	Regenerative and electromechanical
Suspension Type	Active-Trac® (ATX) with Mid-Wheel 6® Technology
High-back Seat Sizes	16"x16"-18", 18"x18"-20", 20"x 20"-22"
Seat-to-Floor Range	21.25" – 33"
Per Charge Range ^{3,5}	Up to 16.6 miles at 220 lbs. Up to 15 miles at 300 lbs.
Standard Electronics	40A, PG VR2
Component Weights ⁷	Base: 252 lbs. w/ U1 batteries 267 lbs. w/ 40 Ah batteries U1 Batteries: 24.5 lbs. 40 Ah Batteries: 32.5 lbs.
Batteries ^{4,6}	12 volt, U-1 (standard), 40 Ah (optional)
Battery Charger	5 amp, off-board
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

HCPCS Code: K0899
POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- ¹ Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ AGM or gel cell type required.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁶ Battery weight may vary based on manufacturer.
- ⁷ Weight includes base, seat and batteries.
- ⁸ WC-19 Occupied Transit.



**PRIDE
MOBILITY®**

**Jazzy®
air₂**

Elevate your life!



FDA Class II
Medical Device¹

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

helpis@alliedmedical.co.nz | www.alliedmedical.co.nz | 0800 31 61 81

Elevate your life!

Height and angle adjustable headrest

High-back or solid seat pan, reclining sport seat featuring 2" of depth adjustment.

Seat sizes:

16"x16"-18"
18"x18"-20"
20"x20"-22"

Angle, height, and depth-adjustable armrests

Convenient "Air" switch on armrest activates elevation

Under seat storage

Lift mechanism elevates up to 12" while driving up to 4 mph

Integrated foot platform

Rear LED lights

Bright LED marker lights

FDA Class II Medical Device¹

Active-Trac® Suspension with Mid-Wheel 6® Technology for a smooth and stable ride.

Bright LED lighting for increased visibility.

Easily accessible "Air" switch activates elevation.

Convenient under seat storage.



AT THE TOUCH OF A SWITCH, the Jazzy Air® 2 elevates 12" in just 11 seconds. With a sleek, more modern look, this power chair safely drives at speeds up to 4 mph while elevated, allowing users to enjoy face-to-face social engagement. With Mid-Wheel 6® Drive and Active-Trac® Suspension, the Jazzy Air 2 offers enhanced performance and stability both indoors and out.

KEY FEATURES

- 12" of power adjustable seat height while driving
- Angle, width, height and depth adjustable armrests
- Integrated foot platform with rubber insert and 2" of adjustment for comfort in both the elevated and non-elevated positions
- Single-sided caster forks for enhanced maneuverability
- WC-19 approved occupied transit

AVAILABLE COLORS



ACCESSORIES

- Cane holder
- Cell phone holder
- Cup holder
- Flag kit
- 50", 60", 70" Lap belt
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- XLR USB Charger

OPTIONS

- 40-amp AGM batteries
- Swing-away joystick
- Stealth® cushions