Jazzy®Air

AU

AGN# 088 609 661

Illustrated Parts Breakdown



Pride®



Your new Pride Illustrated Parts Breakdown (IPB) Manual has been compiled from the latest product information available at the time of publication. Pride® Mobility Products Corporation reserves the right to make changes and updates to this manual and to the product(s) it represents. Any changes, updates, or improvements to the product(s) may result in slight differences between this manual and the product you purchased.

An IPB is an illustrated parts listing of all product components which have maintenance significance. IPB's utilize "exploded" drawings to detail the individual components for each assembly and an itemized list of reference numbers, part numbers, descriptions, and quantities per assembly. Some parts may be included in an assembly but may not be available for purchase separately. Refer to the parts list to locate the part numbers for assemblies and components for ordering what you need.

To achieve proper identification and relationship of detail parts and next higher assemblies, parts listings are presented in a convenient alphanumeric reference system. Reference numbers for parent assemblies identify one part number that can be used to order all of the components listed within it. Components of an assembly will have the same alpha character as its parent assembly.

■ For example, the parent assembly part number using reference number of A1 — A4 is the only number you need to place an order for parts with references of A1, A2, A3, and A4. Instead of ordering 4 separate parts, use the parent assembly part number to get the complete assembly.

Master assembly part numbers have reference numbers that include several subassemblies. By ordering parts using master assembly part numbers, you can save a lot of time and reduce the chance for error.

For example, A1 — G10, a master assembly reference number, includes all components of assemblies A, B, C, D, E, F, and G. Use caution when ordering assemblies. Check to see if a master assembly part number exists for the components that you need. If an entire master assembly is not needed, then the individual assembly part number (i.e. for reference number G1 — G10) will be sufficient to use for ordering.

Asterisks (*) and NOTES are used throughout this service manual to point out important information regarding assembly structure and contents.

■ For example, in some illustrations, only the left OR right side is shown. In other drawings, only the front OR back is depicted. The part numbers for ALL components are listed. An asterisk will also be used to indicate which components are not illustrated. Illustrated Items without callouts are shown for assembly purposes.

Updates to the contents of this IPB will be made on a regular basis. Providers may retrieve updates immediately upon their release by going to the Pride Service Website. The web address is www.pridemobility.com. We welcome any questions or comments regarding the new Pride IPB Manual that you may have.

Please note all pricing included with this manual is subject to change without prior notice from Pride.

Thank you for choosing products from Pride[®] Mobility Products Corporation.


Copyright © 2018 Pride® Mobility Products Corporation www.pridemobility.com INFMANU4765/Rev.E/031318



SECTION I.

Frame Assembly	
Main Frame & Actuator	4
Front Caster Arm with Motor	6
Rear Articulating Beam	8
SECTION II.	
Electronics Assembly	
SECTION III.	
Joystick, Buddy Button & Light Switch	12
SECTION IV.	
Charger	13
SECTION V.	
Foot Platform and Seat Interface	
Mounting Plate for Foot Riggings	
SECTION VI.	
Seat Assembly	18
SECTION VII.	
Shroud Assembly	
SECTION VIII.	
Accessories	
Cane/Crutch Holder	
Oxygen Holder	
Rear Basket	
Walker Holder	_
Docking Device Assembly	
SECTION IX.	
Drive Wheel Assembly	27
SECTION X.	
Owner's Package	
SECTION XI.	
Decal Package	20





I. Frame Assembly — Main Frame & Actuator




REF#	PART#	DESCRIPTION	QTY	11i MSRP (AU\$)
A1	MOT153402	THOMSON ACTUATOR WITH LIMIT ASSEMBLY, LC2000	1	6135.82
B1	MEC155338	IGUS CHAIN ASSEMBLY,W/SUPPORT LINKS,P668765	1	262.9
C1	HDW156648	SCREW,SFHCS,M3 X 0.5 X 8MM,BLACK ZINC,CLASS 10.9	2	6.04
D1	MEC153194	FRAME WELDMENT, JAZZY SELECT AIR	1	N/A
E1	MEC153356	ARTICULATING BUMPER STOP	2	N/A
*F1	MEC154624	BATTERY STRAP	2	37.04
*G1	MEC154626	BACK STRAP,SMALL BATTERY	2	24.66

Notes:

^{*} Not shown.






I. Frame Assembly — Front Caster Arm with Motor



REF#	PART#	DESCRIPTION	QTY	11 MSRP (AU\$)
A1	MEC155351	FRAME, ASSY, FRONT ARM W/CASTER & MOTORS, RIGHT, JAZZY AIR	1	1491.68
B1	MEC155350	FRAME, ASSY, FRONT ARM W/CASTER & MOTORS, LEFT, JAZZY AIR	1	1491.68
C1	HDWASMB1078	ANTI-FLUTTER KIT,ASSY,JAZZYSELECT,B-1060-050	1	68.4
D1	CASASMB7110004	CASTER,ASSY,6" FLAT-FREE WHEEL,SILVER RIM/BLACK TIRE.RIGHT.JAZZYSELECT.VERSION 2.D-1202-034	1	177.3
E1	CASASMB7110002	CASTER,ASSY,6" FLAT-FRÉE WHEEL,SILVER RIM/BLACK TIRE,LEFT,JAZZYSELECT,VERSION 2,D-1202-034	1	177.3
F1	WHLASMB2047	WHEEL, ASSY, CASTER, FLAT FREE, SILVER RIM/BLACK TIRE, 6" X 2", W/BEARINGS & SPACER, E-1202-019	1	180.96
G1	MEC154375	FREEWHEEL LEVER, RIGHT	1	10.88
H1	MEC154374	FREEWHEEL LEVER, LEFT	1	10.88





When applicable, assemblies are grouped by color. All components with that color are included in the assembly.

A1	Complete Articulating Beam Assy (Includes B1 - E1)
B1	Anti Flutter Kit
C1	Right Caster Assy (Includes E1)
D1	Left Caster Assy (Includes E1)
E1	Wheel Assy




I. Frame Assembly — Rear Articulating Beam



REF#	PART#	DESCRIPTION	QTY	11i MSRP (AU\$)
A1	MEC153073	CASTER BEAM ASSEMBLY,6-INCH	1	484.52
B1	HDWASMB1078	ANTI-FLUTTER KIT,ASSY,JAZZYSELECT,B-1060-050	2	68.4
C1	CASASMB7110004	CASTER,ASSY,6" FLAT-FREE WHEEL,SILVER RIM/BLACK TIRE,RIGHT,JAZZYSELECT,VERSION 2,D-1202-034	1	177.3
D1	CASASMB7110002	CASTER,ASSY,6" FLAT-FREE WHEEL,SILVER RIM/BLACK TIRE,LEFT,JAZZYSELECT,VERSION 2,D-1202-034	1	177.3
E1	WHLASMB2047	WHEEL,ASSY,CASTER,FLAT FREE,SILVER RIM/BLACK TIRE,6" X 2",W/BEARINGS & SPACER,E-1202-019	1	180.96






II. Electronics Assembly




REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	ELE156112	ELECTRONIC,ASSY,VR2,JAZZY AIR,OFFBOARD CHARGER	1	1965.48
B1	DWR1010H128	ELECTRONIC,HARNESS,VR2 EXTENSION CABLE,H-1010-128,(P&G: SA77974)	1	54.5
C1	ELE154541	HARNESS,SUSPENSION LOCK,JAZZY AIR	1	38.96
D1	ELE154823	HARNESS,LEFT SUSPENSION SWITCH,JAZZY AIR	1	32.1
E1	ELE154824	HARNESS,RIGHT SUSPENSION SWITCH,JAZZY AIR	1	32.1
F1	ELE155078	JAZZY AIR LIGHTING HARNESS	1	78.34
G1	ELE155079	HARNESS, JAZZY AIR, BUDDY BUTTON EXTENSION HARNESS	1	61.96
H1	SWTMCRO1081	ELECTRONIC,SWITCH,BLACK,EGG,(STEALTH:ESBLK)	1	161.7
I1	DWR1060H001	ELECTRONIC,HARNESS,BATTERY CABLE,GC,JAZZYSELECT,H-1060-001	2	116.74
J1	ELECBRK1001	ELECTRONIC,CIRCUIT BREAKER,25AMP,SC150	1	83.56
K1	ELE155648	HARNESS, LIFT ACTUATOR SWITCH, JAZZY AIR	1	38.8
L1	ELE154604	JAZZY AIR POWER HARNESS	1	147.16
M1	CTL156111	CONTROLLER,VR2,ASSY,40 AMP,1 ACTUATOR,JAZZY AIR PROGRAM	1	1208.38
N1	ELE154790	JAZZY AIR LED BOARD	1	142.64






REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	CTL154825	CONTROLLER, JOYSTICK MODULE, 6 KEY, VR2 2.0 METER CABLE (P&G:D51145)	1	1094.14
B1	MEC156107	JOYSTICK MOUNTING BRACKET,ASSY,GC/VR2 CONTROLLER,INLINE,JAZZY AIR	1	26.44
C1	SWTMCRO1081	ELECTRONIC,SWITCH,BLACK,EGG,(STEALTH:ESBLK)	1	161.7
D1	MEC168058	UNDER ARM BUDDY BUTTON MOUNTING BRACKET ASSY, JAZZY AIR	1	50.1
E1	ELE160094	LIGHT SWITCH,ASSEMBLY,JAZZY AIR	1	61
F1	MEC168060	RETRACT 4 SWING AWAY JOYSTICK ASSY, JAZZY AIR	1	1100





REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	ELEASMB6657	CHARGER,ASSY,OFFBOARD,LOW INHIBIT,100V-240V,3.5 AMP,W/XLR HARNESS,CE/CUL/UL APPROVED,W/AUSTRALIAN POWERCORD	1	N/A
B1	DWR1000H174	ELECTRONIC,HARNESS,AC POWER CORD,INSULATED BLADES,AUSTRALIA,H-1000-174	1	44.98








REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	MEC153429	SELECT AIR FOOTREST	1	180.7
B1	MEC160111	SEAT INTERFACE ASSY, JAZZY AIR	1	231.3
C1	HDW125250	SCREW,SOCKET HEAD,M8X1.25 X 20,ISO 4762,CL 12.9,BZ,PVTP	4	6.04



When applicable, assemblies are grouped by color. All components with that color are included in the assembly.

A1	KIT1801691 - Includes Complete Assy S/A #8, Clover Leaf Seat Frame Assy
B1	KIT1801692 - Includes Compete Assy Taper Pin, Clover Leaf Seat Frame Assy
C1	Seat Mounting Screws




V. Mounting Plate for Foot Riggings



REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	KIT1801691	KIT ASSEMBLY,MOUNTING,SWING A WAY LEGRESTS,JAZZY AIR	1	N/A
B1	KIT1801692	KIT ASSEMBLY,MOUNTING,ELEVATING LEGRESTS,JAZZY AIR	1	N/A
C1	HDW125250	SCREW, SOCKET HEAD, M8X1.25 X 20, ISO 4762, CL 12.9, BZ, PVTP	4	6.04






REF#	PART#	DESCRIPTION	QTY	11 MSRP (AU\$)
A1	SET1701577	SEAT,ASSY,FOLD FLAT CAPTAIN SEAT,BLACK VINYL,16 X 16-18 W/PRIDE LOGO,60" LAP BELT,(VISCOUNT: PR-102)	1	1399.12
B1	SET1701578	SEAT,ASSY,FOLD FLAT CAPTÂIN SEAT,BLACK VIŃYL,18 X 18-20 W/PRIDE LOGO,60" LAP BELT,(VISCOUNT: PR-103)	1	1399.12
C1	SET1701579	SEAT,ASSY,FOLD FLAT CAPTAIN SEAT,BLACK VINYL,20 X 20-22 W/PRIDE LOGO,60" LAP BELT,(VISCOUNT: PR-104)	1	1466.6
D1	SET1701580	SEAT, ASSY, FOLD FLAT SOLID SEAT PAN, BLACK VINYL, 16 X 16-18 W/PRIDE LOGO, 60" LAP BELT, (VISCOUNT: PR-105)	1	1318.16
E1	SET1701581	SEAT, ASSY, FOLD FLAT SOLID SEAT PAN, BLACK VINYL, 18 X 18-20 W/PRIDE LOGO, 60" LAP BELT, (VISCOUNT: PR-106)	1	1318.16
F1	SET1701582	SEAT, ASSY, FOLD FLAT SOLID SEAT PAN, BLACK VINYL, 20 X 20-22 W/PRIDE LOGO, 60" LAP BELT, (VISCOUNT: PR-107)	1	1385.62
G1	ARM156696	ARMREST ASSEMBLY,LEFT/RIGHT,(VISCOUNT PR-15)	2	463.3
H1	SET160081	HEADREST ASSEMBLY,(VISCOUNT:H-PR15)	1	206.46
*I1	ACCBELT1002	LAP BELT,60",C-1010-336	1	111.42
*J1	HDW160283	HARDWARE,ASSY,VISCOUNT SEAT,JAZZY AIR	1	68.82
*K1	HDW125250	SCREW,SOCKET HEAD,M8X1.25 X 20,ISO 4762,CL 12.9,BZ,PVTP	4	6.04

Notes:

Items marked with "*" are not shown.









REF#	PART#	DESCRIPTION	QTY	11i MSRP (AU\$)
A1	DGN160789	TOP SHROUD ASSY,ONYX BLACK,JAZZY AIR	1	254.74
B1	DGN162570	TOP SHROUD ASSY,SAPPHIRE BLUE,JAZZY AIR	1	254.74
C1	DGN167718	TOP SHROUD ASSY,GOLD,JAZZY AIR	1	254.74
D1	DGN167716	TOP SHROUD ASSY,AMBER ORANGE,JAZZY AIR	1	254.74
E1	DGN167715	TOP SHROUD ASSY,ROSE QUARTZ,JAZZY AIR	1	254.74
F1	DGN1611019	TOP SHROUD ASSY,RUBY RED,JAZZY AIR	1	254.74
G1	DGN167719	TOP SHROUD ASSY,SILVER,JAZZY AIR	1	254.74
H1	DGN160790	TOP SHROUD ASSY,PEARL WHITE,JAZZY AIR	1	254.74
I1	DGN167717	TOP SHROUD ASSY,CITRINE YELLOW,JAZZY AIR	1	254.74
J1	DGN1702393	TOP SHROUD ASSY,PINK TOPAZ,JAZZY AIR	1	254.74
K1	DGN153478	CONTROLLER COVER	1	30.28
L1	DGN160107	SHROUD ASSY,REAR,W/LIGHTS,JAZZY AIR	1	146.96
M1	DGN160106	SHROUD ASSY,FRONT,W/LIGHTS,JAZZY AIR	1	217.8
N1	DGN160108	SHROUD,ASSY,BATTERY COVER,LEFT,JAZZY AIR	1	41.7
01	DGN160109	SHROUD,ASSY,BATTERY COVER,RIGHT,JAZZY AIR	1	41.7
P1	PLSKNCP1070	THUMB SCREW,YELLOW,PLASTIC,M8-1.25 X 30MM,CUSTOM,C-1502-069	1	6.02






REF#	PART #	DESCRIPTION	QTY	MSRP (AU\$)
A1	ACC1610775	CANE/CRUTCH HOLDER ASSY,W/ACCESSORY TUBE AND	1	893.92
		HEADREST BRACKET, JAZZY AIR		
B1	MEC165273	ACCESSORY TUBE WELDMENT, JAZZY AIR	1	450.04






REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	ACC1610773	OXYGEN HOLDER ASSY,W/ACCESSORY TUBE AND HEADREST BRACKET,JAZZY AIR	1	913.82
B1	MEC165273	ACCESSORY TUBE WELDMENT, JAZZY AIR	1	450.04









REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	ACC1610772	REAR BASKET ASSY,W/ACCESSORY TUBE,JAZZY AIR	1	696.6
B1	PINCOTD1006	PIN,COTTERLESS DETENT W/RING,.375" DIA X 2.00" GRIP LENGTH,B-1005-169	1	11.94
C1	MEC165273	ACCESSORY TUBE WELDMENT, JAZZY AIR	1	450.04








REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	ACC1610776	WALKER HOLDER ASSY,W/ACCESSORY TUBE,JAZZY AIR	1	877.56
B1	MEC165273	ACCESSORY TUBE WELDMENT, JAZZY AIR	1	450.04






When applicable, assemblies are grouped by color. All components with that color are included in the assembly.

REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	MEC1611900	JAZZY AIR DOCKING ASSY,COMPLETE	1	N/A





When applicable, assemblies are grouped by color. All components with that color are included in the assembly.

A1 Wheel Assy

REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
A1	WHLASMB2045	WHEEL,ASSY,FLAT FREE,SILVER RIM/BLACK TIRE,10",D-1203-065	1	210.18



REF#	PART#	DESCRIPTION	QTY	MSRP (AU\$)
*N/A	INFMANU3992	MANUAL, CONSUMER SAFETY GUIDE, POWERCHAIR	1	5.84
*N/A	INFINFB3072	MANUAL, PRODUCT INSERT SHEET, POWER CHAIR WARRANTY, AUSTRALIA	1	N/A
*N/A	INFMANU4722	MANUAL,OWNER'S,JAZZY AIR,ENGLISH	1	28.66
*N/A	INFINFB3432	MANUAL, PRODUCT INSERT, SPECIFICATIONS, JAZZY AIR, ENGLISH	1	N/A
*N/A	INFINFB3503	MANUAL, PRODUCT INSERT SHEET, VR2, ENGLISH	1	N/A



REF#	PART#	DESCRIPTION	QTY	11 MSRP (AU\$)
*N/A	DWR1234L799	DECAL,ICON,1" X 2.5",WEATHER/MANUAL/LOOSE GARMENT,L- 1234-799	1	6.04
*N/A	DWR1234L848	DECAL,LOGO,3.60" X 0.985",JAZZY,L-1234-848	1	6.14
*N/A	DWR1235L116	DECAL,ICON,TIE DOWN HOOK,L-1235-116	4	5.92
*N/A	DWR1235L515	DECAL,WARNING,BATTERY WARNINGS,L-1235-515	1	7.16
*N/A	DWR1235L520	DECAL,CAUTION,EXTENSION CORD/GROUNDING PRONG,L-1235-520	1	7.16
*N/A	DWR1235L521	DECAL,WARNING,CONTROLLER CORD PLACEMENT,L-1235-521	1	6.6
*N/A	DEC8120018	DECAL,NOT RATED FOR OCCUPIED TRANSPORT,POWER CHAIR	2	17.48
*N/A	LAB148325	DECAL,READ OWNERS MANUAL, 19.1 DIAMETER	1	6.04
*N/A	LAB155762	DECAL,FREEWHEEL, JAZZY AIR	2	6.04
*N/A	LAB155754	DECAL,AESTHETIC, JAZZY AIR	2	6.04
*N/A	LAB155755	DECAL,DATA PLATE, JAZZY AIR(US)	1	6.04

Notes:

Not applicable.



Pride® Mobility Products Corporation

182 Susquehanna Avenue Exeter, PA 18643 USA

Pride® Mobility Products Company

5095 South Service Road Beamsville, ON LOR 1B3 Canada

Pride® Mobility Products Ltd.

32 Wedgwood Road Bicester, Oxfordshire OX26 4UL UK

Pride® Mobility Products Australia Pty. Ltd.

20-24 Apollo Drive Hallam, Victoria 3803 Australia

Pride® Mobility Products New Zealand Ltd.

38 Lansford Crescent Avondale Auckland, New Zealand 1007

Pride® Mobility Products Italia S.r.l.

Via del Porgresso - ang. Via del Lavoro Loc. Prato della Corte 00065 - Fiano Romano (RM)

Pride® Mobility Products Europe B.V. (Authorised EU Representative)

Castricummer Werf 26 1901 RW Castricum The Netherlands

Pride® Mobility Products Spain

Calle Velazquez 80 6D 28001 Madrid Spain

Pride® Mobility Products France

26 rue Monseigneur Ancel 69 800 Saint-Priest

INFMANU4765/Rev.E/031318