

# KARMA

## *Karma 5000 Reclining Wheelchair*



The KM5000 standard lightweight self-propelling recliner is a highly versatile folding wheelchair with an adjustable backrest that allows the carer to continuously adjust the angle from sitting at 90-degrees to a more relaxed 180-degrees making it a great option for stroke patients, the elderly and those recovering from accidents or surgery.



**ALLIED  
MEDICAL  
LIMITED**

## STANDARD FEATURES

1. Reclining backrest
2. Headrest with movable cushion to support head
3. Quick-release, removable rear wheels for ease of storing
4. Fully removable full length armrests, to allow for easy transfers
5. High mounted, push to lock brakes
6. Elevating swingaway legrests with calf pad
7. Reinforced cross bracing - strong but easily folded
8. Multiple-position rear wheels and various castor mounts to optimise centre of gravity and seat-to-floor height adjustments
9. Pump and tool kit for ongoing standard maintenance



## ACCESSORIES

- ◆ Pelvic positioning belt
- ◆ Stump supports
- ◆ 70-degree legrests
- ◆ Angle adjustable footplates
- ◆ Metal footplates (standard on 20")
- ◆ Spokeguards
- ◆ 20" and 22" spoked wheels
- ◆ 24" mag wheels
- ◆ 6" x 1" solid castors
- ◆ 8" x 2" pneumatic castors
- ◆ Heel loops
- ◆ Foot (toe) straps
- ◆ Vinyl coated pushrims
- ◆ Capstan pushrims

## SPECIFICATIONS



**Seat width**

16", 18", 20"



**Seat depth**

16", 18"



**Rear wheel**

24" x 1 3/8" pneumatic with anodised pushrims



**Front castor**

8" x 2" solid



**Seat to floor height**

20"  
(Can be lowered using different sized wheels & castors)



**Wheelchair weight**

Self Propel 20 kg  
(With all standard features)



**User weight limit**

16" - 18" width frames  
115 kg  
20" width frame  
140 kg