

karma

AGILE

Simple but never ordinary, the Agile is built with aircraft-grade aluminium, making it lightweight and durable. With advanced engineering, this chair is easy to push and manoeuvre, along with being easily transportable.



**ALLIED
MEDICAL**
LIMITED



Lightweight, durable and easy to manoeuvre

PERFORMANCE, COMFORT AND STYLE

The Karma Agile Wheelchair is all about performance, comfort, and style. Improved from the tested Karma ECON series, the Soma Agile not only inherited its multi-adjustability but also evolved to provide a wider variety of options to serve individual needs. The Soma Agile is the new standard for transferring wheelchairs and a cost effective option for anyone who is unable to walk but maintains upper body control.

FEATURES

- 8" castors and quick release 24" rear wheels
- Standard flip-back armrests making moving from side to side easy and quick
- 80° removable footrests to ensure a safe transfer into and out of the wheelchair
- The Agile wheelchair is also designed with the carer in mind. With advanced engineering, this chair is very easy to push and manoeuvre
- Fixed backrest
- Compact and lightweight for easy transport
- Quick-release rear wheels allow for even lighter carrying weight and smaller storage space
- Built with aircraft-grade aluminium
- Elevating legrests and anti-tippers available for purchase separately
- Transit kit available for purchase on request
- Available in Pearl Silver with black upholstery
- Three Year Warranty

ACCESSORIES

- Pelvic Belt
- O2 Holder
- Cane Holder
- IV Pole

SPECIFICATIONS



Seat Width
16", 18", 20"



Seat Depth
16", 18"




Wheels
8" Castors
24" Quick Release
Wheels



**Seat to Floor
Height**
19.9"



Overall Width
16" Width: 605 mm
18" Width: 655 mm
20" Width: 705 mm



Folded Width
280 mm



Back Height
16" Depth: 430 mm
18" Depth: 510 mm



Overall Length
16" Depth: 1065 ~
1070 mm



Overall Height
18" Depth: 1115 ~
1120 mm



Total Weight
Approx 13.8 kg



User Weight Limit
125 kg



**ALLIED
MEDICAL
LIMITED**

0800 31 61 81
www.alliedmedical.co.nz

May 2024. This information is correct at the time of printing.
From time to time, product specifications are changed by manufacturers without prior notice.

*Rear wheel size.