

karma



✓ **Lightweight**

✓ **Ergonomic Back and Seat**

Ergo 3

TheNextGeneration of Lightweight Mobility


Karma's new Ergo Series wheelchairs feature the innovative 'Double S-Ergo' seat and backrest system, which offers unparalleled ergonomic support. This patented design ensures a stable, upright posture, keeping users comfortably aligned with their backs straight and heads held high.





World-Renowned S-Ergo Seat

Follows the body's anatomic curves to evenly distribute weight from the hip bone to the lower thigh, maximizing comfort and encouraging proper posture.



Intuitive Design

Designed with both users and caregivers in mind, all movable parts are



Flip-Up Armrests



Detachable Swing-Away Footrests

Relaxed Sitting, Healthy Living

Karma's innovative S-Ergo Backrest enhances both comfort and stability, allowing users—especially those with reduced muscular endurance—to maintain proper, healthy posture with ease.



PATENTED

Ergo Lumbar Support

Karma's Ergo Series lightweight wheelchairs include a set of adjustable curved inserts that provide direct lumbar support and maximize comfort.



Multi-Functional Side Guards

Discreet, easily accessible pockets are elegantly integrated into the inner sides of the side guards for added convenience, while the outer sides feature reflective accents to enhance visibility and safety outdoors.



Lightweight Innovation

Despite its many innovative features, Ergo 3 is one of the lightest wheelchairs on the market, with a total weight of just 11.9 kg (about 26 pounds) in its base configuration.



Comfort in Style



Color Options



Accessories

- Ergo Anti-Tippers
- Basic Anti-Tippers
- Pelvic Belt
- Cane Holder
- Oxygen Tank Holder
- IV Pole



PATENTED

Ergo Anti-Tippers

Ergo 3 Specifications

Seat Width (mm)	410					460					510				
Seat Depth (mm)	435					435					435				
Front Caster Diameter (in)	6		7			6		7			6		7		
Rear Wheel Diameter (in)	14	20	14	20	24	14	20	14	20	24	14	20	14	20	24
Seat Surface Height at Front Edge (mm)	435		485			435		485			435		485		
Backrest Height (mm)	435					435					435				
Footrest to Seat Distance (mm)	385 - 460					385 - 460					385 - 460				
Seat Plane Angle (°)	6					6					6				
Seat to Backrest Angle (°)	104					104					104				
Max. User Weight. (kg)	125*					125*					125*				
Min. Transport Weight (kg)	Q 8.0 F 10.1	Q 8.2 F 11.3	Q 8.4 F 10.5	Q 8.4 F 11.5	Q 8.4 F 12	Q 8.3 F 10.4	Q 8.5 F 11.6	Q 8.7 F 10.8	Q 8.7 F 11.8	Q 8.7 F 12.3	Q 8.6 F 10.7	Q 8.8 F 11.9	Q 9 F 11.1	Q 9 F 12.1	Q 9 F 12.6
Min. Total Mass (kg)	11.9	13.1	12.3	13.3	13.8	12.2	13.4	12.6	13.6	14.1	12.5	13.7	12.9	13.9	14.4
Min. Overall L * W * H (mm)	1050 x 575/595 x 860 1050 x 575/595 x 910					1050 x 625/645 x 860 1050 x 625/645 x 910					1050 x 675/695 x 860 1050 x 675/695 x 910				
Min. Folded L * W * H (mm)	Q 770 x 245 x 660/F14 770x 285 x 700 F20 780x 310 x 700/F24 820x 310 x 700					Q 770 x 245 x 700/F14 770x 285 x 735 F20 780x 310 x 735/F24 820x 310 x 735					Q 770 x 245 x 740/F14 770x 285 x 780 F20 780x 310 x 780/F24 820x 310 x 780				

*The Ergo 3 is rated for users up to 125 kg. However, when equipped with Push & Pull attendant brakes, the maximum user weight is reduced to 115 kg. Karma reserves the right to modify information herein without further notice. " / " means optional : " - " means adjustable ±1cm ± 0.5kg
Rev.01 Sep.2024

