

Customer Details

Client						
Assessor						
Organisation						
	ACC		DSS		Other	
User Weight						

Karma Flexx Adapt Wheelchair

Please note: Scripts are not required for Karma Flexx Adapt wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

Standard Specifications:

- Frame Colour: Apple Green
- Standard Seat Depth Plate
- Seat Plate Angle: 5°
- Front Seat Height: 18"
- Adjustable Back Angle - tool adjustable 90-100°
- Tilt Angle: up to 45°
- Height-adjustable T-Bar Armrests with Desk Armpad
- Swingaway Legrests with Quick-release Knob
- Plastic Flip-up Footplates with Adjustable Base: Size Medium
- Rear Wheels: 16" Quick Release PU Tyres with Drum Brakes
- Height-adjustable Anti-tippers with Foot Release
- Castors: 6" PU Tyres
- Aluminium Height-adjustable Push Handle
- 75kg User Weight Limit

Seat Size (Width x Depth)

	12 x 12		12 x 14
	14 x 14		

Leg Support

	Left	Right	
80° Footrest Hanger			Std
Articulating Elevating Legrest Hanger			\$

Calf Support Pad

	Left	Right	
None			Std
6" x 3" (140 x 75 mm)*			\$
8" x 5" (200 x 120 mm)*			\$

*Required for articulating elevating legrests

Cushions

Varilite Reflex	\$
Vicair Adjuster O2	\$
Vicair Active O2	\$
Vicair Vector O2	\$

Front Seat Height

6" Caster

	17"
	18"

8" Caster

	19"
--	-----

Seat to Floor Height

17" STF - Set in Position 1	\$	18" STF - Set in Position 3	Std
19" STF - Set in Position 3	\$		

Back Support

None	Std
ICON Back Support (please specify size below)	\$
Alternative Back Support (please specify below)	\$

Bodypoint Positioning Belts

2pt	\$	Other (please specify below)
4pt	\$	

Notes