

# Karma Flexx Wheelchair Script



Customer Details					
Client					
Assessor					
Organisation					
	ACC		Whaikaha		Other
User Weight					

## Karma Flexx Wheelchair

Please note: Scripts are not required for Karma Flexx wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

### Standard Specifications:

- 7" Castor Wheels
- 19" Seat to Floor Height
- 1" Seat Dump
- Height Adjustable Ergo Flipback Arms
- 80° Swing-away Leg supports
- Flip-up Foot Supports
- Push to Lock Wheel Locks
- Swing Away Anti-Tippers
- Calf Strap
- Standard Hip Belt
- Seat widths 14-18": 130kg max user weight
- Seat widths 20"+: 170kg max user weight

## Wheelchair Type

Self Propelling

- Quick Release 24" Pneumatic Spoked Rear Wheel
- Ergo Handrim

Transit

- Fixed 14" Solid Tyre Plastic Rim Rear Wheels

## Seat Size (Width x Depth)

14 x 16
16 x 16
16 x 18
16 x 20
18 x 16
18 x 18
18 x 20

20 x 18 <sup>^</sup>
20 x 20 <sup>^</sup>
22 x 18 <sup>^</sup>
22 x 20 <sup>^</sup>
24 x 18 <sup>^</sup>
24 x 20 <sup>^</sup>

<sup>^</sup> Heavy Duty Model

## Front Seat Height

		Seat to Floor Height*					
Castor Size	1) 5"	15.5"	16"	17"	18"		
	6"			17"	18"	19"	
	7"			17.5"	18"	19"	
	8"					19"	20"

<sup>1)</sup> 5" Castor not available on models 20" and wider  
<sup>0.5"</sup> Increments available on fork setting

## Back Cane Height

		Back Height					
Fixed Height	Std			17"			*
Height Adjustable (Adjustable with tools only)	\$	15"	16"	17"	18"	19"	20"

\* Standard on 16 x 20 and 18 x 20 models

## Rear Seat Height

		Rear Seat to Floor Height					
Rear Wheel Size	14"					18"	
	20"	14"	15"	16"	17"	18"	
	22"		15"	16"	17"	18"	19"
	24"			16"	17"	18"	19"

## Anti Tippers

Swing Away Height Adjustable

Rear Seat to Floor Height 17-20".....Std
Rear Seat to Floor Height 14-16".....\$

Pneumatic tyres       Solid tyres

## Centre of Gravity

Rear       Centre       Forward

**Seat Upholstery**

Standard (Padded) .....	Std
Standard with velcro* .....	\$
Omit Seat Upholstery .....	N/C
Cushion (please specify below) .....	\$
Spare Cushion Cover .....	\$
Cushion	

\*Available on most sizes. Some sizes may be indent only.

**Transit Option**

Transit Tie Down Option.....	\$
------------------------------	----

\* Not available on HD

**Leg Support**

	Left	Right
80° Swing-away .....		.....Std
90° Swing-away .....		.....\$
Elevating .....		.....\$
Stump Support .....		.....\$

**Foot Support**

	Left	Right
Standard Flip-up .....		.....Std
Flip-up Angle Adjustable*.....		.....\$
Flip-up Angle & Depth Adj. Small		.....\$
Large (16" or wider chairs only)		.....\$

\* Not available on HD

**Calf Strap/Heel Loops**

Calf Strap .....	Std
Heel Loops* .....	\$

\*Standard on HD with Flip-up Foot Supports only, otherwise an upcharge.

**Attendant Brake**

	Left	Right
Self Propelling 17" Height		.....\$
Self Propelling 20" Height		.....\$
Transit 17" Height		.....\$
Transit 20" Height		.....\$

**Other Parts Required**

**Back Support**

Standard (Padded) .....	Std
Tension Adjustable (not available on height adj back canes) .....	\$
Omit Back Upholstery .....	N/C
Alternative Back Support (please specify below) .....	\$
Alternative Back Support	

**Arm Support**

Flip-back Ergo (desk length) Height Adjustable .....	Std
Flip-back Boat (full length) Height Adjustable .....	\$

**Wheel Locks**

Push to Lock with Brake extension .....	Std
Pull to Lock* .....	\$

\* Available only on 20, 22 and 24" Rear wheels

**Accessories**

Clear Spoke Guards .....	\$
Oxygen Tank Holder .....	\$
Push Handle Rigidiser Bar* .....	\$
IV Pole .....	\$

\* Std on height adj back canes 18" and wider and HD chairs

**Head Support**

Foldable Medium (14 - 17" widths) .....	\$
Foldable Large (17 - 20" widths) .....	\$
Other .....	\$

**Pelvic Positioning Belt**

Lap Belt .....	Std
2 Point .....	\$
4 Point .....	\$

**Tech Notes**