

Karma Flexx Wheelchair Script



Customer Details					
Client					
Assessor					
Organisation					
	ACC		DSS		Other
User Weight					

Karma Flexx Wheelchair

Please note: Scripts are not required for Karma Flexx wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

Standard Specifications:

- 7" Castor Wheels
- 19" Seat to Floor Height
- 1" Seat Dump
- Height Adjustable Ergo Flipback Arms
- 80° Swing-away Leg supports
- Flip-up Foot Supports
- Push to Lock Wheel Locks
- Swing Away Anti-Tippers
- Calf Strap
- Standard Hip Belt
- Seat widths 14-18": 130kg max user weight
- Seat widths 20"+: 170kg max user weight

Wheelchair Type

Self Propelling

- Quick Release 24" Pneumatic Spoked Rear Wheel
- Ergo Handrim

Transit

- Fixed 14" Solid Tyre Plastic Rim Rear Wheels

Seat Size (Width x Depth)

	14 x 16		20 x 18 [^]
	16 x 16		20 x 20 [^]
	16 x 18		22 x 18 [^]
	16 x 20		22 x 20 [^]
	18 x 16		24 x 18 [^]
	18 x 18		24 x 20 [^]
	18 x 20		

[^] Heavy Duty Model

Front Seat Height

		Seat to Floor Height [*]										
		15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	21"
Castor Size	1) 5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
	6"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	7"					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	8"								<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>

^{1) 5"} Castor not available on models 20" and wider

^{2) 20"} Not available on chair widths 20" and wider

Back Cane Height

Fixed Height	Back Height				
	Std				
			<input type="checkbox"/>		<input type="checkbox"/>
					20"
Height Adjustable (Adjustable with tools only)	\$				
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		15"	16"	17"	18"
				19"	20"

* Standard on 16 x 20 and 18 x 20 models

Rear Seat Height

		Rear Seat to Floor Height					
		14"	15"	16"	17"	18"	19"
Rear Wheel Size	14"					<input type="checkbox"/>	
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Anti Tippers

Swing Away Height Adjustable

<input type="checkbox"/>	Rear Seat to Floor Height 17-20".....Std
<input type="checkbox"/>	Rear Seat to Floor Height 14-16".....\$

Pneumatic tyres Solid tyres

Centre of Gravity

Rear Centre Forward

Seat Upholstery

Standard (Padded)	Std
Standard with velcro*	\$
Omit Seat Upholstery	N/C
Cushion (please specify below)	\$
Spare Cushion Cover	\$
Cushion	

*Available on most sizes. Some sizes may be indent only.

Transit Option

Transit Tie Down Option	\$
-------------------------------	----

* Not available on HD

Leg Support

	Left	Right	
80° Swing-away			Std
90° Swing-away			\$
Elevating			\$
Stump Support			\$

Foot Support

	Left	Right	
Standard Flip-up			Std
Flip-up Angle Adjustable*			\$
Flip-up Angle & Depth Adj. Small			\$
Large (16" or wider chairs only)			\$

* Not available on HD

Calf Strap/Heel Loops

Calf Strap	Std
Heel Loops*	\$

*Standard on HD with Flip-up Foot Supports only, otherwise an upcharge.

Attendant Brake

	Left	Right	
Self Propelling 17" Height			\$
Self Propelling 20" Height			\$
Transit 17" Height			\$
Transit 20" Height			\$

Other Parts Required

Back Support

Standard (Padded)	Std
Tension Adjustable (not available on height adj back canes)	\$
Omit Back Upholstery	N/C
Alternative Back Support (please specify below)	\$
Alternative Back Support	

Arm Support

Flip-back Ergo (desk length) Height Adjustable	Std
Flip-back Boat (full length) Height Adjustable	\$

Wheel Locks

Push to Lock with Brake extension	Std
Pull to Lock*	\$

* Available only on 20, 22 and 24" Rear wheels

Accessories

Clear Spoke Guards	\$
Oxygen Tank Holder	\$
Push Handle Rigidiser Bar*	\$
IV Pole	\$

* Std on height adj back canes 18" and wider and HD chairs

Head Support

Foldable Medium (14 - 17" widths)	\$
Foldable Large (17 - 20" widths)	\$
Other	\$

Pelvic Positioning Belt

Lap Belt	Std
2 Point	\$
4 Point	\$

Tech Notes