

karma

KM7520

Standard Lightweight Paediatric Self
Propelling Wheelchair

A durable, functional wheelchair designed to meet a
range of New Zealand kid's needs



**ALLIED
MEDICAL**
LIMITED



Reclining Function
Also
Available

+ FEATURES

- Easily adjustable push handles – to accommodate carers of all heights
- Quick release, removable rear wheels for ease of storing in vehicles
- Fully removable, height adjustable armrests, that allow for easy transfers
- High mounted, push-to-lock brakes, that can be rear mounted for attendant use
- Reinforced cross bracing - strong but easily folded
- Durable upholstery, with hidden screws
- Multiple-position rear wheel and castor mounts to optimise centre of gravity and seat-to-floor height adjustments
- Pump and tool kit for on-going standard maintenance
- Swing-away 90-degree legrests with calf strap



+ OPTIONAL EXTRAS

- Elevating legrests with calf pads
- Pelvic positioning belt
- Standard foam cushion
- 70-degree legrests
- Solid rear wheels
- Heel loops
- Coated pushrims
- Extended brake handles
- Tension adjustable back upholstery
- Side reducer pads

+ SPECIFICATIONS



Seat Width
11.5" (28 cm)
13.5" (34 cm)
*10.5" available if side reducer pads are used



Seat Depth
11.5" (28 cm)
13.5" (34 cm)
*Customisation available



Rear Wheels
20" x 1 3/8"
Pneumatic




Front Wheels
6" x 1"
Pneumatic



Back Height
16" (41 cm)



Seat to Floor Height
18.5" (46 cm)



Overall Width
20" (51 cm)
When open



Overall Length
32" (81 cm)



User Weight Limit
100 kg



Total Weight
14.5 kg
Including legrests and wheels



**ALLIED
MEDICAL
LIMITED**