

Karma Wheelchair

Quick Reference Guide

User Weight Limits

Model	Seat Widths	User Weight Limit
Agile	16" - 20"	125 kg
Ergo Lite	16" or 18"	100 kg
Flexx Junior	10" - 14"	75 kg
Flexx	14" - 18"	130 kg
Flexx Tall	16" - 19"	130 kg
Flexx HD	20" - 24"	170 kg



Flexx



Flexx Junior



Flexx Tall



Flexx HD

Overall Width

Model	Total Width
Agile	16" Width: 605 mm, 18" Width: 655 mm, 20" Width: 705 mm
Ergo Lite	16" Width: 22" (560 mm) 18" Width: 24" (610 mm)
Flexx Junior	Seat width plus 7.5" or 7.9" (190 or 200 mm)
Flexx / Tall / HD Self-Propel	Self-Propel: Seat width plus 9" Transit: Seat width plus 7"



Agile



Ergo Lite

Agile and Ergo Lite Seat To Floor Height

Model	Rear Wheel	6" Castor	8" Castor
Agile	24"	19"	20"
Ergo Lite	14" (Transit)	19"	N/A

Flexx Range Seat To Floor Height

Rear Wheel	5" Castor	6" Castor	7" Castor	8" Castor
14" (Transit)	18"	18"	18"	N/A
20"	15.5" - 17"	17"	N/A	N/A
22"	16" - 18"	17" - 18"	18"	N/A
24"	17" - 18"	17" - 19"	18" - 19"	19" - 21"

(These measurements are a guide only and do not factor in specifics such as seat rake)

Lower Leg Length			
Model	Standard Length Range	Shortest Length Achievable	Longest Length
Agile	16" (406 mm)	15" (381 mm)	19" (483 mm)
Ergo Lite	16" (406 mm)	15" (381 mm)	16" (406 mm)
Flexx Junior	7" - 11" (180 mm - 280 mm)	7" (180 mm)	16" (400 mm)
Flexx / Tall / HD*	15.5" - 19" (393 mm - 483 mm)	15" (381 mm)	19.5" (495 mm)

****Flexx Junior leghangers can be modified to achieve a shorter lower leg length on these models***

Armrest Heights			Armrest Lengths		
Model	Minimum Height	Maximum Height	Model	Desk	Full
Agile	9" (220 mm)		Agile	N/A	12" (305 mm) x 1.75" (45 mm)
Ergo Lite	9" (210 mm)		Ergo Lite	8.75" (223 mm) x 2" (51 mm)	N/A
Flexx Junior	6" (150 mm)	10" (250 mm)	Flexx Junior	9.6" (244.5 mm) x 1.76 (44.7 mm)	
Flexx / Tall / HD	8" (203 mm)	11.5" (292 mm)	Flexx	10" (254 mm) x 1.88" (47 mm)	15" (381 mm) x 1.88" (47 mm)
			Flexx Tall / HD	11.5" (292 mm) x 2" (50 mm)	15" (381 mm) x 2" (50 mm)

Backrest Heights		
Model	Standard Height (Seat to top of back upholstery)	Adjustment Available
Agile	16" Depth: 17" (430 mm) 18" Depth: 20" (510 mm)	N/A
Ergo Lite	16" (420 mm)	N/A
Flexx Junior	14.5" (370 mm)	14.5" - 17" (370 - 430 mm)
Flexx	17" (432 mm)	15" - 21" (381 - 533 mm)
Flexx Tall	20" (508 mm)	17" - 20" (430 - 508 mm)
Flexx HD	17" (432 mm)	15" - 20" (381 - 508 mm)

Tube Diameter	
Armrest Tube Size	3/4"
Back Cane Tube Size	7/8"

