

Account Number:

Date:

Address:

City:

Phone:

Email:

Quote: PO Number:

Associated Quote/Order#:

Effective: 1/20/2025

Catalyst 5 - Transit Included (K0005) _____ \$

Catalyst 5 - Transit Not Included (K0005) _____ \$

300 lb weight capacity

Frame

Seat Width

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2201)	\$
21" (E2201)	\$
22" (E2201)	\$

Seat Depth

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	\$
20" (E2203)	\$

Transit Option

Transit System (Weight Limit 300lbs)	\$
Not available with 350 lb Weight Capacity	
Omit	N/C

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinerge wheels or Height Adjustable Flip Back T-Arm.

e-motion Ready	N/C
twion Ready	N/C

Titanium Option

Titanium Frame	\$
Includes hangers and base frame. Does not include Cross Braces.	

Weight Capacity

300 lbs	N/C
350 lbs (Heavy Duty)	
Requires minimum 6" caster height	\$

Cross Brace

Standard Cross Braces	N/C
Seating System Cross Braces	N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.	
Titanium Cross Braces	\$
Not available with 21" and 22" Seat Widths	

Seat Pan

Removable Aluminum Seat Pan E2231	
1/2" Shorter than Depth	\$

Rear Frame Type

No selection needed with Titanium Frame

Standard (Curved)	N/C
Select from Standard or Angle Adjustable Back Post on chart on the next page	
Depth Adjustable (K0108)	\$
Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.	
Reclining Frame	\$
Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.	

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0011.43 Pricing and specifications subject to change without notice.

Backrest

Depth Adjustable not available with Titanium Frame

Back Post	Price	Height Range	Back Height				
Standard	N/C	Short ^{3,4}	8.5"	9.5"	10.5"	11.5"	
		Medium ⁴	12.5" ⁶	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	
Angle Adjustable w/ Push Handle	\$	Short ⁴	14"	15"	16"		
		Medium ⁴	16"	17"	18"		
		Tall	18"	19"	20"		
Stroller Handle Back Post (Depth Adjustable)	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable) ^{2,9}	\$	Short ¹	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" ⁸	21" ⁵	22" ⁵	23" ⁵	24" ⁵
Reclining w/ Adjustable Handle (E1225) ^{1,7}	N/C	Short	22"				
		Tall	26"				

¹ Not available with Dynamic Back

² Not available with Standard Back Upholstery

³ Not available with Tension Adjustable Back Upholstery

⁴ Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths

⁵ Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

⁶ Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

⁷ Not available with O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

⁸ Not available with Angle Adjustable Locking Extendable Flip Up Armrest

⁹ Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

Reclining Back Options

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

Standard Spring	N/C
Heavy Duty Spring	N/C

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

Back Options

Straight with Push Handle	N/C
Straight without Push Handle	N/C
8 Degree Bend with Push Handle	\$
Not available with Short Back Height Range	
Straight with Fold Down Push Handle	\$
Not available with 350 lb Weight Capacity	

Stroller Handle

Caregiver Adjustable Height	\$
Not available with Swing away armrests. Available only with straight without push handle backrest.	

Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity

Low Profile Push Handle	\$
2" shorter than standard push handles	

Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

Stroller Handle Extensions	\$
Not available with Height Adjustable Straight with Push Handle Back Post	
Dynamic Back (E2398)	\$
Available only with T-Arms.	

Back Upholstery Options

Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard	N/C
Tension Adjustable with Standard Cover (E2611)	\$
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$
Omit	\$

U.S. Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$

Frame Type / Footrest

Frame Type		Footrest Hangers							
Swing Away	N/C	Swing Away	Extension Tube Mount	60°	70°	80°	90°		Retail
Super Low - Swing Away	N/C		Front Mount ²			1	1		N/C
70° Fixed Front	N/C		Elevating Leg Rests ^{3,4}						\$
80° Fixed Front	N/C		Pro Elevating Leg Rests E0990 ^{3,4}						\$
			Omit Hangers						\$
		Fixed	Extension Tube Mount			1			N/C

¹ Not available with 7" or 8" casters
² Not available with Titanium Frame
³ Available only with Classic Hanger Release
⁴ Not available in combination with 24" wheels and Push to Lock Flush Mount Wheel Locks.

ELR and PRO ELR Cover

Standard ELR Foam Knee Pad (Pair) K0108	\$
Only available with ELR	
Pro ELR Gel Pad Cover (Pair) K0108	\$
Only available with Pro ELR	

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"			
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"			
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"			
Front Mount ^{1,3} Hangers		Standard Selections									Requires 3" Extension		
	Short Front Mount Hanger (70°, 80°, 90°)	4" ²	4.5" ²	5"	5.5"	6"	6.5"				10"	10.5"	11"
		7"	7.5"	8"	8.5"	9"	9.5"				11.5"	12"	12.5"
	Medium Front Mount Hanger (70°, 80°, 90°)	4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"	7.5"		12"	12.5"	13"
		8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"		13.5"	14"	14.5"
	Long Front Mount Hanger (70°, 80°)	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	14"	14.5"	15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"
ELR	Standard ELR Extension Tube Short (12" - 17")												
	Standard ELR Extension Tube Medium (14" - 19")												
	Standard ELR Extension Tube Long (16.5" - 21.5")												
PRO ELR Low Pivot Setting	Pro ELR Low Pivot Setting Short (12" - 17.5")												
	Pro ELR Low Pivot Setting Long (15" - 21")												
PRO ELR High Pivot Setting	Pro ELR High Pivot Setting Short (10.5" - 15.5")												
	Pro ELR High Pivot Setting Long (13" - 19")												
Fixed Front Frame (70°, 80°)	Short (11" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"			
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"			

¹ Not available with Titanium Frame
² Available with 90° only
³ See Minimum Front Seat Height chart below

Minimum Front Seat Height with Front Mount Hangers					
Hanger Length	Short		Medium		Long
	70°	80/90°	70°	80/90°	70° 80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16" 17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5" 20"

3" Extension

3" Extension	\$
3" Extension (Provided as Parts)	\$

Footplates

Composite	N/C
Not available with Front Mount Hanger	
Composite Angle Adjustable (K0040)	\$
Not available with Front Mount Hanger or 14" Seat Width	
Aluminum Angle Adjustable (K0040)	\$
Aluminum Locking Multi-Angle Adjustable Footplate (K0108)	\$
Available only with Front Mount Hanger	
One Piece Flip Up Angle Adjustable (K0040)	\$
Not available with Front Mount Hanger, ELR, or PRO ELR	

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

Calf Straps

Velcro Adjustable (K0038)	\$
Padded Velcro Adjustable (K0038)	\$
Padded Airmesh Calf Strap (K0038)	\$

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$
Adjustable Heel Loops (E0951)	\$
Omit Heel Loops	N/C

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size				
Left	Left Limb Support (E1020)	\$	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5x8 Pad Left Side	7x10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side
Right	Right Limb Support (E1020)	\$	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5x8 Pad Right Side	7x10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Casters

		Caster Size					
		Price	4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C					
	.75" Lighted Rollerblade	\$					
	1.25" Soft Roll	N/C					
	1" Poly Aluminum	\$					
	1" Pneumatic	\$					
	1.5" Soft Roll	\$					
	1.5" Soft Roll Aluminum	\$					
	2" Poly	\$					
	2" Pneumatic w/Foam Insert	\$					

Caster Housing

Standard Housing	N/C
Performance Housing	\$
<i>Not available with 350 lb Weight Capacity</i>	

Caster Pin Locks

Caster Pin Locks (K0073)	\$
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Wheel Locks

Wheel Lock

Push to Lock	N/C
Push to Lock Flush Mount	N/C
Pull to Lock	N/C
Low Profile Scissor Lock	\$
<i>Not available with low poly or full poly tires</i>	
Drum Brake (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	

Fork Options

Best Fit	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem	N/C
3/4" Longer Fork Stem	N/C
1 1/2" Longer Fork Stem	N/C
Frog Legs Suspension Fork	\$
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
Frog Legs II Suspension System	\$
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	

Frog Legs Options

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C
User Weight 160-220 lbs	N/C
User Weight 220-260 lbs	N/C
User Weight 260-300 lbs	N/C

Wheel Lock Options

Extension Handles (E0961)	\$
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	
Anti-Rollback System	\$
<i>PRATIKO Anti-Rollback System. Installed by Ki. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg). Available with 22" and 24" wheel sizes of Spoke, Maxx Spoke and Maxx Mag only. Not available with transit, drum brake, seat pan, OAD, heavy duty, not available on 21" or wider seat widths.</i>	

Rear Seat Height and Wheels

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height									
20"	Low	11.5" ³	12" ¹	12.5" ³	13" ¹	13.5" ³	14"	14.5" ³	15"	15.5" ³	16" ³
	Full	12" ³	12.5" ¹	13" ³	13.5" ¹	14" ³	14.5"	15" ³	15.5"	16" ³	16.5" ³
22"	Low	12.5" ³	13" ¹	13.5" ³	14" ¹	14.5" ³	15"	15.5" ³	16"	16.5" ³	17" ³
	Full	13" ³	13.5" ¹	14" ³	14.5" ¹	15" ³	15.5"	16" ³	16.5"	17" ³	17.5" ³
24"	Low	13.5" ³	14" ¹	14.5" ³	15" ¹	15.5" ³	16"	16.5" ³	17"	17.5" ³	18" ³
	Full	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
25"	Low	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
26"	Low		15" ^{1,2}	15.5" ³	16" ¹	16.5" ³	17"	17.5" ³	18"	18.5" ³	19" ³
	Full	15" ^{2,3}	15.5" ^{1,2}	16" ³	16.5" ¹	17" ³	17.5"	18" ³	18.5"	19" ³	19.5" ³

¹ Not available with Amputee Axle Plate

² Not available with Height Adjustable Flip Back T-Arm

³ Requires Universal Axle Plate

One Arm Drive

One Arm Drive - Quick Release (E0958) \$
 Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.
 Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.

Drive Side

Left Side Drive N/C
 Right Side Drive N/C

Drive Handrim

Aluminum Drive Handrim N/C
 Poly Coated Drive Handrim \$

Axle Plates

No selection needed with One Arm Drive

Universal N/C
 Amputee (E0959) \$
 Not available with Swing Away Armrests or Wheelchair Drive Ready Option.
 Required with Reclining Backrest.

Handrims

Aluminum Anodized N/C
 Superlight Aluminum Anodized \$
 Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinerig wheels. Not available with One Arm Drive.
 Plastic Coated \$
 Projection (8 Vertical Projections) \$
 Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels. Not available with One Arm Drive.

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used.

Close Mount \$
 Reduces overall width by .75"

Rear Wheel Type

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	\$					
	Maxx Spoke	N/C					
	Mag	N/C					
	Colored Mag ³	\$					
	Superlight Spoke ¹	\$				²	
	Spinerig Spox	\$					
	Spinerig LX	\$					

¹ Not available with 350 lb Weight Capacity

² Available only with Superlight and Plastic Coated Handrims

³ Color selection on page 8

Spinerig Options

Select spoke color for Spinerig Wheels

Black Spokes Red Spokes Blue Spokes White Spokes N/C

Ergonomic Handrims

Not available with 20" or 26" wheels

Lite Regular Grip \$
 Standard Regular Grip \$
 Lite Extra Grip \$
 Standard Extra Grip \$
 Lite Omit Grip \$
 Standard Omit Grip \$

3M Micro Grip \$
 Not available on Plastic Coated or Projection Handrims

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" (E2211, E2212)	\$
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	
Schwalbe Downtown 1 3/8" (E2211, E2212)	\$

Axles

No selection needed with One Arm Drive or Drum Brake

Quick Release	N/C
Quad Quick Release	\$
<i>Not available with Mag or Spinergy wheels</i>	

Low Profile

Low Poly 1" (E2218)	N/C
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.</i>	
Sentinel Solid Tire 1" (E2218)	\$
<i>Not available on Superlight Wheels</i>	
High Pressure 1" (E2211, E2212)	\$
<i>Not available on Mag wheel</i>	
Schwalbe One 1" (E2211, E2212)	\$
<i>Not available on Mag wheel, 20", 22" or 26" wheel sizes</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	

Axle Position

Rear Wheel Camber

2° and 4° Not available with Heavy Duty option

0°	2°	4°	N/C
----	----	----	-----

Center of Gravity

Rear	Center	Forward	N/C
------	--------	---------	-----

No selection needed with Amputee Axle Plate

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$
	Full	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 ²	Desk	\$
	Full	
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 ²	Desk	\$
	Full	

1 If requested configuration is unavailable, Standard Height T-Arm will be used

2 Available only with Curved Rear Frame

3 Not available with Seat Depths less than 16"

4 Available only with Standard Back Posts

5 Available only with Depth Adjustable and Reclining Back Posts

6 Not available with Arm Pads or Half-Folding back

7 Not Available with 15" Rear Seat Height and 24" wheels

Armrest Type	Length	Price
Flip Back (10.5") ^{3,4,7}	Desk	N/C
	Full	
Angle Adjustable Locking Extendable Flip Up K0020 ⁵	Desk	\$
	Full	
Swing Away (10.5"-12.5") E0973 ^{4,6}		\$
Swing Away with Side Guard (10.5"-12.5") E0973 ^{4,6}		\$
Omit Armrests		N/C

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guard

Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests

Adult Composite Side Guards (K0108)	\$
Pediatric Composite Side Guards (K0108)	\$
Aluminum Straight - Black Matte (K0108)	\$
Aluminum Straight - Match Frame Color (K0108)	\$
<i>Not available with Hydrographics or Titanium Frame</i>	
Aluminum Adjustable Fender w/Foam Pad - Black Matte (K0108)	\$
Aluminum Adjustable Fender w/Foam Pad - Match Frame Color (K0108)	\$
<i>Not available with Hydrographics or Titanium Frame</i>	
Carbon Fiber Adjustable Fender (K0108)	\$

Anti-tips

Rear Anti-Tips (E0971)	\$
Swing-In Anti-Tip Single - Left (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with O2 Holder.</i>	
Omit	N/C

Accessories

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
Small Flex-Mount Positioning Belt (16-24") (E0978)	\$
Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$
Large Flex-Mount Positioning Belt (25-37") (E0978)	\$
Belt Mounting Band Clamp	\$
Clamp includes hardware to mount Belt	

Additional Accessories

Tool Kit	\$
Removable Underseat Pouch	\$
Pediatric Ki Mobility Backpack	\$
Ki Mobility Backpack	\$
Luggage Carrier	\$
Not available with ELR or Pro ELR	
Black Neoprene Impact Guards	\$
Clear Spoke Guards (K0065)	\$
Black Spoke Guards (K0065)	\$
Cane and Crutch Holder	\$
IV Pole	\$

O2 Holder

Not Available with Depth Adjustable Rear Frame

O2 Holder - Left Side (E2208)	\$
O2 Holder - Right Side (E2208)	\$

Colour

Touch-up paint is available for purchase through the parts manual

Options	Price	Colours																			
		Gloss										Matte									
		Black	Ruby Red	Radiant Orange	Yellow	Electric Green	Rich Teal	Glacier Blue	Blue Steel	Deep Sapphire Blue	Purple Plum	Opalescent Violet	Dark Pink	Light Pink	Champagne	Shadow Black	Iridescent Navy	Bright Red	Tactical Green	Dry Earth Tan	Platinum
Frame	N/C																				
Mag Colour																					
Colour Accent Package ¹	N/C																				
Colour Anodizing Package ²	Black	Blue	BrightGreen	Orange	Pewter	Pink	Purple	Red	Gold												
	N/C	\$	\$	\$	\$	\$	\$	\$	\$												

¹ Includes axle plates and forks. Universal axle plate only. Not available with amputee axle plate, or one-arm drive.

² Select desired color. If no selection is made, black will be used. Includes standard caster housing, caster fork, and axle plates. Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System. Not available with amputee axle plate. Extended lead time may apply.

Hydrographics Options

Not available with Frame Color. Hydrographics Base Color required with Hydrographics. Zebra not available with black.

Options	Price	Hydrographics						Hydrographics Base Colour					
		AOR1 (Desert Digital Camo)	NEXT Camo	Pink Camo	Carbon Fiber	Muddy Girl	Zebra	Black	Blue	Granny Smith	Pink	Purple	Red
Frame	\$												
Base Colour	N/C												

Catalyst Front Seat Height Matrix																															
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster															
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork												
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"												
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4													
15.5 ¹	14 ¹			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5			20.5	n/a	1 1/2													
16	14.5	3/4		16	14.5	Std	5"	17	15.5	Std	6"	19	18	1 1/2		19	17.5	Std	7"												
16.5 ¹	15 ¹			16.5	15			17.5	16			19.5	18																		
17	15.5			16.5	15.5			18	16.5			20	18.5			20	18.5														
15.5	14	Std	17	16	3/4	17.5		16.5	3/4	18.5		17	Std	20.5	19																
16	14.5		17.5	16		18		17		19		17.5		19.5	18.5																
16	15	3/4	18	16.5		1 1/2		18.5		17.5		Std		19.5	18	20	19														
16.5	15.5		16.5	15	17			15.5	19	17.5				19.5	18	20.5	n/a														
17	15.5	1 1/2	17	16	Std	6"		18.5	17	1 1/2		7"	19	18	3/4	21	n/a	3/4													
17.5	16		17.5	16				19.5	18				20.5	n/a																	
16.5	15	17.5	16	19			17.5	20	18.5		21		n/a																		
16.5	15	Std	17.5	16	3/4		18	16.5	Std	20	19		1 1/2	21	n/a	1 1/2															
17	15.5		18	16.5			18.5	17		20.5	n/a																				
17	16	3/4	18.5	17	1 1/2		7"	19	17.5	Std	7"			20	18.5	1 1/2															
17.5	16.5		18	16.5				19.5	18				20.5	19																	
18	16.5	18.5	17	18.5				17.5	21				n/a																		
18.5	17	1 1/2	19	17.5	3/4			7"	19	18			3/4	7"																	
18	16.5		18	16.5		19.5			18.5																						
18.5	17	18.5	17	20		19																									
18.5	17.5	3/4	19	17.5	1 1/2	7"			19.5	18		1 1/2	7"																		
19	18		18.5	17.5					20	18.5																					
19.5	18.5	19	18	20.5					19																						
20	19	1 1/2	19.5	18.5	3/4		7"		21	n/a	1 1/2	7"																			
				19.5					18	6" x 2" caster not available with 5" fork																					
				20					18.5																						
				20.5	19																										

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Cushions and Backrests

AXIOM

Ki Mobility

ADI
A Stealth Product

VARILITE

ROE

Cushion:

Backrest:

Other Accessories: