

Account Number:

Date:

Address:

Quote: PO Number:

City:

Phone:

Email:

Associated Quote/Order#:

Effective: 6/17/2024

Focus CR - Transit Included (E1161) \$

Focus CR - Transit Not Included (E1161) \$

300 lb weight capacity Tilt range: -5° anterior to 50° posterior

Frame

Frame Width

Select Frame Width within desired Frame Width Range

Frame Width	Frame Width				
Small (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2201) \$
Large (18" - 22")	18" N/C	19" N/C	20" (E2201) \$	21" (E2201) \$	22" (E2201) \$
X-Large ¹ (20" - 24")	20" (E2201) \$	21" (E2201) \$	22" (E2201) \$	23" (E2201) \$	24" (E2202) \$

¹ 400 lb Weight Capacity is required

Seating Setup

Make one selection from either Option 1 at right or Option 2 below based on desired seating setup

Setup By Ki Mobility (Option 1) _____ N/C

Select Seat Depth corresponding to desired Frame Depth Range row and Backrest Offset column. Note: All Seat Depths in a column and within the same frame depth range are achievable with future adjustment.

Setup by Frame Depth (Option 2) _____ N/C


Select Frame Depth within desired Frame Depth Range

Option 1: Setup By Ki Mobility

Frame Depth Range	Backrest Offset A						
	-2"	-1"	Flush	+1"	+2"	+3"	+4"
Seat Depth							
Short	16"	15"	14"				
	17"	16"	15"	14"			
	18"	17"	16"	15"	14"		
	19"	18"	17"	16"	15"	14"	
Medium	20" ¹	19"	18"	17"	16"	15"	14"
	21" ¹	20" ¹	19"	18"	17"	16"	15"
	22" ²	21" ¹	20" ¹	19"	18"	17"	16"
	23" ¹	22" ²	21" ¹	20" ¹	19"	18"	17"
Long	24" ¹	23" ¹	22" ¹	21" ¹	20" ¹	19"	18"
	25" ¹	24" ¹	23" ¹	22" ¹	21" ¹	20" ¹	19"
	26" ²	25" ¹	24" ¹	23" ¹	22" ¹	21" ¹	20" ¹
	27" ²	26" ²	25" ¹	24" ¹	23" ¹	22" ¹	21" ¹

¹ Price is \$630 E2203

² Price is \$1056 E2204



Option 2: Setup By Frame Depth

Frame Depth Range	Frame Depth				
Short	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Medium	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$
Long	18" N/C	19" N/C	20" (E2203) \$	21" (E2203) \$	22" (E2204) \$

Seat Options

Growing Adjustable Aluminum Seat Pan _____ \$

Omit Seat Pan _____ N/C

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0031.43 Pricing and specifications subject to change without notice.

Frame Options

Weight Capacity	Offset Option
300 lbs _____ N/C	Offset Hanger with Tilt Stops _____ \$
400 lbs (Heavy Duty) _____ \$	<i>Tilt Stop included. No anterior tilt.</i>
<i>Available only with Long Frame Depth Range and 8" x 2" or 6" x 2" Casters</i>	
Transit	Vent and Battery Tray
Transit System _____ \$	<i>Requires minimum 16" Seat Height. Not available on 14" or 15" Frame Widths, Attendant Foot Lock, 12" wheels, One Arm Drive, 5" casters</i>
Omit Transit _____ N/C	Vent Tray (E1029) _____ \$
	<i>Not available with Self Propel Axle Plate or Foot Tilt. Requires Anti-Tip.</i>
	Battery Box Holder (E1029) _____ \$
	<i>Available only with Dual Hand Tilt</i>
	Vent Tray and Battery Box Holder (E1029) _____ \$
	<i>Requires Anti-Tip. Available only with Dual Hand Tilt. Not available with Self Propel Axle Plate.</i>

Tilt Mechanism

Tilt Mechanism	Power Tilt Options
<i>Power Tilt is not available with Vent Tray, Battery Box Holder, 23" and 24" Wide, or 400 lb Weight Capacity. Power Tilt includes Anterior and Posterior Tilt Stops.</i>	<i>Power Tilt selection required</i>
Hand Tilt - Right _____ N/C	Toggle Switch _____ N/C
Hand Tilt - Left _____ N/C	Pendant Control _____ N/C
Hand Tilt - Dual _____ N/C	<i>Sent as accessory</i>
Foot Tilt - Right _____ N/C	Interface for Adaptive Switching _____ \$
Foot Tilt - Left _____ N/C	<i>Sent as accessory</i>
Power Tilt - Right (E1002) _____ \$	
Power Tilt - Left (E1002) _____ \$	

Backrest

Backrest	Stroller Handle Back Post Option
Stroller Handle Back Post (Short 20") (Folding) _____ N/C	<i>Available only with Stroller Handle Back Post</i>
Stroller Handle Back Post (Tall 24") (Folding) _____ N/C	Removable Stroller Handle Extension _____ \$
Fixed Height Back with Adjustable Handle (Short 22") (Folding) _____ N/C	
<i>Not available with Dynamic Back</i>	
Fixed Height Back with Adjustable Handle (Tall 26") (Folding) _____ N/C	Other Backrest Options
<i>Not available with Dynamic Back</i>	Non Folding Back _____ N/C
Adjustable Height Back with Adjustable Handle (Short 18" - 25") (Folding) _____ N/C	<i>Standard with Reclining Back and Dynamic Back</i>
<i>Not available with O2 Holder or Dynamic Back</i>	Dynamic Back (E2398) _____ \$
Adjustable Height Back with Adjustable Handle (Tall 20" - 27") (Folding) _____ N/C	<i>Available only with Non Folding Back. Not available with Angle Adjustable Locking Extendable Flip Up Armrest or 400 lb Weight Capacity.</i>
<i>20"-22" Heights not available with Dynamic Back</i>	Adjustable Height Rigidizer Bar _____ \$
Reclining Back with Adjustable Handle (Short 22") (Non-Folding) _____ \$	<i>Not available with 23" or 24" frame widths</i>
<i>Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.</i>	
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding) _____ \$	Reclining Back Options
<i>Not available with O2 Holder or Dynamic Back. If no Recline Angle and Gas Spring selections are made, 0° - 65° and Standard Spring will be used.</i>	<i>Required when Reclining Back is selected</i>
	Recline Angle
	0-65° Recline Angle (E1225) _____ N/C
	10-75° Recline Angle (E1225) _____ N/C
	20-85° Recline Angle (E1226) _____ N/C
	Gas Spring
	Standard Spring _____ N/C
	Heavy Duty Spring _____ N/C
Backrest Height	
<i>Required when Adjustable Height Back with Adjustable Handle Back Post is selected</i>	
18" ¹ 19" ¹ 20" 21" 22" 23" 24" 25" 26" ² _____ N/C	
27" ² _____ N/C	
<i>1 Available only with Short</i>	
<i>2 Available only with Tall</i>	

Footrest

Hanger / Front Frame Bend						
Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount			1			N/C
Front Mount			1	2		N/C
Pro Elevating Leg Rests E0990 ⁴						\$
Contracture Footrest Bilateral Mount ³						\$
Contracture Footrest Center Mount ³						\$
Omit Hangers						(\$)

- 1 7" and 8" Casters may interfere with less than 5° tilt
- 2 Not available with 7" or 8" Casters
- 3 Aluminum Angle Adjustable footplates included. See image to right.
- 4 Available only with classic hanger release

Pro ELR Cover


Available only on Pro ELR

Pro ELR Gel Pad Cover (PAIR) _____ \$

3" Extension

Available only with Front Mount Hangers

3" Extension _____ N/C



Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length										
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"		
Extension Tube Hangers (60°,70°,80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"		
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"		
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"		
Front Mount Hangers ¹	Short Front Mount Hanger (70°,80°,90°)	Standard Selections								Requires 3" Extension for \$35		
		4" ²	4.5" ²	5"	5.5"	6"				8.5"	9"	9.5"
	Medium Front Mount Hanger (70°,80°,90°)	6.5"	7"	7.5"	8"					10"	10.5"	11"
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"		10.5"	11"	11.5"
	Long Front Mount Hanger (70°,80°)	7.5"	8"	8.5"	9"	9.5"	10"			12"	12.5"	13"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	12.5"	13"	13.5"
		9"	9.5"	10"	10.5"	11"	11.5"	12"	14"	14.5"	15"	
Contracture Footrest	Short (10" - 13.5")	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"		
	Medium (11.5" - 15.5")	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"		
	Long (14" - 18")	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
Pro ELR Low Pivot Setting	Pro ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"										
	Pro ELR Low Pivot Setting Long (13" - 19")	13" - 19"										
Pro ELR High Pivot Setting	Pro ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"										
	Pro ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"										

- 1 See Minimum Front Seat Height chart below
- 2 Available with 90° only

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend						
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Footplates

No selection needed with Contracture Footrest

Composite	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable (K0040)	\$
<i>Not available with Front Mount Hanger or 14" Seat Width without Offset Hanger</i>	
Aluminum Angle Adjustable (K0040)	\$
Aluminum Locking Multi-Angle Adjustable (K0108)	\$
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable (K0040)	\$
<i>Not available with Front Mount Hanger or PRO ELR</i>	

Hanger Release

No selection needed with Contracture Footrest. If no selection is made 4-Way - In-Line Hanger Release will be used.

Classic	N/C
4-Way - In-Line	N/C
<i>Not available with Pro ELR</i>	
4-Way - Out	N/C
<i>Not available with Pro ELR</i>	

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$
Adjustable Heel Loops (E0951)	\$
Omit Heel Loops	N/C

Calf Straps

Velcro Adjustable (K0038)	\$
Padded Velcro Adjustable (K0038)	\$
Padded Airmesh Calf Strap (K0038)	\$

Residual Limb Support

Not available with Offset Hanger

	Limb Support	Price	Omit	Cover		Pad Size					
				Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Left	Left Limb Support (E1020)	\$	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
	Seat Height					
12" with Pneumatic or Pneumatic with Airless Insert Tire ³	13.5"					
	14"					
	14.5"	14.5"	14.5"		17"	
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
12" with Full Poly Tire ³	16.5"	16.5"	16.5"			
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"		
12" with Drum Brake ³	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20" ²		20" ²	20" ²	20" ²	20" ²
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire ⁴	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5" ¹	16.5" ¹	16.5" ¹	16.5" ¹		16.5" ¹
	19"	19"	19"	19" ²	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
16" with Low Poly Tire ⁴	13"					
	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
			20.5"	20.5"	20.5"	

Rear Wheel Size	Caster Size					
	5"	6"	6" x 2"	7"	8"	8" x 2"
	Seat Height					
16" with Drum Brake ⁴	13.5"					
	14"	14"	14"			
	14.5"	14.5"	14.5"			
	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"			16"
		18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20"
20"	15"	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"		
	16"	16"	16"	16"	16"	16"
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5" ²	18.5" ²	18.5" ²	18.5" ²		
22"	16"	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
24"	16.5" ⁵	16.5" ^{1.5}	16.5" ^{1.5}			
	17" ⁵	17" ⁵	17" ⁵	17" ⁵	17" ⁵	17" ⁵
	17.5" ⁶	17.5" ⁶	17.5" ⁶	17.5"	17.5"	17.5" ⁶
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"	19.5"
	20"	20"	20"	20"	20"	20" ²

¹ Available only with Full Poly tire

² Not available with Full Poly tire

³ Not available with Vent and Battery Tray combination

⁴ Seat heights higher than 16.5" with 16" wheel size not available with Vent and Battery Tray combination

⁵ Not available with Drum Brake

⁶ Not available with Drum Brake in combination with Pneumatic tires

Casters and Rear Wheels

One Arm Drive

One Arm Drive - Quick Release (E0958) _____ \$
*Available only with 24" Spoke wheels. Spoke wheel style is exclusive to OAD.
 Not available with Vent Tray, Battery Box Holder, Self Propel Axle Plate or
 Angle Adjustable Locking Extendable Flip Up. Standard Handrim will be
 Aluminum Anodized. Axle is Quick Release.*

Left Side Drive _____ N/C
 Right Side Drive _____ N/C
 Aluminum Drive Handrim _____ N/C
 Poly Coated Drive Handrim _____ \$

Casters

		Caster Size				
		Price	5"	6"	7"	8"
Caster Type	1" Poly	N/C				
	1" Poly Aluminum	\$118				
	1" Pneumatic	\$118				
	1.5" Poly	\$95				
	1.5" Soft Roll Aluminum	\$190				
	2" Poly	\$97				
	2" Pneumatic w/Foam Insert	\$246				

Axle

No selection needed with One Arm Drive or Drum Brake

Fixed _____ N/C
 Quick Release _____ N/C

Self Propel Axle Plate

Self Propel Axle Plate _____ \$
Not available with Vent Tray or 12" or 16" wheel sizes

Wheel Locks

Wheel locks

Push to Lock _____ N/C
 Pull to Lock _____ N/C
Not available with 12" or 16" Full Poly tires
 Attendant Foot Lock _____ \$
 Combination Attendant Foot Lock _____ \$
*Available only with 12" or 16" wheel sizes. Includes Extension Handles. Push
 to Lock only.*
 Drum Brake (E2228) _____ \$
*Available only with 12" Mag, 16" Mag, or Spoke wheels. Not available with
 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate.
 Axle is Quick Release.*
 Dual - Drum Brake with Push to Lock (E2228) _____ \$
*Available only with 12" Mag, 16" Mag, or Spoke wheels. Not available with
 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate.
 Axle is Quick Release.*
 Dual - Drum Brake with Pull to Lock (E2228) _____ \$
*Available only with 12" Mag, 16" Mag, or Spoke wheels. Not available with
 16" Low Poly tire, Reclining Backrest, One Arm Drive or Self Propel Axle Plate.
 Axle is Quick Release.*
 Omit Wheel Locks _____ N/C

Rear Wheel

No selection needed with One Arm Drive

		Rear Wheel Size				
		Price	12"	16"	20"	22"
Rear Wheel Type	Mag	N/C				
	Spoke	\$205				

Rear Tire Type

Full Poly (E2218) _____ N/C
 Pneumatic (E2211,E2212) _____ N/C
 Pneumatic with Airless Insert (E2211, E2213) _____ \$
 Low Poly (E2218) _____ N/C
Available only with 16" wheel size

Handrims

Handrims not available on 12" or 16" wheel sizes

Standard Handrims

Aluminum Anodized _____ N/C
 Plastic Coated _____ \$
 Projection (8 Vertical Projections) _____ \$
*Available only with full profile tires in combination with 22", 24", and 26"
 MAXX Spoke wheels, or 22" and 24" Mag wheels. No selection needed with
 One Arm Drive.*
 Omit Handrim _____ N/C

Handrim Options

Not available on Plastic Coated or Projection Handrims

3M Micro Grip _____ \$

Wheel Lock Options

Extension Handles _____ \$
Available only with Push to Lock or Pull to Lock Wheel Locks

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$
	Full	

¹ If requested configuration is unavailable, Standard Height T-Arm will be used. Not available with 14" seat depth and Offset Hangers"

² Not available with IV Pole

³ Not available with Dynamic Back

⁴ Minimum armrest height with Reclining Backrest is 12"

Armrest Type	Length	Price
Height Adj. Dual Post Flip Back (9.5"-14.5") E0973 ²	Desk	\$
	Full	
Angle Adjustable Locking Extendable Flip Up K0020 ^{3,4}	Desk	\$
	Full	
Omit Armrests		N/C

Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guard

Composite Side Guards (K0108)	\$
<i>Not available with T-Arm or Dual Post Flip Back Armrests</i>	

Anti-Tip

Rear Anti-Tips (E0971)	\$
<i>Required with Reclining Backrest or Vent Tray</i>	
Omit	N/C

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$
Small Flex-Mount Positioning Belt (16-24") (E0978)	\$
Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$
Large Flex-Mount Positioning Belt (25-37") (E0978)	\$

Belt Mounting Kit

Ki Mobility Belt Mounting Brackets	\$
Belt Mounting Band Clamp	\$

U.S. Military Service Patches

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$

O2 Holder

O2 Holder (E2208)	\$
<i>Not available with Reclining Backrest or Dynamic Back</i>	

Other Accessories

Tool Kit	\$
Removable Underseat Pouch	\$
Ki Mobility Backpack	\$
Black Neoprene Impact Guards	\$
Tilt Stop (Pair)	\$
<i>Included with Offset Hanger</i>	
IV Pole	\$
Canopy with Hardware (K0108)	\$
<i>Not available on frame widths greater than 18"</i>	

Frame Color

Gloss Finish

Black	N/C
Candy Blue	N/C
Candy Purple	N/C
Candy Red	N/C
Granny Smith	N/C
Metallic Orange	N/C
Yellow	N/C
Opalescent Violet	N/C
Radiant Orange	N/C
Electric Green	N/C
Glacier Blue	N/C
Dark Pink	N/C
Light Pink	N/C
Purple Plum	N/C
Rich Teal	N/C

Matte Finish

Charcoal Gray	N/C
Iridescent Burgundy	N/C
Iridescent Navy	N/C
Royal Blue	N/C
Shadow Black	N/C
Bright Red	N/C
White	N/C
Dry Earth Tan	N/C
Tactical Green	N/C

Metallic Finish

Blue Steel	N/C
------------	-----

FOCUS CR

Cushions and Backrests

AXIOM

Ki Mobility

ADI
A Stealth Product

VARILITE

RIE

Cushion:

Backrest:

Other Accessories: