

Ki Mobility Liberty FT Express Script



Customer Details

Client				
Assessor				
Organisation				
User Weight	ACC		Whaikaha	Other

Ki Mobility Liberty FT Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

Choose either A or B		Seat Width			
		14"	16"	18"	20"
A	14"	<input type="checkbox"/>	<input type="checkbox"/> * (.ULLIB161415D)	<input type="checkbox"/> * (.ULLIB181415D)	<input type="checkbox"/> < 56kg <input type="checkbox"/> 56 - 68kg <input type="checkbox"/> 68 - 79kg
	15"	<input type="checkbox"/>	<input type="checkbox"/> * (.ULLIB161415D)	<input type="checkbox"/> * (.ULLIB181415D)	
	16"	<input type="checkbox"/> * (.ULLIB1416D)	<input type="checkbox"/> * (.ULLIB1616D)	<input type="checkbox"/> * (.ULLIB1816D)	
B	18"	<input type="checkbox"/> 68 - 79kg	<input type="checkbox"/>	<input type="checkbox"/> * (.ULLIB181819D)	<input type="checkbox"/> * (.ULLIB201819D)
	19"	<input type="checkbox"/> 79 - 90kg	<input type="checkbox"/>	<input type="checkbox"/> * (.ULLIB181819D)	<input type="checkbox"/> * (.ULLIB201819D)
	20"	<input type="checkbox"/> > 90kg	<input type="checkbox"/> * (.ULLIB1620D)	<input type="checkbox"/> * (.ULLIB1820D)	<input type="checkbox"/> * (.ULLIB2020D)

Code for all Express Packages is

KILIB1414STDNAVY / RED

KILIB1618STDNAVY / RED

* Please enter code in brackets, when shown, as well as the Standard Express Package code, above.

A User Weight Range < 80kg

B User Weight Range 68-90+kg

A 16W needs Back Upholstery (KI004523)

AB 18W needs Back Upholstery (KI004525)

B 20W needs Back Upholstery (KI004527)

Express Package: (Select optional specifications and alternatives on reverse)

- Standard seat sling
- Transit tie downs
- Height adjustable straight back cane with push handle
- Fold down back (angle adjustable)
- Standard back upholstery
- 80° legrest hanger / front frame bend
- Aluminium angle adjustable footplates
- Long extension tube hangers
- 4-Way in line hanger release
- 6" castors - 1.5" soft roll aluminium
- 24in Maxx Spoke Performance Pneumatic Black Wheels
- Push to lock wheel locks
- Aluminium anodised handrims
- Quick release axles
- Height adjustable T-arm armrests with full arm pads
- Rear anti-tips
- Iridescent Navy Frame

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height

<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"			

Seat to Footrest Length

<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"
<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	

Front Seat Height

A	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"
B	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"

Other specifications are available on reverse as a non-standard express package

Optional Specifications

Back Upholstery 1

<input type="checkbox"/>	Omit back upholstery
--------------------------	----------------------

Leg Hanger / Font Frame Bend 2

<input type="checkbox"/>	70° (.ULLIBLEG70NAVY)	<input type="checkbox"/>	Omit Hanger
--------------------------	-----------------------	--------------------------	-------------

Front Seat Height with 8 x 1" Poly + 7" Forks 3

.ULLIBFSH21 (for 14" to 18" deep)

.ULLIBFSH1920 (for 19" to 20" deep)

<input type="checkbox"/>	18.5"	<input type="checkbox"/>	19"	<input type="checkbox"/>	19.5"	<input type="checkbox"/>	20"
--------------------------	-------	--------------------------	-----	--------------------------	-------	--------------------------	-----

Armpad 4

<input type="checkbox"/>	Desk arm pad (.ULLIBARPD)
--------------------------	---------------------------

Adjustable Heel Loops 5

<input type="checkbox"/>	14-15W (.ULLIBHL1415)
<input type="checkbox"/>	16-22W (.ULLIBHL1622)

Calf Strap 6

Velcro adjustable

<input type="checkbox"/>	15-17W (.ULLIBCSVA1517)
<input type="checkbox"/>	18-22W (.ULLIBCSVA1822)

Padded velcro adjustable

<input type="checkbox"/>	15-17W (.ULLIBCSPPVA1517)
<input type="checkbox"/>	18-22W (.ULLIBCSPPVA1822)

Bodypoint padded aeromesh

<input type="checkbox"/>	14-15W (.ULLIBHL1415)
<input type="checkbox"/>	15-17W (.ULLIBCSAM1517)
<input type="checkbox"/>	16-22W (.ULLIBHL1622)
<input type="checkbox"/>	18-22W (.ULLIBCSAM1822)

Notes

Wheel Lock Options (Only on Std. Push to lock) 7

<input type="checkbox"/>	Extension handles (.ULLIBWLKEXT)
--------------------------	----------------------------------

Push Rim Options 8

<input type="checkbox"/>	Plastic coated push rim (KI112821)
<input type="checkbox"/>	Natural fit omit thumb grip (KI100793)
<input type="checkbox"/>	Thumb grip regular (KI100795)

Rear Wheel Transit Kit 9

<input type="checkbox"/>	Rear wheel transit kit (.ULLIBWHL02) Includes the following: • 16"Mag wheel • Street Pneumatic
--------------------------	---

Rear Wheel Kit - Low profile tyres 10

<input type="checkbox"/>	24" Maxx Spoke Low Profile Schwalbe Marathon Plus natural fit LT omit thumb grip (KI119391)
<input type="checkbox"/>	Plastic coated push rim (KI112821)
<input type="checkbox"/>	Std aluminum push rim (KI200349)
<input type="checkbox"/>	Thumb grip regular (KI100795)
<input type="checkbox"/>	Natural fit standard (KI100793)

PRO Elevated Legrest - Angle Adj. FP Long 11

<input type="checkbox"/>	16/17 LHS (504526)	<input type="checkbox"/>	16/17 RHS (504534)
<input type="checkbox"/>	18-25 LHS (504527)	<input type="checkbox"/>	18-25 RHS (504543)

Residual Limb Support 12

5" x 8" with smooth cover

<input type="checkbox"/>	Left (KI502355)	<input type="checkbox"/>	Right (KI502354)
--------------------------	-----------------	--------------------------	------------------

7" x 10" with smooth cover

<input type="checkbox"/>	Left (KI502357)	<input type="checkbox"/>	Right (KI502356)
--------------------------	-----------------	--------------------------	------------------

Accessories 13

<input type="checkbox"/>	Composite side guards (not used with armrests) (.ULLIBSDGC)
<input type="checkbox"/>	Black neoprene impact guards (.ULLIBIMG)
<input type="checkbox"/>	Clear spoke guards (.ULLIBSPGCL24)
<input type="checkbox"/>	Tray

Notes

A large, empty rectangular area with a light gray background, intended for taking notes. It is bounded by a thin black line on the left, right, and bottom, and a dark blue header bar at the top.