

Account Number:

Date:

Address:

Quote: PO Number:

City:

Associated Quote/Order#:

Phone:

Email/Fax:

Effective: 6/17/2024

|   |    |
|---|----|
| Liberty - Transit Included (E1161)            | \$ |
| Liberty - Transit Not Included (E1161)        | \$ |
| Liberty - Rigid, Transit Included (E1161)     | \$ |
| Liberty - Rigid, Transit Not Included (E1161) | \$ |

250 lb weight capacity Tilt range: 0° - 20°

## Frame

### Seat Width

Select Seat Width within desired Seat Width Range

| Frame Width Range  | Frame Width |         |                   |                   |                   |
|--------------------|-------------|---------|-------------------|-------------------|-------------------|
| Small (14" - 18")  | 14" N/C     | 15" N/C | 16" N/C           | 17" N/C           | 18" N/C           |
| Medium (16" - 20") | 16" N/C     | 17" N/C | 18" N/C           | 19" N/C           | 20" (E2201)<br>\$ |
| Large (18" - 22")  | 18" N/C     | 19" N/C | 20" (E2201)<br>\$ | 21" (E2201)<br>\$ | 22" (E2201)<br>\$ |

### Seat Depth

Select Seat Depth within desired Seat Depth Range

| Seat Depth Range   | Seat Depth |         |                   |
|--------------------|------------|---------|-------------------|
| Short (14" - 16")  | 14" N/C    | 15" N/C | 16" N/C           |
| Medium (16" - 18") | 16" N/C    | 17" N/C | 18" N/C           |
| Long (18" - 20")   | 18" N/C    | 19" N/C | 20" (E2203)<br>\$ |

### Assistive Gas Strut Set-Up

|  |     |
|--|-----|
| Client Weight > 200 lbs (Setting 1, 450N Spring)   | N/C |
| Client Weight 176-200 lbs (Setting 2, 450N Spring) | N/C |
| Client Weight 151-175 lbs (Setting 3, 450N Spring) | N/C |
| <i>Standard on chairs less than 16" deep</i>       |     |
| Client Weight 151-175 lbs (Setting 1, 300N Spring) | N/C |
| <i>Standard on chairs 16" deep or larger</i>       |     |
| Client Weight 125-150 lbs (Setting 2, 300N Spring) | N/C |
| Client Weight < 125 lbs (Setting 3, 300N Spring)   | N/C |

### Seat Option

|                                     |     |
|-------------------------------------|-----|
| Standard Seat Sling                 | N/C |
| Growing Adjustable Seat Pan (E2231) | \$  |
| Omit Seat                           | N/C |

### Transit Option

|  |     |
|--|-----|
| Transit System   | \$  |
| <i>No charge when selected with Liberty - Transit Included</i> |     |
| Omit   | N/C |

## Backrest

### Back Post / Height Range / Back Height

| Back Post   | Price | Height Range | Back Height |     |     |     |     |
|---|-------|--------------|-------------|-----|-----|-----|-----|
| Stroller Handle Back Post <sup>1</sup>                  | N/C   | Short        | 20"         |     |     |     |     |
|   |       | Tall         | 24"         |     |     |     |     |
| Height Adjustable Straight with Push Handle             | \$26  | Short        | 14"         | 15" | 16" | 17" | 18" |
|   |       | Medium       | 16"         | 17" | 18" | 19" | 20" |
|   |       | Tall         | 18"         | 19" | 20" | 21" | 22" |
| Height Adjustable 8° Bend with Push Handle <sup>2</sup> | \$26  | Short        | 14"         | 15" | 16" | 17" | 18" |
|   |       | Medium       | 16"         | 17" | 18" | 19" | 20" |
|   |       | Tall         | 18"         | 19" | 20" | 21" | 22" |

<sup>1</sup> Not available with Back Upholstery

<sup>2</sup> Not available with IV Pole

### Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension \_\_\_\_\_ \$

### Backrest Options

Fixed (Non-Folding) \_\_\_\_\_ N/C

Angle Adjustable (Folding) \_\_\_\_\_ \$

### Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard \_\_\_\_\_ N/C

Tension Adjustable (E2611) \_\_\_\_\_ \$

*Not available on Seat Widths less than 16" combined with 21" and 22" Backrest Heights*

Omit \_\_\_\_\_ (\$ )

### U.S. Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

Army \_\_\_\_\_ Navy \_\_\_\_\_ Marines \_\_\_\_\_ Air Force \_\_\_\_\_  
 Coast Guard \_\_\_\_\_ US Purple Heart \_\_\_\_\_ \$

Combo Army/Purple Heart \_\_\_\_\_ Combo Navy/Purple Heart \_\_\_\_\_  
 Combo Marines/Purple Heart \_\_\_\_\_ Combo Air Force/Purple Heart \_\_\_\_\_  
 Combo Coast Guard/Purple Heart \_\_\_\_\_ \$

## Footrest

### Hanger

| Hanger                        | 60° | 70° | 80° | 90° |  | Retail |
|-------------------------------|-----|-----|-----|-----|--|--------|
| Extension Tube Mount          |     |     | 1   |     |  | N/C    |
| Front Mount                   |     |     | 1   | 2   |  | N/C    |
| Pro Elevating Leg Rests E0990 |     |     |     |     |  | \$     |
| Omit Hangers                  |     |     |     |     |  | (\$ )  |

<sup>1</sup> Not available with 8" casters

<sup>2</sup> Not available with 8" casters in combination with 3" extension

### PRO ELR Cover

Available only on PRO ELR

Pro ELR Gel Pad Cover (PAIR) \_\_\_\_\_ \$

### Footplates

Composite \_\_\_\_\_ N/C

*Not available with Front Mount Hanger*

Composite Angle Adjustable (K0040) \_\_\_\_\_ \$

*Not available with Front Mount Hanger or 14" Seat Width*

Aluminum Angle Adjustable (K0040) \_\_\_\_\_ \$

Aluminum Locking Multi-Angle Adjustable (K0108) \_\_\_\_\_ \$

*Available only with Front Mount Hanger*

One Piece Flip Up Angle Adjustable (K0040) \_\_\_\_\_ \$

*Not available with Front Mount Hanger or PRO ELR*

## Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

| Footrest                                  | Extension  | Seat to Footrest Length |                   |       |       |       |       |       |       |                       |       |       |
|---|--|-------------------------|-------------------|-------|-------|-------|-------|-------|-------|-----------------------|-------|-------|
|   |  | 10"                     | 10.5"             | 11"   | 11.5" | 12"   | 12.5" | 13"   | 13.5" | 15.5"                 | 18"   |       |
| Extension Tube Hangers<br>(60°, 70°, 80°) | Short (10" - 13.5")                                | 10"                     | 10.5"             | 11"   | 11.5" | 12"   | 12.5" | 13"   | 13.5" |                       |       |       |
|   | Medium (11.5" - 15.5")                             | 11.5"                   | 12"               | 12.5" | 13"   | 13.5" | 14"   | 14.5" | 15"   | 15.5"                 |       |       |
|   | Long (14" - 18")                                   | 14"                     | 14.5"             | 15"   | 15.5" | 16"   | 16.5" | 17"   | 17.5" | 18"                   |       |       |
| Front Mount Hangers <sup>1</sup>          | Short Front Mount Hanger<br>(70°, 80°, 90°)        | Standard Selections     |                   |       |       |       |       |       |       | Requires 3" Extension |       |       |
|   |  | 4" <sup>2</sup>         | 4.5" <sup>2</sup> | 5"    | 5.5"  | 6"    |       |       |       | 8.5"                  | 9"    | 9.5"  |
|   | Medium Front Mount Hanger<br>(70°, 80°, 90°)       | 6.5"                    | 7"                | 7.5"  | 8"    |       |       |       |       | 10"                   | 10.5" | 11"   |
|   |  | 4" <sup>2</sup>         | 4.5" <sup>2</sup> | 5"    | 5.5"  | 6"    | 6.5"  | 7"    |       | 10.5"                 | 11"   | 11.5" |
|   | Long Front Mount Hanger<br>(70°, 80°)              | 7.5"                    | 8"                | 8.5"  | 9"    | 9.5"  | 10"   |       |       | 12"                   | 12.5" | 13"   |
|   |  | 5"                      | 5.5"              | 6"    | 6.5"  | 7"    | 7.5"  | 8"    | 8.5"  | 12.5"                 | 13"   | 13.5" |
|   |  | 9"                      | 9.5"              | 10"   | 10.5" | 11"   | 11.5" | 12"   | 14"   | 14.5"                 | 15"   |       |
| PRO ELR Low Pivot<br>Setting              | PRO ELR Low Pivot Setting Short<br>(10.5" - 15.5") | 10.5" - 15.5"           |                   |       |       |       |       |       |       |                       |       |       |
|   | PRO ELR Low Pivot Setting Long<br>(13" - 19")      | 13" - 19"               |                   |       |       |       |       |       |       |                       |       |       |
| PRO ELR High Pivot<br>Setting             | PRO ELR High Pivot Setting Short<br>(8.5" - 17")   | 8.5" - 17"              |                   |       |       |       |       |       |       |                       |       |       |
|   | PRO ELR High Pivot Setting Long<br>(11.5" - 17.5") | 11.5" - 17.5"           |                   |       |       |       |       |       |       |                       |       |       |

<sup>1</sup> See Minimum Seat Height chart below

<sup>2</sup> Available with 90° only

| Minimum Seat Height with Front Mount Hangers |       |        |        |        |       |       |
|--|-------|--------|--------|--------|-------|-------|
| Hanger Length                                | Short |        | Medium |        | Long  |       |
|  | 70°   | 80/90° | 70°    | 80/90° | 70°   | 80°   |
| Minimum Seat Height with Standard Selection  | 13"   | 13"    | 13"    | 13.5"  | 14.5" | 15.5" |
| Minimum Seat Height with 3" Extension        | 13.5" | 14"    | 15.5"  | 16"    | 17"   | 18"   |

### Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951) \_\_\_\_\_ \$  
Omit Heel Loops \_\_\_\_\_ N/C

### Calf Straps

Velcro Adjustable (K0038) \_\_\_\_\_ \$  
Padded Velcro Adjustable (K0038) \_\_\_\_\_ \$  
Padded Airmesh Calf Strap (K0038) \_\_\_\_\_ \$

### Hanger Extension Tubes

Available only with Front Mount Hangers

3" Extension \_\_\_\_\_ N/C  
3" Extension (Provided as Parts) \_\_\_\_\_ N/C

### Residual Limb Support

|       | Limb Support               | Price | Omit              | Cover              |                      | Pad Size           |                     |                     |                     |                      |                      |
|-------|----------------------------|-------|-------------------|--------------------|----------------------|--------------------|---------------------|---------------------|---------------------|----------------------|----------------------|
|       |                            |       |                   | Smooth Cover Left  | Air Mesh Cover Left  | 5X8 Pad Left Side  | 7X10 Pad Left Side  | 9x10 Pad Left Side  | 9x14 Pad Left Side  | 11x10 Pad Left Side  | 11x14 Pad Left Side  |
| Left  | Left Limb Support (E1020)  | \$431 | Omit Left Hanger  | Smooth Cover Left  | Air Mesh Cover Left  | 5X8 Pad Left Side  | 7X10 Pad Left Side  | 9x10 Pad Left Side  | 9x14 Pad Left Side  | 11x10 Pad Left Side  | 11x14 Pad Left Side  |
| Right | Right Limb Support (E1020) | \$431 | Omit Right Hanger | Smooth Cover Right | Air Mesh Cover Right | 5X8 Pad Right Side | 7X10 Pad Right Side | 9x10 Pad Right Side | 9x14 Pad Right Side | 11x10 Pad Right Side | 11x14 Pad Right Side |

### Hanger Release

4-Way - In-Line \_\_\_\_\_ N/C  
Not available with PRO ELR  
4-Way - Out \_\_\_\_\_ N/C

## Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

| Rear Wheel Size | Caster Size        |                  |                  |                  |
|-----------------|--------------------|------------------|------------------|------------------|
|                 | 4"                 | 5"               | 6"               | 8"               |
|                 | Seat Height        |                  |                  |                  |
| 12"             | 13"                | 14"              | 15"              |                  |
|                 | 13.5"              | 14.5"            | 15.5"            | 17"              |
|                 | 14"                | 15"              | 16"              | 17.5"            |
|                 | 14.5"              | 15.5"            | 16.5"            | 18"              |
|                 | 15"                | 16"              | 17"              | 18.5"            |
|                 | 15.5"              | 16.5"            | 17.5"            | 19" <sup>2</sup> |
|                 | 16"                | 17"              | 18"              |                  |
|                 |                    | 17.5"            | 18.5"            |                  |
|                 |                    | 18"              | 19" <sup>2</sup> |                  |
|                 |                    | 18.5"            |                  |                  |
|                 | 19" <sup>2</sup>   |                  |                  |                  |
| 16"             | 13.5" <sup>1</sup> | 14" <sup>1</sup> | 15"              |                  |
|                 | 14" <sup>1</sup>   | 14.5"            | 15.5"            | 17"              |
|                 | 14.5"              | 15"              | 16"              | 17.5"            |
|                 | 15"                | 15.5"            | 16.5"            | 18"              |
|                 | 15.5"              | 16"              | 17"              | 18.5"            |
|                 | 16"                | 16.5"            | 17.5"            | 19"              |
|                 |                    | 17"              | 18"              | 19.5"            |
|                 |                    | 17.5"            | 18.5"            | 20"              |
|                 |                    | 18"              | 19"              |                  |
|                 |                    | 18.5"            |                  |                  |
|                 | 19"                |                  |                  |                  |
| 20"             | 13"                | 14"              | 15"              |                  |
|                 | 13.5"              | 14.5"            | 15.5"            | 17"              |
|                 | 14"                | 15"              | 16"              | 17.5"            |
|                 | 14.5"              | 15.5"            | 16.5"            | 18"              |
|                 | 15"                | 16"              | 17"              | 18.5"            |
|                 | 15.5"              | 16.5"            | 17.5"            | 19"              |
|                 | 16"                | 17"              | 18"              | 19.5"            |
|                 |                    | 17.5"            | 18.5"            | 20"              |
|                 |                    | 18"              | 19"              |                  |
|                 |                    | 18.5"            |                  |                  |
|                 | 19"                |                  |                  |                  |

<sup>1</sup> Available only with Low Profile tires

<sup>2</sup> Not available with Full Poly tire

| Rear Wheel Size | Caster Size        |                    |                  |       |
|-----------------|--------------------|--------------------|------------------|-------|
|                 | 4"                 | 5"                 | 6"               | 8"    |
|                 | Seat Height        |                    |                  |       |
| 22"             | 13.5" <sup>1</sup> | 14"                | 15"              |       |
|                 | 14"                | 14.5"              | 15.5"            | 17"   |
|                 | 14.5"              | 15"                | 16"              | 17.5" |
|                 | 15"                | 15.5"              | 16.5"            | 18"   |
|                 | 15.5"              | 16"                | 17"              | 18.5" |
|                 | 16"                | 16.5"              | 17.5"            | 19"   |
|                 |                    | 17"                | 18"              | 19.5" |
|                 |                    | 17.5"              | 18.5"            | 20"   |
|                 |                    | 18"                | 19" <sup>2</sup> |       |
|                 |                    | 18.5"              |                  |       |
|                 | 19" <sup>2</sup>   |                    |                  |       |
| 24"             | 14.5" <sup>1</sup> | 14.5" <sup>1</sup> | 15"              |       |
|                 | 15"                | 15"                | 15.5"            | 17"   |
|                 | 15.5"              | 15.5"              | 16"              | 17.5" |
|                 | 16"                | 16"                | 16.5"            | 18"   |
|                 |                    | 16.5"              | 17"              | 18.5" |
|                 |                    | 17"                | 17.5"            | 19"   |
|                 |                    | 17.5"              | 18"              | 19.5" |
|                 |                    | 18"                | 18.5"            | 20"   |
|                 |                    | 18.5"              | 19"              |       |
|                 |                    | 19"                |                  |       |

## Casters

|             |                            | Caster Size |    |    |    |    |
|-------------|----------------------------|-------------|----|----|----|----|
|             |                            | Price       | 4" | 5" | 6" | 8" |
| Caster Type | 1" Poly                    | N/C         |    |    |    |    |
|             | 1" Poly Aluminum           | \$          |    |    |    |    |
|             | 1" Pneumatic               | \$          |    |    |    |    |
|             | 1.5" Poly                  | \$          |    |    |    |    |
|             | 1.5" Soft Roll Aluminum    | \$          |    |    |    |    |
|             | 2" Poly                    | \$          |    |    |    |    |
|             | 2" Pneumatic w/Foam Insert | \$          |    |    |    |    |

## Rear Wheel

|                 |       | Rear Wheel Size |     |     |     |     |
|-----------------|-------|-----------------|-----|-----|-----|-----|
|                 |       | Price           | 12" | 16" | 20" | 22" |
| Rear Wheel Type | Mag   | N/C             |     |     |     |     |
|                 | Spoke | N/C             |     |     |     |     |

## Rear Tire Type

### Full Profile

|  |     |
|--|-----|
| Street Pneumatic (E2211, E2212)              | N/C |
| Full Poly 1 3/8" (E2218)                     | N/C |
| Pneumatic with Airless Insert (E2211, E2213) | \$  |

### Low Profile

|  |     |
|--|-----|
| Low Poly 1" (E2218)                                | N/C |
| <i>Available only with 16" wheel size</i>          |     |
| Sentinel Solid Tire 1" (E2218)                     | \$  |
| Schwalbe Marathon Plus Evolution 1" (E2211, E2212) | \$  |
| <i>Available only on 24" Spoke</i>                 |     |

## Wheel Locks

### Wheel locks

|  |     |
|--|-----|
| Push to Lock   | N/C |
| Pull to Lock   | N/C |
| <i>Not available with 12" or 16" Full Poly tires</i> |     |
| Hemi Wheel Lock - Left (K0108)                       | \$  |
| <i>Extension Handle included</i>                     |     |
| Hemi Wheel Lock - Right (K0108)                      | \$  |
| <i>Extension Handle included</i>                     |     |

## Handrims

### Standard Handrims

*Handrims not available on 12" or 16" wheel sizes*

|   |     |
|---|-----|
| Aluminum Anodized   | N/C |
| Plastic Coated  | \$  |
| Projection (8 Vertical Projections)   | \$  |
| <i>Available only with full profile tires in combination with 22" and 24" Spoke wheels, or 22" and 24" Mag wheels</i> |     |
| Omit Handrim  | N/C |

### Standard Handrim Mount

*If no selection is made, Standard will be used*

|                                      |     |
|--------------------------------------|-----|
| Close Mount                          | N/C |
| <i>Reduces overall width by .75"</i> |     |

### Handrim Options

*Not available on Plastic Coated or Projection Handrims*

|               |    |
|---------------|----|
| 3M Micro Grip | \$ |
|---------------|----|

### Axles

*No selection needed with Drum Brake*

|               |     |
|---------------|-----|
| Fixed         | N/C |
| Quick Release | N/C |

### Wheel Lock Options

|  |    |
|--|----|
| Extension Handles (E0961)  | \$ |
| <i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>                                      |    |
| Grade Aids (E0974)   | \$ |
| <i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i> |    |

## Armrests

| Armrest Type   | Length | Price |
|--|--------|-------|
| Height Adjustable T-Arm Tall (11.5"-15") E0973                 | Desk   | \$    |
|  | Full   |       |
| Height Adjustable T-Arm (8.5"-13.5") E0973                     | Desk   | \$    |
|  | Full   |       |
| Angle Adjustable Locking Extendable Flip Up K0020 <sup>1</sup> | Desk   | \$    |
|  | Full   |       |
| Omit Armrests  |        | N/C   |

<sup>1</sup> Available only with 12" and 16" wheel sizes and Stroller Handle Back Post

## Accessories

### Removable Side Guard

Composite Side Guards (K0108) \$  
*Not available with T-Arms*

### Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978) \$  
 1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978) \$  
 2" Airline Buckle Pelvic Positioning Belt (E0978) \$  
 Small Flex-Mount Positioning Belt (16-24") (E0978) \$  
 Medium Flex-Mount Positioning Belt (19-29") (E0978) \$  
 Large Flex-Mount Positioning Belt (25-37") (E0978) \$

### Belt Mounting Kit

Seat Belt Mounting Hardware \$  
 Belt Mounting Band Clamp \$

### O2 Holder

O2 Holder (E2208) \$  
*Mounted to the left side. Not available with Fold Down Back.*

## Frame Color

|                     |     |
|---------------------|-----|
| Black               | N/C |
| Candy Blue          | N/C |
| Charcoal Gray       | N/C |
| Candy Purple        | N/C |
| Iridescent Burgundy | N/C |
| Candy Red           | N/C |
| Iridescent Navy     | N/C |
| Bright Red          | N/C |
| Rich Teal           | N/C |
| Dry Earth Tan       | N/C |
| Tactical Green      | N/C |
| Radiant Orange      | N/C |
| Dark Pink           | N/C |

### Arm Pads

Standard Arm Pad N/C  
 Waterfall Arm Pad \$  
*Not available with Angle Adjustable Locking Extendable Flip Up Armrests*  
 Foam Grip N/C  
*Available only with Angle Adjustable Locking Extendable Flip Up Armrest*

### Anti-Tip

Rear Anti-Tips (E0971) \$  
 Omit N/C

### Transit Wheels

Transit Wheels (K0108) \$  
*Not available with Long Depth Seat Frame, Seat Heights less than 14" or greater than 18", 12" or 16" rear wheels, Cane and Crutch Holder, and O2 Holder*

### Other Accessories

Tool Kit \$  
 Removable Underseat Pouch \$  
 Ki Mobility Backpack \$  
 Black Neoprene Impact Guards \$  
 Cane and Crutch Holder \$  
 IV Pole - Mounted Right Side \$  
*Not available with Transit Option or 8° Bend with Push Handle Back Post*  
 IV Pole - Mounted Left Side \$  
*Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O2 Holder.*

Cushions and Backrests

AXIOM



Cushion:

Backrest: