

Account Number:

Date:

Address:

Quote: PO Number:

City:

Phone:

Email:

Effective: 6/17/2024

Flip - Rigid Base - Transit Included (E1233, E1161)	\$
Flip - Rigid Base - Transit Not Included (E1233, E1161)	\$
Flip - Folding Base - Transit Included (E1234, E1161)	\$
Flip - Folding Base - Transit Not Included (E1234, E1161)	\$

165 lb weight capacity Tilt range: 0° - 45° Maximum Seat Width of 18"

## Frame

### Frame Width

Select box corresponding to desired Frame Width

Frame Width									
10"	11"	12"	13"	14"	15"	16"	17"	18"	

### Frame Depth

Select Frame Depth within desired Frame Depth Range

Frame Depth Range	Frame Depth				
Short (12" - 16")	12" N/C	13" N/C	14" N/C	15" N/C	16" N/C
Medium (14" - 18")	14" N/C	15" N/C	16" N/C	17" N/C	18" N/C
Long (16" - 20")	16" N/C	17" N/C	18" N/C	19" N/C	20" (E2203) \$

### Transit Option

Transit System	\$
<i>No charge when selected with Base Options with Transit Included</i>	
Omit Transit	N/C

### Offset Hanger Option

Offset Hanger	\$
<i>Available only with Rigid Base</i>	

### Tilt Mechanism

Hand Tilt	N/C
Foot Tilt	N/C

### Seat Pan

Growing Solid Seat Pan (E2231)	\$
Omit Seat Pan	N/C

## Backrest

### Back Post / Height Range / Back Height

Stroller Handle Back Post (Short 20")	N/C
Stroller Handle Back Post (Tall 24")	N/C
Fixed Height Back with Adjustable Handle (Short 22")	N/C
<i>Not available with Dynamic Back</i>	
Fixed Height Back with Adjustable Handle (Tall 26")	N/C
<i>Not available with Dynamic Back</i>	
Reclining Back with Adjustable Handle (Short 22") (Non-Folding) (E1225)	\$
<i>Not available with Dynamic Back or Rigidizer Bar. Available only with Attendant Foot Lock. Minimum Seat Height of 15" with 12" Mag wheel and 16" with 16" Mag wheel.</i>	
Reclining Back with Adjustable Handle (Tall 26") (Non-Folding) (E1225)	\$
<i>Not available with Dynamic Back or Rigidizer Bar. Available only with Attendant Foot Lock. Minimum Seat Height of 15" with 12" Mag wheel and 16" with 16" Mag wheel.</i>	
Adjustable Height Back with Adjustable Handle (Short 18" - 25")	N/C
<i>Not available with Dynamic Back</i>	
Adjustable Height Back with Adjustable Handle (Tall 20" - 27")	N/C
<i>20"-22" Heights not available with Dynamic Back</i>	

### Adjustable Height Back Post Height

18" Back Height	19" Back Height	20" Back Height	21" Back Height
22" Back Height	23" Back Height	24" Back Height	25" Back Height
26" Back Height	27" Back Height		N/C

### Stroller Handle Back Post Option

*Available only with Stroller Handle Back Post*

Removable Stroller Handle Extension	\$
-------------------------------------	----

### Other Backrest Options

Non Folding Back	N/C
<i>Standard with Dynamic Back</i>	
Dynamic Back (Non-Folding) (E2398)	
<i>Not available with Flip Up Armrests. Available only with Non Folding Back.</i>	
Adjustable Height Rigidizer Bar	\$

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN129.25 Pricing and specifications subject to change without notice.

## Footrest

Footrest Hanger Type					Pro ELR Cover	
Hanger	70°	80°	90°		Retail	
Extension Tube Mount		1			N/C	Available only on Pro ELR Pro ELR Gel Pad Cover (PAIR) _____ \$
Front Mount		1	1		N/C	
Pro Elevating Leg Rests E0990 <sup>2</sup>					\$	
Pediatric Pro Elevating Leg Rests E0990					\$	
Omit Hangers					(\$ )	

<sup>1</sup> Not available with 8" casters

<sup>2</sup> Available only with Classic hanger release. Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers.

Hanger Extension Tubes	
Available only with Front Mount Hangers	
3" Extension	\$
3" Extension (Provided as Parts)	\$

## Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length												
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"			
Extension Tube Hangers (70°,80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"				
Front Mount Hangers <sup>1</sup>	Short Front Mount Hanger (70°,80°,90°)	Standard Selections										Requires 3" Extension for \$55		
		4" <sup>4</sup>	4.5" <sup>4</sup>	5"	5.5"	6"						8.5"	9"	9.5"
	Medium Front Mount Hanger (70°,80°,90°)	6.5"	7"	7.5"	8"							10"	10.5"	11"
		4" <sup>4</sup>	4.5" <sup>4</sup>	5"	5.5"	6"	6.5"	7"				10.5"	11"	11.5"
	Long Front Mount Hanger (70°,80°)	7.5"	8"	8.5"	9"	9.5"	10"					12"	12.5"	13"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"			12.5"	13"	13.5"
Peds PRO ELR	Low Pivot Setting	6" <sup>3</sup>	6.5" <sup>3</sup>	7" <sup>3</sup>	7.5"	8"	8.5"	9" <sup>2</sup>	9.5" <sup>2</sup>	10" <sup>2</sup>	10.5" <sup>2</sup>	11" <sup>2</sup>	11.5" <sup>2</sup>	
	High Pivot Setting	6" <sup>3</sup>	6.5" <sup>3</sup>	7" <sup>3</sup>	7.5" <sup>2</sup>	8" <sup>2</sup>	8.5" <sup>2</sup>	9" <sup>2</sup>	9.5" <sup>2</sup>	10" <sup>2</sup>				
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"												
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"												
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"												

<sup>1</sup> See Minimum Front Seat Height chart to the right

<sup>2</sup> Comes standard with 3" Extension

<sup>3</sup> Not available on chairs 14" wide or more

<sup>4</sup> Available only with 90° Front Mount Hangers

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Front Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

### Footplates

*No selection needed with Contracture Footrest*

Composite	N/C
<i>Not available with Front Mount Hanger, Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hanger.</i>	
Composite Angle Adjustable (K0040)	\$
<i>Not available with Front Mount Hanger or 10"-12" Frame Widths. 13" and 14" Frame Widths require Offset Hanger</i>	
Aluminum Angle Adjustable (K0040)	\$
Aluminum Locking Multi-Angle Adjustable (K0108)	\$
<i>Available only with Front Mount Hanger</i>	

### Hanger Release

*No selection needed with Contracture Footrest. If no selection is made 4-Way - In line will be used.*

Classic	N/C
4-Way - In-Line	N/C
<i>Allows chair to be folded without removing footrest hanger. Not available with Pro ELR or Peds Pro ELR.</i>	
4-Way - Out	N/C
<i>Not available with Pro ELR or Peds Pro ELR</i>	

### Heel Loops

*Not available with One Piece Flip Up Angle Adjustable Footplate*

Heel Loops (E0951)	\$
<i>Not available with Frame Widths less than 14" without Offset Hangers or Frame Widths less than 12" with Offset Hangers.</i>	
Adjustable Heel Loops (E0951)	\$
Omit Heel Loops	N/C

### Calf Straps

Velcro Adjustable (K0038)	\$
Padded Velcro Adjustable (K0038)	\$
Padded Airmesh Calf Strap (K0038)	\$

### Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
						5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Left	Left Limb Support (E1020)	\$	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5X8 Pad Left Side	7X10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5X8 Pad Right Side	7X10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

## Seat Height, Casters, and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size				
	5"	6"	6" x 2"	7"	8"
	Seat Height				
12" with Pneumatic or Pneumatic with Airless Insert Tire	13.5"				
	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17" <sup>3</sup>	17" <sup>3</sup>
	17.5"	17.5"	17.5"	17.5"	17.5"
12" with Full Poly Tire	13.5"				
	14"				
	14.5"	14.5"	14.5"		
	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16" <sup>4</sup>
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17" <sup>3</sup>	17" <sup>3</sup>
	17.5"				
12" with Drum Brake	13.5"				
	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16" <sup>4</sup>
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
16" with Pneumatic, Pneumatic with Airless Insert or Full Poly Tire	15"	15"	15"	15" <sup>4</sup>	
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
19.5" <sup>2</sup>	19.5" <sup>4</sup>	19.5" <sup>4</sup>	19.5" <sup>4</sup>		
16" With Low Poly Tire	14"	14" <sup>4</sup>	14" <sup>4</sup>		
	14.5"	14.5"	14.5"		
	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"

1 Self Propel Axle Plate selections limit tilt range

2 Not available with 12" Frame Depth

3 Not available with 12-16" Frame Depth

4 Not available with 17-20" Frame Depth

Rear Wheel Size	Caster Size				
	5"	6"	6" x 2"	7"	8"
	Seat Height				
16" with Drum Brake	15"	15"	15"		
	15.5"	15.5"	15.5"	15.5"	
	16"	16"	16"	16"	16"
	16.5"	16.5"	16.5"	16.5"	16.5"
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5" <sup>4</sup>	18.5" <sup>4</sup>	18.5"
	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"
20" with Pneumatic, Pneumatic with Airless Insert, Low Poly, Full Poly Tire Or Drum Brake	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
	19.5"	19.5"	19.5"	19.5"	19.5"
	20" <sup>4</sup>	20" <sup>4</sup>	20" <sup>4</sup>	20"	20"
	20.5"	20.5"	20.5"	20.5"	20.5"
20" Self Propel or Self Propel with Drum Brake	15.5" <sup>10</sup>	15.5" <sup>10</sup>	15.5" <sup>10</sup>	15.5" <sup>10</sup>	
	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>	16" <sup>10</sup>
	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>
	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5" <sup>4</sup>	18.5"	18.5"	18.5" <sup>4</sup>	18.5"
	19"	19"	19"	19"	19"
22" Self Propel or Self Propel with Drum Brake	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>	16.5" <sup>10</sup>
	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>	17" <sup>10</sup>
	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>	17.5" <sup>10</sup>
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
	19.5" <sup>4</sup>	19.5" <sup>2</sup>	19.5" <sup>4</sup>	19.5" <sup>4</sup>	
24" Self Propel <sup>1</sup>	17"	17"	17"	17"	17"
	17.5"	17.5"	17.5"	17.5"	17.5"
	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"
	19.5" <sup>7,8</sup>	19.5"	19.5"	19.5" <sup>9</sup>	19.5" <sup>7,8</sup>
	20" <sup>5</sup>	20" <sup>6</sup>	20" <sup>6</sup>	20"	20"
	20.5" <sup>4,11</sup>	20.5" <sup>11</sup>	20.5" <sup>11</sup>	20.5" <sup>11</sup>	20.5" <sup>11</sup>
24" Self Propel with <sup>1</sup> Drum Brake	19"	19"	19"	19"	19"
	19.5" <sup>11</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>	19.5" <sup>9,12</sup>
	20" <sup>4</sup>	20" <sup>6</sup>	20"	20"	20"

5 Not available in combination with spoke wheel and Full Poly tire

6 Not available in combination with spoke wheel, Full Poly tire and 17-20" Frame Depth

7 Not available in combination with Mag wheel and 17-20" Frame Depth

8 Not available in combination with spoke wheel, Pneumatic tires and 17-20" Frame Depth

9 17-20" Frame Depth requires spoke wheel and Full Poly tire

10 Not available with Drum Brake

11 Available only with Full Poly tire 12

12" Frame Depth requires spoke wheel and Full Poly tire

### Casters

		Caster Size				
		Price	5"	6"	7"	8"
Caster Type	.75" Lighted Rollerblade	\$				
	1" Poly	N/C				
	1" Poly Aluminum	\$				
	1" Pneumatic	\$				
	1.5" Poly	\$				
	1.5" Soft Roll Aluminum	\$				
	2" Poly	\$				

### Axles

No selection needed with Drum Brake

Fixed \_\_\_\_\_ N/C  
*Not available with Folding Base Frame and Standard Wheel Locks*  
 Quick Release \_\_\_\_\_ N/C

### Axle Plates

Standard \_\_\_\_\_ N/C  
 Self Propel Axle Plate \_\_\_\_\_ N/C

### Wheels

		Rear Wheel Size				
		12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C				
	Colored Mag <sup>1</sup>	\$				
	MAXX Spoke	\$				

<sup>1</sup> Mag color selection is on page 6

### Rear Tire Type

Low Poly (E2218) \_\_\_\_\_ N/C  
*Available only with 16" wheel size*  
 Full Poly (E2218) \_\_\_\_\_ N/C  
 Pneumatic (E2211,E2212) \_\_\_\_\_ \$  
 Pneumatic with Airless Insert (E2211, E2213) \_\_\_\_\_ \$

### Handrims

Handrims not available on 12" or 16" wheel sizes

Aluminum Anodized \_\_\_\_\_ N/C  
 Plastic Coated \_\_\_\_\_ \$  
 Projection (8 Vertical Projections) \_\_\_\_\_ \$  
*Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels*  
 Omit Handrim \_\_\_\_\_ N/C

### Handrim Option

Not available on Plastic Coated or Projection Handrims

3M Micro Grip \_\_\_\_\_ \$

## Wheel Locks

### Wheel Lock

Rear wheels must be removed to fold chair with Push to Lock, Pull to Lock and Dual wheel lock selections

Push to Lock \_\_\_\_\_ N/C  
 Pull to Lock \_\_\_\_\_ N/C  
*Not available with 12" or 16" Full Poly tires*  
 Attendant Foot Lock \_\_\_\_\_ N/C  
*Available only with 12" and 16" Mag wheels*  
 Drum Brake (E2228) \_\_\_\_\_ \$  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire or Reclining Backrest. Axle is Quick Release.*  
 Dual - Attendant Foot Lock with Push to Lock \_\_\_\_\_ \$  
*Available only with 12" and 16" Mag wheels*  
 Dual - Attendant Foot Lock with Pull to Lock \_\_\_\_\_ \$  
*Available only with 12" and 16" Mag wheels*  
 Dual - Drum Brake with Push to Lock (E2228) \_\_\_\_\_ \$  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire or Reclining Backrest. Axle is Quick Release.*  
 Dual - Drum Brake with Pull to Lock (E2228) \_\_\_\_\_ \$  
*Available only with 12" Mag, 16" Mag, or MAXX Spoke wheels. Not available with 16" Low Poly tire or Reclining Backrest. Axle is Quick Release.*  
 Omit \_\_\_\_\_ N/C

### Wheel Lock Option

Extension Handles (E0961) \_\_\_\_\_ \$  
*Available only with Push to Lock or Pull to Lock Wheel Locks*

## Armrests

Armrest Type	Length	Price
Height Adjustable Pediatric E0973 <sup>1</sup>	Desk(6.5"-9")	\$
	Desk Tall(8"-10.5")	
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973 <sup>2</sup>	Desk	\$
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>2,3</sup>	Desk	\$
	Full	

<sup>1</sup> Best fit with Frame Depths less than 14"

<sup>2</sup> Not available with 12" Frame Depth

<sup>3</sup> If requested configuration is unavailable, Standard Height T-Arm will be used

<sup>4</sup> Not available with Dynamic Back. Available only with 12" and 16" Rear Wheels.

<sup>5</sup> Not available with Dynamic Back, Reclining Back or Arm Pads

Armrest Type	Length	Price
Angle Adjustable Locking Extendable Flip Up K0020 <sup>4</sup>	Desk	\$
	Full	
Tubular Flip Up <sup>5</sup>		N/C
Omit Armrests		N/C

### Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrest</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

## Accessories

### Removable Side Guards

*Not available with Height Adjustable Armrests*

Adult Composite Side Guards (14" - 16") (K0108)	\$
Pediatric Composite Side Guard (10" - 14") (K0108)	\$

### Anti-Tips

Rear Anti-Tips (E0971)	\$
Omit	N/C

### Pelvic Positioning Belts

1" Auto Buckle Pelvic Positioning Belt (E0978)	\$
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$
Small Flex-Mount Positioning Belt (16-24") (E0978)	\$
Medium Flex-Mount Positioning Belt (19-29") (E0978)	\$
Large Flex-Mount Positioning Belt (25-37") (E0978)	\$

### Belt Mounting Kit

Belt Mounting Band Clamp	\$
Seat Belt Mounting Hardware	\$

### Spoke Guards

Clear Spoke Guards (K0065)	\$
<i>Not available with Mag wheel</i>	
Black Spoke Guards (K0065)	\$
<i>Available only with 22" and 24" Spoke wheels</i>	
Blue Monster Spoke Guards (K0065) <sup>1</sup>	\$
Robo Fun Spoke Guards (K0065) <sup>1</sup>	\$
Tutu Cute Spoke Guards (K0065) <sup>1</sup>	\$
Tuxedo Time Spoke Guards (K0065) <sup>1</sup>	\$

<sup>1</sup> Available only with 20" and 22" wheel sizes. Not available with Mag wheel

### Other Accesories

Tool Kit	\$
Removable Underseat Pouch	\$
Pediatric Ki Mobility Backpack	\$
Ki Mobility Backpack	\$
Black Neoprene Impact Guards	\$
Tilt Stop Posterior (Pair) (K0108)	\$
Tilt Stop Anterior (Pair) (K0108)	\$
<i>Available only with Rigid Base</i>	
IV Pole	\$
Canopy with Hardware (K0108)	\$

## Color

Select desired frame and accent colors. Accent color is used on side panel; if no accent color is selected, frame color will be used.

Options	Price	Colors																							
		Gloss													Matte						Metallic				
		Black	Yellow	Light Pink	Dark Pink	Glacier Blue	Electric Green	Radiant Orange	Opalescent Violet	Purple Plum	Rich Teal	Metallic Orange	Granny Smith	Candy Red	Candy Blue	Candy Purple	Tactical Green	Dry Earth Tan	White	Shadow Black	Bright Red	Royal Blue	Charcoal Gray	Iridescent Navy	Iridescent Burgundy
Frame	N/C																								
Accent	N/C																								
Mag Color																									

Cushions and Backrests

AXIOM

Ki Mobility

ADI  
A Stealth Product

VARILITE

R I O E

Cushion:

Backrest:

Other Accessories: