

Customer Details					
Client					
Assessor					
Organisation					
	ACC		Whaikaha		Other
User Weight					

Ki Mobility Rogue2 Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

		Seat Width		
		16"	18"	20"
Seat Depth	16"	<input type="checkbox"/>	<input type="checkbox"/>	
	17"	<input type="checkbox"/>	<input type="checkbox"/>	
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Express Codes

- KIROGUE21616STD
- KIROGUE21818STD
- KIROGUE21618STD
- KIROGUE22018STD
- KIROGUE21816STD

Add ROGUEEXPPACKAGE to all orders

Use KIROGUEXXXSTD where available

Express Package: (Select optional specifications and alternatives on reverse)

- Standard Seat Sling
- Omit Transit
- 2in Depth Adjustment
- 85 degree Front Frame Bend
- Sling position 0°
- 1in Tapered Frame
- Height adjustable Medium straight back posts with no Push Handle
- Back Angle set to Standard
- Height adjustable Rigidising Bar
- Tension adjust. back upholstery
- Standard Fork
- Angle Adjustable Flip Under Footplates
- 24in Maxx Spoke Performance Pneumatic Black Wheels
- 5x1.5 Soft Roll Aluminium Castors
- Ergo Lite Handrims omit thumb guard
- Rear Wheel Camber set to 2 degrees
- COG set to 0.5"
- Caster Optimised
- Rear Wheel spacing set to 1.5" to accommodate T-Arms
- Quick Release Axles
- Push to Lock Flush Wheel Locks
- T-Arms with Standard Desk Armpads
- Rear Anti Tipppers
- Iridescent Navy Frame

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height

<input type="checkbox"/>	13"	<input type="checkbox"/>	14"	(For 13" - 14" use .ULROGBUPTA1 if back upholstery required)
<input type="checkbox"/>	15"	<input type="checkbox"/>	16"	
				<input type="checkbox"/> Refer to box #1

Rear Seat Height

16W:	<input type="checkbox"/>	15"	<input type="checkbox"/>	15.5"	<input type="checkbox"/>	16"	<input type="checkbox"/>	16.5"
	<input type="checkbox"/>	17"						
18W:	<input type="checkbox"/>	15"	<input type="checkbox"/>	15.5"	<input type="checkbox"/>	16"	<input type="checkbox"/>	16.5"
	<input type="checkbox"/>	17"	<input type="checkbox"/>	17.5"	<input type="checkbox"/>	18"		
20W:	<input type="checkbox"/>	17"	<input type="checkbox"/>	17.5"	<input type="checkbox"/>	18"	<input type="checkbox"/>	18.5"
	<input type="checkbox"/>	19"	<input type="checkbox"/>	19.5"	<input type="checkbox"/>	20"		

Backrest Angle

<input type="checkbox"/>	Std	<input type="checkbox"/>	+2°	<input type="checkbox"/>	+4°	<input type="checkbox"/>	+6°
<input type="checkbox"/>	+8°	<input type="checkbox"/>	+10°	<input type="checkbox"/>	-2°	<input type="checkbox"/>	-4°
<input type="checkbox"/>	-6°	<input type="checkbox"/>	-8°	<input type="checkbox"/>	-10°		

Relaxed Back 0° +4° +9°

Seat to Footrest Length

16W:	<input type="checkbox"/>	13.5"	<input type="checkbox"/>	14"	<input type="checkbox"/>	14.5"	<input type="checkbox"/>	15"
18W:	<input type="checkbox"/>	14"	<input type="checkbox"/>	14.5"	<input type="checkbox"/>	15"	<input type="checkbox"/>	15.5"
	<input type="checkbox"/>	16"						
20W:	<input type="checkbox"/>	15"	<input type="checkbox"/>	15.5"	<input type="checkbox"/>	16"	<input type="checkbox"/>	16.5"
	<input type="checkbox"/>	17"	<input type="checkbox"/>	17.5"				

Front Seat Height

16W:	<input type="checkbox"/>	17.5"	<input type="checkbox"/>	18"	<input type="checkbox"/>	18.5"
18W:	<input type="checkbox"/>	18.5"	<input type="checkbox"/>	19"	<input type="checkbox"/>	19.5"
20W:	<input type="checkbox"/>	19.5"	<input type="checkbox"/>	20"	<input type="checkbox"/>	20.5"
<input type="checkbox"/> Refer to box #7						

Refer to box #4

Optional Specifications

Backrest Height 1

Swap to Tall (use .ULROG2BOPTL) Only use .ULROGBUPTA - if upholstery required

17" 18" 19" 20"

Backrest Upholstery 2

Omit Back Upholstery (add narrative)

Backrest Options 3

<input type="checkbox"/>	Standard Push Handle Med (.ULROGBCKSTDMED)
<input type="checkbox"/>	Standard Push Handle Med (.ULROG2BCKSTDMED)
<input type="checkbox"/>	Standard Push Handle Tall (.ULROGBCKSTD TALL)
<input type="checkbox"/>	Fold Down Push Handle Med (.ULROGBCKFLDMED)
<input type="checkbox"/>	Fold Down Push Handle Tall (.ULROGBCKFLD TALL)

Footrest Options 4

High Mount Alum. Ang. Adj. Footplates A

	5"	6"	7"	8"	9"	10"	11"	12"	13"	14"
16W										
18W										
20W										

(.ULROGFPLHM)

Alum. Ang. Adj. Footplates B

	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"
16W									
18W									
20W									

(.ULROGFPLAL or .ULROG2FPLAL1)

Side Guards (not available with T-Arms) 5

Clamp On Side Guards (.ULROGSDGDC)

Tall Armpad Short (14"- 17" deep)

Short side shld and std pad desk (.ULEXPTARSD)

Short side shld and std pad full (.ULEXPTARSF)

Camber 6

0° 4° Rogue2
(.ULROG2CAM0) (.ULROG2CAM4)

0° 4° Rogue
(.ULROGCAM0) (.ULROGCAM4)

Castor Options 7

4 x 1.5 Soft Roll Castors (.ULROGCASSR4)

For 16" Width, Front Seat Height 17"
For 18" Width, Front Seat Height 18"
For 20" Width, Front Seat Height 19"

Push Rim Options 8

Plastic coated push rim (.ULEXPHRMPC)

Standard aluminum push rim (.ULEXPHRMSALI)

Natural Fit:

Thumb grip regular (.ULEXPHRMTGREG)

Ergonomic Fit:

Ergo Lite Regular (.ULEXPHRMELREG)

Short Tabs (.ULEXPHRMSTAB)

Rear Wheel Kit - Low profile tyres 9

Rear wheel kit (KIEXPWHL02)

Includes the following:

- Spinery spox black 24"
- Low profile - Schwalbe Marathon Plus Evolution 1"
- Natural fit LT omit thumb grip

OR

24" Maxx Spoke Low Profile Schwalbe Marathon Plus Ergo Lite Omit thumb grip (.ULEXP24MAXLPFEL)

Plastic coated push rim (.ULEXPHRMPC)

Standard aluminum push rim (.ULEXPHRMSALI)

Natural Fit:

Standard gr handrim for MAXX SP (.ULCATHRMNF02)

Ergonomic Fit:

Ergo Lite Regular (.ULEXPHRMELREG)

Short Tabs (.ULEXPHRMSTAB)

Wheel Locks 10

Low Profile Scissor Locks (.ULROGWLKLPS)

Push To Lock (.ULROGWLKPTL)

Wheel Lock Extension (.ULROGWLKEXT) Only with Push to Lock

Rogue 2 Options 11

User Activated Anti Tips (No more than 2" dump .ULROG2UANTTP)

Accessories 12

Clear Spoke Guards (.ULROGSPGCL24)

Impact Guards (.ULROGIMG)

Bodypoint Belt Options:

Bodypoint 2pt belt – Centre Pull Flat Mount (BPHB205-M46-A1) + (BPHW200-90-2)

Bodypoint 4pt belt Centre Pull Flat Mount (BPHB405-M46-A1) + (BPHW200-90-2)

Bodypoint Evoflex Pelvic Stabiliser (BPEB205-M46) + (BPHW200-90-2)

Bodypoint Aeromesh Calf Strap (BPSP103M)

Blank area for notes.