

Ki Mobility Rogue Express Script



Customer Details					
Client					
Assessor					
Organisation					
	ACC		DSS		Other
User Weight					

Ki Mobility Rogue2 Express Script

If you only require specifications that fall under the Standard Express Package, you do not need to fill in the reverse page.

		Seat Width		
		16"	18"	20"
Seat Depth	16"	☐	☐	☐
	17"	☐	☐	☐
	18"	☐	☐	☐
	19"	☐	☐	☐
	20"	☐	☐	☐

Express Codes

☐ KIROGUE21616STD	☐ KIROGUE21818STD
☐ KIROGUE21618STD	☐ KIROGUE22018STD
☐ KIROGUE21816STD	

Add ROGUEEXPPACKAGE to all orders

☐ Can use Rogue 10 where available

Express Package: (Select optional specifications and alternatives on reverse)

- Standard Seat Sling
- Omit Transit
- 2in Depth Adjustment
- 85 degree Front Frame Bend
- Sling position 0°
- 1in Tapered Frame
- Height adjustable Medium straight back posts with no Push Handle
- Back Angle set to Standard
- Height adjustable Rigidising Bar
- Tension adjust. back upholstery
- Standard Fork
- Angle Adjustable Flip Under Footplates
- 24in Maxx Spoke Performance Pneumatic Black Wheels
- 5x1.5 Soft Roll Aluminium Castors
- Ergo Lite Handrims omit thumb guard
- Rear Wheel Camber set to 2 degrees
- COG set to 0.5"
- Caster Optimised
- Rear Wheel spacing set to 1.5" to accommodate T-Arms
- Quick Release Axles
- Push to Lock Flush Wheel Locks
- T-Arms with Standard Desk Armpads
- Rear Anti Tippers
- Iridescent Navy Frame

Please note: Express Package includes some upcharge products.

Choose Your Specifications (No additional codes are required on the below parts)

Backrest Height

☐ 13"	☐ 14"	(For 13" - 14" use .ULROGBUPTA1 if back upholstery required)
☐ 15"	☐ 16"	☐ Refer to box #1

Backrest Angle

☐ Std	☐ +2°	☐ +4°	☐ +6°
☐ +8°	☐ +10°	☐ -2°	☐ -4°
☐ -6°	☐ -8°	☐ -10°	

Relaxed Back	☐ 0°	☐ +4°	☐ +9°
--------------	-----------------------------	------------------------------	------------------------------

*Relaxed back only available on Rogue2

Front Seat Height

16W:	☐ 17.5"	☐ 18"	☐ 18.5"
18W:	☐ 18.5"	☐ 19"	☐ 19.5"
20W:	☐ 19.5"	☐ 20"	☐ 20.5"

☐ Refer to box #7
--

Rear Seat Height

16W:	☐ 15"	☐ 15.5"	☐ 16"	☐ 16.5"
	☐ 17"			
18W:	☐ 15"	☐ 15.5"	☐ 16"	☐ 16.5"
	☐ 17"	☐ 17.5"	☐ 18"	
20W:	☐ 17"	☐ 17.5"	☐ 18"	☐ 18.5"
	☐ 19"	☐ 19.5"	☐ 20"	

Seat to Footrest Length

16W:	☐ 13.5"	☐ 14"	☐ 14.5"	☐ 15"
18W:	☐ 14"	☐ 14.5"	☐ 15"	☐ 15.5"
	☐ 16"			
20W:	☐ 15"	☐ 15.5"	☐ 16"	☐ 16.5"
	☐ 17"	☐ 17.5"		

☐ Refer to box #4
--

Optional Specifications																																													
Backrest Height 1	Push Rim Options (Standard Wheels) 8																																												
Swap to Tall (use .ULROG2BOPTL) Only use .ULROGBUPTA - if upholstery required 17" 18" 19" 20"	Plastic coated push rim (.ULEXPHRMPC) Standard aluminum push rim (.ULEXPHRMSALI)																																												
Backrest Upholstery 2	Natural Fit: 9																																												
Omit Back Upholstery (add narrative)	Thumb grip regular (.ULEXPHRMTGREG)																																												
Backrest Options 3	Ergonomic Fit: 10																																												
Standard Push Handle Med (.ULROGBCKSTDMED) Standard Push Handle Med (.ULROG2BCKSTDMED) Standard Push Handle Tall (.ULROGBCKSTDALL) Fold Down Push Handle Med (.ULROGBCKFLDMED) Fold Down Push Handle Tall (.ULROGBCKFLDTALL)	Ergo Lite Regular (.ULEXPHRMELREG) Short Tabs (.ULEXPHRMSTAB)																																												
Footrest Options 4	Rear Wheel Kit - Low profile tyres 9																																												
☐ High Mount Alum. Ang. Adj. Footplates A <table border="1" style="border-collapse: collapse; text-align: center; width: 100%;"> <tr> <th></th><th>5"</th><th>6"</th><th>7"</th><th>8"</th><th>9"</th><th>10"</th><th>11"</th><th>12"</th><th>13"</th><th>14"</th></tr> <tr> <td>16W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr> <td>18W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr> <td>20W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> </table> (.ULROGFPLHM)		5"	6"	7"	8"	9"	10"	11"	12"	13"	14"	16W											18W											20W											Rear wheel kit (KIEXPWHL02) Includes the following: • Spinerge spox black 24" • Low profile - Schwalbe Marathon Plus Evolution 1" • Natural fit LT omit thumb grip
	5"	6"	7"	8"	9"	10"	11"	12"	13"	14"																																			
16W																																													
18W																																													
20W																																													
☐ Alum. Ang. Adj. Footplates B <table border="1" style="border-collapse: collapse; text-align: center; width: 100%;"> <tr> <th></th><th>13.5"</th><th>14"</th><th>14.5"</th><th>15"</th><th>15.5"</th><th>16"</th><th>16.5"</th><th>17"</th><th>17.5"</th></tr> <tr> <td>16W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr> <td>18W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> <tr> <td>20W</td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr> </table> Rogue .ULROGFPLAL Rogue2 .ULROG2FPLAL1		13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	16W										18W										20W										24" Maxx Spoke Low Profile Schwalbe Marathon Plus Ergo Lite Omit thumb grip (.ULEXP24MAXLPFEL) Plastic coated push rim (.ULEXPHRMPC) Standard aluminum push rim (.ULEXPHRMSALI)				
	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"																																				
16W																																													
18W																																													
20W																																													
☐ OR	Natural Fit: Standard gr handrim for MAXX SP (.ULCATHRMNF02)																																												
Side Guards (not available with T-Arms) 5	Ergonomic Fit: 11																																												
Clamp On Side Guards (.ULROGSDGDC)	Ergo Lite Regular (.ULEXPHRMELREG) Short Tabs (.ULEXPHRMSTAB)																																												
Camber 6	Wheel Locks 10																																												
0° (.ULROG2CAM0) 4° (.ULROG2CAM4) Rogue2	Low Profile Scissor Locks (.ULROGWKLKPS) Push To Lock (.ULROGWLKPTL) Wheel Lock Extension (.ULEXPWLKEXT) Only with Push to Lock																																												
0° (.ULROGCAM0) 4° (.ULROGCAM4) Rogue	Anti Tips 11																																												
☐ Tall Armpad Short (14"- 17" deep) Short side shld and std pad desk (.ULEXPTARSD) Short side shld and std pad full (.ULEXPTARSF)	User Activated Anti Tips (No more than 2" dump .ULROG2UANTTP)																																												
Castor Options 7	Accessories 12																																												
☐ 4 x 1.5 Soft Roll Castors (.ULROGCASSR4) For 16" Width, Front Seat Height 17" For 18" Width, Front Seat Height 18" For 20" Width, Front Seat Height 19"	Clear Spoke Guards (.ULROGSPGCL24) Impact Guards (.ULROGIMG)																																												
☐ Bodypoint Belt Options: Bodypoint 2pt belt – Centre Pull Flat Mount (BPHB205-M46-A1) + (BPHW200-90-2) Bodypoint 4pt belt Centre Pull Flat Mount (BPHB405-M46-A1) + (BPHW200-90-2) Bodypoint Evoflex Pelvic Stabiliser (BPB205-M46) + (BPHW200-90-2) Bodypoint Aeromesh Calf Strap (BPSP103M)	☐ Side Guards (not available with T-Arms) 5																																												

