




LECKEY BeMe^{CE}

Name	
Date	
Email	
Contact Number	
Purchase Order	
Invoice Address	
Delivery Address	

Packages



Option Code	Description	Price
Size		
151BF40	BeMe Size 1	
151BF41	BeMeSize 2	
151BF42	BeMe Size 3	
151BF43	BeMe Size 4	
Colour - Note, all trial equipment will come in blue. Other colours available for purchase only.		
151YJ18	Blue	
151YJ19	Purple	
151YJ20	Red	
151YJ21	Grey	
Package		
151QE38	Includes - Basic Seat Basic Backrest 1 Piece Footplate Armrests Hip guides 4 point belt	
151QE39	Includes - Contoured Seat Moderate Backrest Split Footplates & Sandals Flip Away Slim Laterals Butterfly Harness Armrests BeMe Headrest 4 point belt	
151QE40	Includes - Split Seat Advanced Backrest Split Footplates & Sandals BeMe Headrest Flip Away Slim Laterals Hip guides Butterfly Harness Armrests 4 point belt 4 x Gables	

Chassis		
		
151NB118	Multi height and tilt Indoor - Manual Mechanism	
151NB119	Multi height and tilt Indoor/Outdoor - Manual Mechanism	
151NB120	Multi height and tilt Indoor - Electric Mechanism	
151NB121	Multi height and tilt Indoor/Outdoor - Electric Mechanism	
Push Handles		
151PH1	Push Handles	
Interfaces		
		
151NB122	Interface for Size 1 BeMe on to Zippie Iris, Quickie Iris, Quickie SR45	
151NB123	Interface for Size 1 BeMe on to Q300 M -std	
151NB124	Interface for Size 1 BeMe on to Q300 M - lift & tilt	
151NB125	User control mount bracket (Joystick Mount)	
Chest Harness		
		
151NH121	Dynamic Butterfly Harness	
Tray		
		
1517A11	Tray	
1517A12	Padded Tray Cover	
1517A13	Tray Elbow Blockers	
Leg Support		
		
151F138	Cushion for One-Piece Footplate (Only Available on Basic Package)	

Foot Support		
  		
151DI86	Ankle Huggers	
150-1852	Extra Small Booties	
151DI88	Small Booties	
151DI89	Sandals Small	
151DI90	Sandals Medium	
151DI91	Sandals Large	
151DI92	Small Ankle Wraps (Only with Sandals)	
151DI93	Large Ankle Wraps (Only with Sandals)	
Headrests		
 		
151HE52	BeMe Headrest	
151HE53	Leckey Contoured Headrest	
151HE54	Leckey Contoured Occipital Headrest	
151HE55	Leckey Flat Headrest	
151HE56	Leckey Flat Headrest Laterals	
Miscellaneous		
1518A129	Communication mount	
1518A132	25mm Spacer to Reduce Seat Depth	
1518A133	25mm Spacer to Reduce Seat Depth	
Assemble		
151HA7	Assemble at Sunrise	
Technical Information		

	Size 1	Size 2	Size 3	Size 4
Seat Depth [mm] *	120-270	225-350	320-430	400-540
Seat Width Maximum (No Hip Guides) [mm]	280	325	410	475
Femoral Width on Split Seat [mm]		90-125	110-160	135-170
Maximum Weight Limit [kg] **	35	50	70	80
Back Height [mm]	300-400	350-470	400-530	475-650
Standard Hip Guide Width [mm]	145-215	175-260	245-350	290-405
Slim Hip Guide Width [mm]	170-240	200-285	270-375	315-430
Lateral Width [mm]	130-230	200-300	225-350	350-450
Inset Lateral Width [mm]	80-180	150-250	175-300	300-400
Backrest Angle [degrees]			+10 to -25	
Maximum Weight Limit [kg] **	35	50	70	80