







# Leckey Horizon

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote  Order




Purchase Order	
----------------	--



Item Code	Description	Quantity
<b>Chassis</b>		
		
LHZ/1/177-EXP	Size 1 Horizon Chassis	
LHZ/2/177-EXP	Size 2 Horizon Chassis	
LHZ/3/177-EXP	Size 3 Horizon Chassis	
LKS1127	Horizon Stander Charger	
<b>Spine</b>		
		
LHZ/1/198-GRB-EUR	Horizon Size 1 Extrusion with Chest Pad Hip Pad Chest Laterals - Grey	
LHZ/1/198-GB-EUR	Horizon Size 1 Extrusion with Chest Pad Hip Pad Chest Laterals - Green	
LHZ/1/198-BB-EUR	Horizon Size 1 Extrusion with Chest Pad Hip Pad Chest Laterals - Blue	
LHZ/2/198-GRB-EUR	Horizon Size 2 Extrusion with Chest Pad Hip Pad Chest Laterals - Grey	
LHZ/2/198-GB-EUR	Horizon Size 2 Extrusion with Chest Pad Hip Pad Chest Laterals - Green	
LHZ/2/198-BB-EUR	Horizon Size 2 Extrusion with Chest Pad Hip Pad Chest Laterals - Blue	
LHZ/3/198-GRB-EUR	Horizon Size 3 Extrusion with Chest Pad Hip Pad Chest Laterals - Grey	
LHZ/3/198-GB-EUR	Horizon Size 3 Extrusion with Chest Pad Hip Pad Chest Laterals - Green	
LHZ/3/198-BB-EUR	Horizon Size 3 Extrusion with Chest Pad Hip Pad Chest Laterals - Blue	
<b>Pelvic Support</b>		
		
	Hip Laterals	Pelvic Belt
LHZ/1/164	Horizon Size 1 Hip Laterals	
LHZ/2/164	Horizon Size 2 Hip Laterals	
LHZ/3/164	Horizon Size 3 Hip Laterals	
LHZ/1/181	Horizon Size 1 SUPINE Pelvic Belt	
LHZ/2/181	Horizon Size 2 SUPINE Pelvic Belt	
LHZ/3/181	Horizon Size 3 SUPINE Pelvic Belt	
<b>Tray</b>		
		
	Supine Stander tray	Prone Stander tray
	Padded tray insert	
LHZ/1/158	Horizon Size 1 PRONE TRAY	
LHZ/1/143	Horizon Size 1 SUPINE TRAY	
LHZ/2/158	Horizon Size 2-3 PRONE TRAY	
LHZ/2/143	Horizon Size 2-3 SUPINE TRAY	
LHZ/1/109	Padded Tray Insert for Size 1 Tray	
LHZ/2/109	Padded Tray Insert for Size 2/3 Tray	
101-3690	Horizon Size 1 Mygo Stander Tray for Horizon	
101-3691	Horizon Size 2-3 Mygo Stander Tray for Horizon	

Chest Support (Laterals Included with Extrusion)			
			
	101-2120	Hyperextension Pads for Laterals	
Knees			
			
		Knee Pads	Knee Extension Brackets
	LHZ/1/167	Horizon Size 1 Knees (pair) Short adjustment = - 10 to + 45mm or long adjustment = - 10 to + 85mm	
	LHZ/2/167	Horizon Size 2-3 Knees (pair) Short adjustment = - 10 to + 45mm or long adjustment = - 10 to + 85mm	
	LHZ/4/167	Horizon Size 4 Large Knees (pair) Short adjustment = - 10 to + 45mm or long adjustment = - 10 to + 85mm	
	LHZ/1/104	Horizon Size 1 Wide Knees (pair)	
	LHZ/2/104	Horizon Size 2-3 Wide Knees (pair)	
	101-2224	Medium Knee Adjustment Bracket (125mm-175mm adjustment)	
	101-2225	Long Knee Adjustment Bracket (175mm-225mm adjustment)	
Horizon Headrests			
			
		Flat Headrest	Head lateral supports
			Extended Head Rest Brackets
	LHZ/1/169	Horizon Size 1 Headsupport complete	
	LHZ/2/169	Horizon Size 2-3 Headsupport complete	
	LHZ/2/117	Horizon Size 2-3 Head Lateral Supports	
	101-1959	Headrest 4 Inch Extension Bracket	
Whitmyer Plush			
			
	PRO-P1W	Plush Headrest, 10 Inch with 8 Inch PRO Bar	
	PRO-P2W	Plush Headrest, 14 Inch with 8 Inch PRO Bar	
	PRO-P4W	Plush Headrest, 19 Inch with 8 Inch PRO Bar	
	PRO-AP1L	Adjust-A-Plush Narrow Headrest with 8 Inch PRO Bar	
	PRO-AP2L	Adjust-A-Plush Standard Headrest with 8 Inch PRO Bar	
Whitmyer Contoured Cradle			
			
	PRO-CCS2	Small Contoured Cradle Head on 8 Inch PRO Bar	
	PRO-CCS3	Medium Contoured Cradle Head on Pro 8 Inch PRO Bar	
	PRO-CCS4	Large Contoured Cradle Head on 8 Inch PRO Bar	
Whitmyer Heads Up			
			

	PRO-HU-ISSCL-05	Heads Up Short Infant (Size 1 - 8.3 Cm) including chin prompt	
	PRO-HU-SSSCL-06	Heads Up Short Child (Size 2 - 10.1 Cm) including chin prompt	
	PRO-HU-MSSCL-08	Heads Up Short Teenager (Size 3 - 11.6 Cm) including chin prompt	
	PRO-HU-LSSCL-010	Heads Up Short Adult (Size 4 - 14.7 Cm) including chin prompt	
	PRO-HU-ITSCL	Heads Up Tall Infant (Size 1 - 8.3 Cm) including chin prompt	
	PRO-HU-STSCl	Heads Up Tall Child (Size 2 - 10.1 Cm) including chin prompt	
	PRO-HU-MTSCl	Heads Up Tall Teenager (Size 3 - 11.6 Cm) including chin prompt	
	PRO-HU-LTSCl	Heads Up Tall Adult (Size 4 - 14.7 Cm) including chin prompt	
<b>Whitmyer Complex</b>			
			
	PRO-3W	Paediatric S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
	PRO-4W	Small Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
	PRO-1W	Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
<b>Whitmyer Forehead Support</b>			
			
	D1010	Whitmyer Bracket for Cradle headrest to attach forehead support	
	DFS3	DYNAMIC FOREHEAD ASSEMBLY - Paediatric	
	DFS1	DYNAMIC FOREHEAD ASSEMBLY - Adult	
	D1017	DYNAMIC FOREHEAD STRAP (NEOPRENE AGAINST SKIN)	
	D1018	DYNAMIC FOREHEAD STRAP (TOWELLING AGAINST SKIN)	
<b>Whitmyer Bars and Accessories</b>			
	PRC7030	Whit 8" vertical/horizontal PRO mount	
	PRC7070	PRO 63mm Offset Bar	
	M9450	Dynamic Headrest Kit	
	LFC25L	Single Lateral Facial Support for 3 Pad System	
	LFC25-CSL	Single Lateral Facial Support for Contoured Cradle	
	LFC25-PLL	Single Lateral Facial Support for Plush	
	M5200	Swingaway System for headrests	
<b>Foot Support</b>			
 			
<div style="display: flex; justify-content: space-around; width: 100%;"> <span>Flip Up Footplate</span> <span>3D Footplate</span> </div>			
    			
<div style="display: flex; justify-content: space-around; width: 100%;"> <span>Sandals</span> <span>Open-back sandals</span> <span>Booties</span> <span>Ankle Huggers</span> <span>Ankle Wraps</span> </div>			
	LK101-2675	Horizon Basic Footplate Assembly Size 2 (Limited stock available, discontinued)	
	LHZ/1/161	Horizon Size 1 Flip-up Footplate Assy	
	LHZ/2/161	Horizon Size 2-3 Flip-up Footplate Assy	
	LHZ/3/160	Horizon Size 1-3 3D Footplate Assy INC. SANDALS	
	152-1600	Sandals Size 1	
	152-2600	Sandals Size 2	
	152-3600	Sandals Size 3	
	LHZ/1/100	Size 1 open-back sandals	
	LHZ/2/100	Size 2 open-back sandals	
	LHZ/3/100	Size 3 open-back sandals	
	LHZ/4/100	Size 4 open-back sandals	
	AH01	Ankle Huggers	
	173-2305	Small Ankle Wraps (use with sandals)	
	173-3305	Large Ankle Wraps (use with sandals)	

Accessories			
			
Spine Extension		Lateral Spacer	
	LHZ/1/166	Horizon Bracket to extend outer shaft/spine by 150mm	
	RAF-4-1723	Horizon 1 Inch Lateral Spacer	
Harness			
			
	155-N001-L600	Small Butterfly Harness & attachment straps	
	155-N001-L601	Medium Butterfly Harness & attachment straps	
	155-N001-L602	Large Butterfly Harness & attachment straps	
	155-N001-L603	X Large Butterfly Harness & attachment straps	

	Size 1	Size 2	Size 3
Age	4-10	8-14	14-Adult
User Height	1000mm-1400mm (39.3"-55.1")	1400mm-1800mm (55.1"-70.8")	1550mm-1800mm (61"-70.8")
Max. User Weight	50kg (110lbs)	80kg (176lbs)	100kg (220lbs)
Shoulder Height (footplate to top of pad)	820mm-1150mm (32.2"-45.2")	1000mm-1450mm (39.3"-57")	1000mm-1450mm (39.3"-57")
Hip Pad Height (footplate to top of pad)	450mm-750mm (17.7"-29.5")	750mm- 950mm (29.5"-37.4")	750mm-1150mm (29.5"-45.2")
Hip Pad Width	175mm-270mm (6.8"-10.6")	210mm-390mm (8.2"-15.3")	310mm-470mm (12.2"-18.5")
Knee Height (footplate to middle of pad)	250mm- 450mm (9.8"-17.7")	400mm-550mm (15.7"-21.6")	400mm-650mm (15.7"-25.5")
Width Between Chest Laterals	200mm-300mm (7.8"-11.8")	240mm-380mm (9.4"-14.9")	340mm-460mm (13.3"-18.1")
Footplate Angle Adjustment +/-	10°	10°	10°
Angle Adjustment	0°-90°- Prone 0°-90°- Supine	0°-90°- Prone 0°-90°- Supine	0°-90°- Prone 0°-90°- Supine
Tray Angle Adjustment	0°-40°- Prone +/- 20°/30°- Supine	0°-40°- Prone +/- 20°/30°- Supine	0°-40°- Prone +/- 20°/30°- Supine

**Comments**