



**ALLIED
MEDICAL**

Leckey MyWay+ and MyWay Pedal

Name:	
Date:	
Email:	
Contact Number:	
Invoice Address:	
Delivery Address:	

Purchase Order	
-----------------------	--

☐ Quote

☐ Order



Item Code		Description	Quantity
Frame			
	LK209-1600	MyWay+ S1 MyWay Frame	
	LK209-2600	MyWay+ S2 MyWay Frame	
	LK209-3600	MyWay+ S3 MyWay Frame	
Harness			
			
	LK173-4100-04	MyWay S1 Purple MyWay Harness	
	LK173-4200-07	MyWay S2 Green MyWay Harness	
	LK173-4300-02	MyWay S3 Blue MyWay Harness	
	LK173-4400-09	MyWay S4 Red MyWay Harness	
	LK173-4500-06	MyWay S5 Orange MyWay Harness	
	LK173-4600-08	MyWay S6 Grey MyWay Harness	
	LK173-1800	MyWay Hoist Straps (one size)	
Headrest			
			
	LK173-1680	MyWay Headrest	
	LK209-1680	MyWay+ Complex Headrest	
User Handles			
			
	LK209-1667	MyWay+ Size 1 User Handles (Pair)	
	LK209-2667	MyWay+ Size 2 User Handles (Pair)	
	LK209-3667	MyWay+ Size 3 User Handles (Pair)	
	LK209-1740	MyWay+ Size 1 Complex User Handles incl. Ball Grip (Pair)	
	LK209-2740	MyWay+ Size 2 Complex User Handles incl. Ball Grip (Pair)	
	LK209-3740	MyWay+ Size 3 Complex User Handles incl. Ball Grip (Pair)	
	LK209-1319	MyWay+ Cylinder Grip (Pair)	
Adult Handles			
			
	LK209-1691	MyWay+ S1 Adult Handles	
	LK209-2691	MyWay+ S2 Adult Handles	
Ankle Guides			
			
	LK209-1700	MyWay+ S1 Ankle Guides (L&R)	
	LK209-2700	MyWay+ S2 Ankle Guides (L&R)	
	LK209-3700	MyWay+ S3 Ankle Guides (L&R)	

Saddle			
			
	LK173-1671	MyWay S1 Saddle (recommended for MyWay S1 & S2 frames)	
	LK173-2671	MyWay S2 Saddle (recommended for MyWay S2 & S3 frames)	
	LK209-1671	MyWay+ S1 Narrow Saddle (recommended for MyWay S1 & S2 frames)	
	LK209-2671	MyWay+ S2 Narrow Saddle (recommended for MyWay S2 & S3 frames)	
	LK209-1679	MyWay+ S1 Drop (fits with MyWay S1 saddle only)	
	LK209-2679	MyWay+ S2 Drop (fits with MyWay S2 saddle only)	
	LK209-3679	MyWay+ S1 Drop (fits with MyWay S2 saddle only)	
MyWay Pedal (Size 2 and 3 MyWay only)			
			
	LK173-1750	Smart Pedal Walker (one size) requires sandals	
	LK152-1600	Sandals Size 1	
	LK152-2600	Sandals Size 2	
	LK152-3600	Sandals Size 3	
	LK173-2305	Small Ankle Wraps (use with sandals)	
	LK173-3305	Large Ankle Wraps (use with sandals)	

Important Note - MyWay accessories not shown on this order form are no longer compatible or available to order with MyWay+

Frame Size	1	2	3			
Age (approx.)	1-5	4-11	9-16			
User Shoulder Height (to floor) (mm)	590-870	870-1190	1070-1490			
User Saddle Height (inside leg) (mm)	240-410	400-610	550-810			
Max. User Weight (kg)	25	50	80			
Frame Weight (kg)	14	7.5	11			
Upper Support Unit Weight (kg)	Size 1 frame & upper support unit do not separate	10	11			
Frame Width (mm)	570	670	720			
Frame Length (mm)	840	930	1180			
Harness Size	1	2	3	4	5	6
Age (approx.)	1-4	2-6	5-9	8-12	12-15	Up to 16
Waist at Navel Level (cm)	46-56	50-60	55-65	60-70	70-80	80-90
Max. User Weight (kg)	50	50	50	80	80	80
Colour	Purple	Green	Blue	Red	Orange	Grey

Comments