



# Leckey Pronestander

<b>Name</b>	
<b>Date:</b>	
<b>Email</b>	
<b>Contact Number</b>	
<b>Invoice Address</b>	
<b>Delivery Address</b>	

Quote       Order

Purchase Order	
----------------	--



Item Code	Description	Quantity
<b>Frame</b>		
LPS/1-MW	Size 1 Pronestander - Meadow	
LPS/2-MB	Size 2 Pronestander - Mondeo Blue	
LPS/3-MB	Size 3 Pronestander - Mondeo Blue	
<b>Pelvic Band</b>		
LPS/1/B-MW	Pronestander Pelvic Band - size 1 - Meadow	
LPS/2/B-MB	Pronestander Pelvic Band - size 2 - Mondeo Blue	
LPS/3/B-MB	Pronestander Pelvic Band - size 3 - Mondeo Blue	
LFS/1/B-MW	Freestander Pelvic Band - size 1 - Meadow	
LFS/2/B-MB	Freestander Pelvic Band - size 2 - Mondeo Blue	
LFS/3/B-MB	Freestander Pelvic Band - size 3 - Mondeo Blue	
<b>Soft Knees</b>		
LPS/1/K-MW	Freestander/Pronestander Upholstered Kneecups - size 1 - Meadow	
LPS/2/K-MB	Freestander/Pronestander Upholstered Kneecups - size 2 - Mondeo Blue	
LPS/3/K-MB	Freestander/Pronestander Upholstered Kneecups - size 3 - Mondeo Blue	
LPS/4/K-MB	Freestander/Pronestander Upholstered Kneecups - size 4 - Mondeo Blue	

### Foot Supports



Sandals



Open-back sandals



Booties



Ankle Huggers



Ankle Wraps

	152-1600	Sandals Size 1	
	152-2600	Sandals Size 2	
	152-3600	Sandals Size 3	
	LHZ/1/100	Size 1 open-back sandals	
	LHZ/2/100	Size 2 open-back sandals	
	LHZ/3/100	Size 3 open-back sandals	
	LHZ/4/100	Size 4 open-back sandals	
	150-1852	Extra Small Booties	
	AH01	Ankle Huggers	
	173-2305	Small Ankle Wraps (use with sandals)	
	173-3305	Large Ankle Wraps (use with sandals)	
	117-787	1Inch Footplate Raiser for Split Footplate	

Technical data	Size 1		Size 2		Size 3	
Age (years)	1.5 – 5		4 – 10		9 – 18	
Body height (cm/inches)	Min.	80/31.5	Min.	110/43.3	Min.	140/55.1
	Max.	115/45.3	Max.	145/57	Max.	185/72.8
User weight (kg/lbs)	60/132		80/176		90/198	
Chest pad height (cm/inches)	Min.	47/18.5	Min.	68/26.8	Min.	97/38.4
	Max.	68/26.8	Max.	98/38.4	Max.	114/53.1
Width between chest supports (cm/inches)	Min.	15/6	Min.	20/8	Min.	25/10
	Max.	25/10	Max.	31/12	Max.	36/14

### Comments