






Leckey Squiggles + / TT

Name	
Date	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



Item Code	Description	Quantity
Frame and Chassis		
		
184-764	Squiggle + Pivot Chassis (gas operated - can be adjusted with child in frame)	
184-600	Squiggle + Stander Support Shell (central support includes hip & chest support with sternum prompt chest & hip lateral supports kneecups footplate and tray attachment)	
Cover Pack		
Includes 2 x support harnesses, 2 x support cushions, 2 x lateral covers, 1 x pair of kneepads		
		
184-1800-09VIN	Cover Pack (Size 1 Knees)- Red Vinyl	
184-1800-09	Cover Pack (Size 1 Knees) - Red	
184-1800-04	Cover Pack (Size 1 Knees)- Purple	
184-1800-02	Cover Pack (Size 1 Knees)- Blue	
184-2800-09VIN	Cover Pack (Size 2 Knees) - Red Vinyl	
184-2800-09	Cover Pack (Size 2 Knees)- Red	
184-2800-04	Cover Pack (Size 2 Knees)- Purple	
184-2800-02	Cover Pack (Size 2 Knees)- Blue	
Tray		
		
118-619	Squiggles Stander Activity Tray	
188-1854-05	Squiggles Black Padded Tray Cover	

Foot Supports



Sandals



Booties



Ankle Huggers



Ankle Wraps

152-1600

Sandals Size 1

AH01

Ankle Huggers

150-1852

Extra Small Booties

173-2305

Small Ankle Wraps (use with sandals)

117-787

1Inch Footplate Raisers for Split Footplate

Head Supports



Contoured



Contoured Occipital



Flat



Laterals (for use with flat)

137-699-05

Contoured headrest Assembly and Cushion

137-699-05VIN

Contoured headrest Assembly and Cushion - Vinyl

117-699-05

Contoured Occipital Headrest

117-699-05VIN

Contoured Occipital Headrest - Vinyl

117-851-05

Flat Headrest Assembly and Cushion

117-851-05VIN

Flat Headrest Assembly and Cushion - Vinyl

120-893-05

Flat Headrest Laterals with covers

120-893-05VIN

Flat Headrest Laterals with covers Vinyl

Whitmyer Plush



PRO-P1W

Plush Headrest, 10 Inch with 8 Inch PRO Bar

PRO-P2W

Plush Headrest, 14 Inch with 8 Inch PRO Bar

PRO-P4W

Plush Headrest, 19 Inch with 8 Inch PRO Bar

PRO-AP1L

Adjust-A-Plush Narrow Headrest with 8 Inch PRO Bar

PRO-AP2L

Adjust-A-Plush Standard Headrest with 8 Inch PRO Bar

Whitmyer Contoured Cradle



PRO-CCS2

Small Contoured Cradle Head on 8 Inch PRO Bar

PRO-CCS3

Medium Contoured Cradle Head on Pro 8 Inch PRO Bar

PRO-CCS4

Large Contoured Cradle Head on 8 Inch PRO Bar

Whitmyer Heads Up



PRO-HU-ISSCL-05

Heads Up Short Infant (Size 1 - 8.3 Cm) including chin prompt

PRO-HU-SSSCL-06

Heads Up Short Child (Size 2 - 10.1 Cm) including chin prompt

PRO-HU-MSSCL-08




Heads Up Short Teenager (Size 3 - 11.6 Cm) including chin prompt

PRO-HU-LSSCL-010

Heads Up Short Adult (Size 4 - 14.7 Cm) including chin prompt

PRO-HU-ITSCL

Heads Up Tall Infant (Size 1 - 8.3 Cm) including chin prompt

	PRO-HU-STSC	Heads Up Tall Child (Size 2 - 10.1 Cm) including chin prompt	
	PRO-HU-MTSC	Heads Up Tall Teenager (Size 3 - 11.6 Cm) including chin prompt	
	PRO-HU-LTSC	Heads Up Tall Adult (Size 4 - 14.7 Cm) including chin prompt	
Whitmyer Complex			
			
	PRO-3W	Paediatric S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
	PRO-4W	Small Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
	PRO-1W	Adult S.O.F.T. Three Pad Head Support with 8" PROvertical mounting bar	
Whitmyer Forehead Support			
			
	D1010	Whitmyer Bracket for Cradle headrest to attach forehead support	
	DFS/3	DYNAMIC FOREHEAD ASSEMBLY - Paediatric	
	DFS/1	DYNAMIC FOREHEAD ASSEMBLY - Adult	
	D1017	DYNAMIC FOREHEAD STRAP (NEOPRENE AGAINST SKIN)	
	D1018	DYNAMIC FOREHEAD STRAP (TOWELLING AGAINST SKIN)	
Whitmyer Bars and Accessories			
	PRC7030	Whit 8" vertical/horizontal PRO mount	
	PRC7070	PRO 63mm Offset Bar	
	M9450	Dynamic Headrest Kit	
	LFC25L	Single Lateral Facial Support for 3 Pad System	
	LFC25-CSL	Single Lateral Facial Support for Contoured Cradle	
	LFC25-PLL	Single Lateral Facial Support for Plush	
	M5200	Swingaway System for headrests	
Targetted Training Conversion Kit (includes dynamic base and wraparound supports)			
			
	188-680-02	Squiggle + to Squiggle TT conversion kit blue	
	188-680-04	Squiggle + to Squiggle TT conversion kit purple	
	188-680-09	Squiggle + to Squiggle TT conversion kit red	
	188-680-09VIN	Squiggle + to Squiggle TT conversion kit red vinyl	

Age (approx)	1-5
Max User Weight	22kg/48.4lbs
User Height	Min 730mm / 28.5 inches Max 1110mm / 43.5 inches
Chest Support	Angle adjustment +/- 30° Depth adjustment 50mm / 2 inches
Chest Support Height (footplate to armpit)	Min 460mm / 18.1 inches Max 820mm / 32 inches
Chest Width (distance between laterals)	Min 160mm / 6.3 inches Max 230mm / 9 inches
Hip Support (height from footplate to mid buttock)	Min 300mm / 11.8 inches Max 530mm / 20.8 inches
Distance between midline of knees	Min 140mm / 5.5 inches Max 215mm / 8.4 inches
Footplate Angle	Plantarflexion 10° / Dorsiflexion 10°

Comments