







Leckey Squiggles Seat

Name	
Date:	
Email	
Contact Number	
Invoice Address	
Delivery Address	

Quote Order

Purchase Order	
----------------	--



Item Code	Description	Quantity
Seat		
		
120-610-MIN	Squiggles Seat Squiggles Shell-Configured in Minimum Seat Depth Set-up	
Interface		
		
194-600	Squiggles Seat Interface for Sunrise bases	
High Low Base		
 		
120-800	Squiggles Seat Hi low chassis - foot pedal	
122-710	Mygo/Squiggles Standard Black Push handle	
122-2710	Mygo/Squiggles Height Adjustable Push Handle (Grey Metal)	

Mygo Base



Indoor



Indoor/Outdoor

	120-818	Interface to Attach Squiggles Seat to Mygo Chassis	
	102-2753-G	Hi-Low Mygo Grey Indoor Manual Chassis - With Height Adjustable Handle	
	102-2754-G	Hi-Low Mygo Grey Indoor Powered Chassis - With Height Adjustable Handle	
	102-2766-G	Hi-Low Mygo Grey Indoor Hydraulic Chassis - With Height Adjustable Handle	
	102-2767-G	Hi-Low Mygo Grey Indoor/Outdoor Manual Chassis - With Height Adjustable Handle	
	102-2768-G	Hi-Low Mygo Grey Indoor/Outdoor Hydraulic Chassis - With Height Adjustable Handle	
	102-2769-G	Hi-Low Mygo Grey Indoor/Outdoor Powered Chassis - With Height Adjustable Handle	
	102-2753	Hi-Low Mygo Black Indoor Manual Chassis - With Push Handle	
	102-2754	Hi-Low Mygo Black Indoor Powered Chassis - With Push Handle	
	102-2766	Hi-Low Mygo Black Indoor Hydraulic Chassis - With Push Handle	
	102-2767	Hi-Low Mygo Black Indoor/Outdoor Manual Chassis - With Push Handle	
	102-2768	Hi-Low Mygo Black Indoor/Outdoor Hydraulic Chassis - With Push Handle	
	102-2769	Hi-Low Mygo Black Indoor/Outdoor Powered Chassis - With Push Handle	

Tray



Tray




Padded Tray Cover



Grab Rail

	120-672	Squiggles Seat tray	
	120-688	Squiggles Seat tray 50 mm extended height	
	117-769	MYGO/Squiggles Grab Rails	
	188-1854-05	Squiggles Black Padded Tray Cover	

Covers



Blue

Green

Pink

	120-710-04	Squiggles Seat covers - pink (Includes femorals seat and backrest)	
	120-710-03	Squiggles Seat covers - blue (Includes femorals seat and backrest)	
	120-710-01	Squiggles Seat covers - green (Includes femorals seat and backrest)	

Shoulder Support



	120-702	Squiggles Seat Shoulder Support Assembly	
	120-828-04	Squiggles Seat Pink Shoulder Cushion	
	120-828-03	Squiggles Seat Blue Shoulder Cushion	
	120-828-01	Squiggles Seat Green Shoulder Cushion	

Leg Guide			
			
	120-856-04	Squiggles Seat Leg Guide Pad Assembly - pink	
	120-856-03	Squiggles Seat Leg Guide Pad Assembly - blue	
	120-856-01	Squiggles Seat Leg Guide Pad Assembly - green	
Footrest and Pommel			
 			
	120-641	Squiggles Seat pommel	
	120-639	Squiggles Seat Footrest	
Back Spacer			
			
	120-779-04	Squiggles Seat 1" back spacer pad & flipaway lateral support - Pink	
	120-779-03	Squiggles Seat 1" back spacer pad & flipaway lateral support - Blue	
	120-779-01	Squiggles Seat 1" back spacer pad & flipaway lateral support - Green	
	120-777-04	Squiggles Seat 1" back spacer pad & rigid lateral support - Pink	
	120-777-03	Squiggles Seat 1" back spacer pad & rigid lateral support - Blue	
	120-777-01	Squiggles Seat 1" back spacer pad & rigid lateral support - Green	
Chest Support			
 			
	137-808-05	Oval Chest Harness - Black	
	137-808-05VIN	Oval Chest Harness - Black Vinyl	
	155-N001-L600	Small Butterfly Harness & attachment straps	
	155-N001-L601	Medium Butterfly Harness & attachment straps	
Foot Support			
   			
		Sandals	
		Booties	
		Ankle Huggers	
		Ankle Wraps	
	152-1600	Sandals Size 1	
	150-1852	Extra Small Booties	
	AH01	Ankle Huggers	
	173-2305	Small Ankle Wraps (use with sandals)	

	120-816	1Inch Footplate Raiser for Single Footplate Head Support	
			
		Contoured	
		Contoured Occipital	
		Flat	
		Laterals (for use with flat)	
	137-699-05	Contoured headrest Assembly and Cushion	
	137-699-05VIN	Contoured headrest Assembly and Cushion - Vinyl	
	117-699-05	Contoured Occipital Headrest	
	117-699-05VIN	Contoured Occipital Headrest - Vinyl	
	117-851-05	Flat Headrest Assembly and Cushion	
	117-851-05VIN	Flat Headrest Assembly and Cushion - Vinyl	
	120-893-05	Flat Headrest Laterals with covers	
	120-893-05VIN	Flat Headrest Laterals with covers Vinyl	
Pelvic Support			
			
		Pelvic Harness	
		Pelvic Spacer	
	120-844-03	Squiggles pelvic harness small - blue	
	120-844-01	Squiggles pelvic harness small - green	
	120-844-04	Squiggles pelvic harness small - pink	
	117-848-03	Squiggles pelvic harness medium - blue	
	117-848-01	Squiggles pelvic harness medium - green	
	117-848-04	Squiggles pelvic harness medium - pink	
	120-789	MYGO/Squiggles Pelvic Harness spacer pads (pair)	
Laterals			
			
		Rigid Laterals	
		Flipaway Laterals	
		Flipaway Laterals with Protraction Pads	
	120-656	Squiggles Seat rigid lateral support - Hardware	
	120-855	Squiggles Seat flipaway lateral support Hardware	
	120-838-03	Squiggles Seat lateral support covers - blue	
	120-838-01	Squiggles Seat lateral support covers - green	
	120-838-04	Squiggles Seat lateral support covers - pink	
	117-644-05	Flipaway Laterals with Protraction Pads	

Product Sizing			
		mm	Inches
Age (Approx)	1-5 Years		
Stature	Min	750	29.5
	Max	1100	43.3
Max User Weight (kg/Lbs)		22Kg	48.4Lbs
Seat Depth (Backrest to back of knee)	Min	190	7.5
	Max	290	11.4
Seat width	Min	160	6.3
	Max	260	10.2
Hip Width (Distance between hip guides)	Min	160	6.3
	Max	220	8.7
Backrest Height (Sitting shoulder height)	Min	320	12.6
	Max	380	15
Backrest Height (with shoulder support)	Min	420	16.5
	Max	480	18.9
Chest Width (Distance between laterals)	Min	160	6.3
	Max	90	3.5
Distance betweenommel & side of seat	Min	190	7.5
	Min	122.5	4.8
Seat to Footplate (Sole of foot to back of knee)	Min	270	10
	Max	650	25
Seat to Floor (Top of Seat to Floor (Hi-Low Chassis))	Prone	10°	10°
	Recline	25°	25°
Backrest Angle	Plantarflexion	10°	10°
	Dorsiflexion	10°	10°
Footplate Angle		460/560	18.1/22
Tray Size			

Comments