

Let's Dream

Indoor rollator

Black TCLD502

Orange TCLD501



Let's Dream, the worldwide success from TrustCare, is a very practical walking aid for people who need a light and flexible rolling aid for indoor use. It also has a beautiful and elegant design. Available in two different colors. It is easy to operate and easy to use in narrow spaces such as in the bathroom. It has a removable tray for transporting, for example, a plate or a cup of coffee. You can use it as a food trolley, a small table to eat from, and as a support if you need to get up from bed in the middle of the night. A fabric basket is available as an accessory and can be used together with the tray.

The handles are clad with soft foam rubber, adjustable in height to fit different statures. The brake can easily be operated with only one hand and is active on both rear wheels. The rollator also has of course a parking brake. It can be easily folded and can stand alone when it is folded.



Black Swan Art. No. LD502



Sunset Orange Art. No. LD501



Transparent Tray



Fabric Basket

Technical specifications:

Tested according to ISO 11199-2:2005

Max. user weight 100 kg

Width 55 cm

Weight 6.5 kg

Handle height is adjustable from 83 – 95 cm

Max. tray load 5 kg

Max. basket load 3 kg

Black tray included

