

Let's Enjoy

Shower chair

Black TCLEN101
White TCLEN201



Let's Enjoy is a new shower chair in a modern design intended for use in the shower or the bathroom. The chair features a beautiful and stable construction for safe and secure use. The legs are easy to adjust to desired height with one easy motion. The chair is intended for intimate cleaning.

Let's Enjoy is adjustable from 42 to 53 cm. The height is adjusted by pressing the buttons on the legs and pulling the lower legs to the desired positions, after which the buttons will be locked in place in the adjustment holes.

Important: Always adjust all the legs to the same height to avoid falling accidents.



White. Art. No. LEN201



Ice Blue. Art. No. LEN301



Copper. Art. No. LEN401

Technical specifications:

Tested according to ISO 12182:2012

Max. user weight 150 kg
Adjustable height 42 – 53 cm
Weight 2.4 kg



Material:

Seat ABS
Tube frame Aluminum
Feet TPE

