

Let's Fly

Outdoor rollator

White	TC0220
Black	TC0230
Red	TC0250
Beige	TC0260
Grey	TC0270



Let's Fly has soft suspension of the rear wheels for smooth run.

Let's Fly is one of the leading rollators in the world today in term of design and functionality. This outdoor rollator is light and steady and its weight is only 6.8 kg. Width only 57 cm. It is simple to fold and can stand alone when folded. Let's Fly is made of cast aluminum and its design is comparable with the most extremely designed products. It has large elegant front wheels that makes it easier to cross curbs. Soft suspension of the rear wheels gives a smooth and soft run. Let's Fly has ergonomically designed rubber handles, adjustable from 78.5 to 92 cm to fit different body lengths. Of course it has an effective brake and parking brake that are easy to operate. Accessories are also available for Let's Fly such as bag, backrest band, crutch holder, foot lever, bag holder and lamps.

When Let's Fly is folded the width is only 20 cm and it is easy to bring with you in a car or a bus.



Crutch holder. Art. No. LF315



Black backrest band. Art. No. 0123
Brown backrest band. Art. No. 0153



Bag. Art. No. LF303

Technical specifications:

Tested according to ISO 11199-2:2005

Max. user weight	130 kg
Width	57 cm
Length	65.5 cm
Weight	6.8 kg
Handle height is adjustable from	78.5 – 92 cm
Seat height	60 cm
Front wheel diameter	24.5 cm
Rear wheel diameter	20 cm
Max. bag load	12 kg

