



**ALLIED  
MEDICAL**  
LIMITED

# USER MANUAL

Pedal Exerciser

MP30401



## Parts of the product



1	Pedal
2	Adjusting knob
3	Steel Frame
4	Tip

## Warranty

1 year for main frame only, not including any consuming parts( wheels, brakes, and so on...)

## Intended Use

Designed for low impact rehabilitation and exercise this Pedal Exerciser is NOT a suitable alternative to a full size gym exercise bike or similar. It is ideal for those with reduced mobility or those who wish to complete a little gentle exercise.

## Caution

Do not stand on the mini exercise bike! Do not overload the pedals! It has risk of damage to the mini exercise bike.

- Do not use the device before carefully reading the user manual.
- Use of the following device for purposes other than those defined in this manual is prohibited.
- For correct use of the device, always consult your doctor or therapist.
- Keep the product packaged away from any heat source; the packaging is made of cardboard.
- The lifespan of the device is determined by the wear of non-repairable and/or replaceable parts.
- Always pay attention to the presence of children.

## Specification Sheet

<b>Product name</b>	PEDAL EXERCISER
<b>Product item number</b>	HE30100010
<b>Product N.W.</b>	2.2KG
<b>Weight capacity</b>	100KG
<b>Product Size L x W x H</b>	578x1850x402mm
<b>Frame material</b>	Steel
<b>Package2</b>	PCS/CTN