

in.basic is a mobility base designed for indoor use. It is compatible with the seats of our strollers: MOUSE and MEWA. It is perfect for home, school and rehabilitation centres.

The **in.basic** indoor base has an adjustable seat mounting height, so it can be adjusted to fit most home tabletops. This gives the child the opportunity to be active while playing, learning or eating under the supervision of a guardian. Each wheel has independent brakes to guarantee stability and safety.

Seats mounted on the **in.basic** frame have all the functionalities of strollers, such as backrest angle adjustment, tilt in space seat adjustment, footrest angle adjustment, lying position and the option to install the seat backwards or forwards.



indications and contraindications.

INDICATIONS

- paralysis and paresis of lower limbs of various origins
- paralysis and paresis of the trunk
- weakening of postural muscles responsible for maintaining the correct body position
- cerebral palsy
- myelomeningocele
- muscular dystrophies

CONTRAINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.

mouse

Mewa



construction.

The **in.basic** mobility base is a light and stable device that changes the image of MyWam strollers. It takes a few moments to put the seat on and there is no need to adjust it.

The seat of the MOUSE and MEWA strollers is attached in the same way as in a regular stroller. Thanks to this solution, one mobility base can be used with three seats and will serve the user for many years.

Small size after folding (107 x 57 x 22 cm) makes it easy to store and transport the **in.basic** mobility base.

**This is a medical device.
Use it according to the instruction manual or label.**

adjustment and features.

Once one of the two MyWam seats – MOUSE or MEWA– is mounted on the **in.basic** mobility base, it is possible to freely modify the parameters of the complete home version stroller.

In order to adjust the parameters of the seat to the individual needs of the user, follow the INSTRUCTIONS for individual strollers.

RECLINING BACKREST



TILT IN SPACE SEAT



SEAT HEIGHT ADJUSTMENT



colour options.

The **in.basic** indoor base is available in WHITE. However, it is possible to fulfill individual color orders. You may order **in.basic** indoor base in your ideal colour design and we will be more than happy to support you in making this dream come true.

☐ **COVERING**
white glossy

parameters.

TECHNICAL INFORMATION

- indoor base length: 95 cm
- indoor base width: 58 cm
- indoor base height: 104 cm
- total weight: 7,5 kg
- seat height regulation: 53 – 71 cm
- length / width / height of the base folded for transport: 107 / 57 / 22 cm

USER PARAMETERS

-  max. patient weight: 40 kg
-  height of the patient: 80-125 cm

contact.



helpis@alliedmedical.co.nz | www.alliedmedical.co.nz | 0800 31 61 81