


Meywalk® MK4



The Meywalk 4 comes standard with:

Dynamic Saddle Seat, Parking & drag brakes, swivel locks, choice of springs and a choice of Trunk support - Universal



Meywalk has a 2 Year Warranty

TGA-ARTG#108757 - Walking aid, two arms, frame, wheeled, Class 1.

DIN EN 12182:1999 - Technical aids for disabled persons.

Description	Item No.	
Select Meywalk Size		
Meywalk 4® Small Green - seat height 50-67cm	013000	
Meywalk 4® Medium Red - seat height 60-84cm	011200	
Meywalk 4® Large Black - seat height 80-104cm	011250	

Select Spring size		
Meywalk 4® Orange Springs - from 10 - 15 Kg - for Micro, Mini	104296	
Meywalk 4® Purple Springs - from 15 - 20 Kg - for Micro, Mini	104297	
Meywalk 4® White Springs - from 20 - 50 Kg - for Micro, Mini, Small	104298	
Meywalk 4® Yellow Springs - from 45-65 Kg - for Mini, Small	104299	
Meywalk 4® Blue Springs - from 65-80 Kg - for Mini, Small	104300	
Meywalk 4® Yellow Springs - from 45 - 65 Kg - for Medium & Large	104292	
Meywalk 4® Blue Springs - 65 - 80 Kg - for Medium & Large	104293	
Meywalk 4® Red Springs - from 80 - 100 Kg - for Medium & Large	104294	
Meywalk 4® Silver Springs - 100 - 130 Kg - for Medium & Large	104295	

Select Trunk Support solution		
Traditional Thoracic Ring - 70cm circumference	011101	
Traditional Thoracic Ring - 90cm circumference	011105	
Traditional Thoracic Ring - 105cm circumference	011109	
Traditional Thoracic Ring - 120cm circumference	011113	
Traditional Thoracic Ring - 130cm circumference	011118	
Universal Thoracic Support Mini / Small - adjusts from 70 - 85cm	100501	 
Universal Thoracic Support Medium - adjusts from 80 - 105cm	100502	
Universal Thoracic Support Large - adjusts from 105 - 135cm	100503	
		Not available for trial
		Not available for trial
		Not available for trial


	VELA Meywalk 4 Micro	VELA Meywalk 4 Small	VELA Meywalk 4 Medium	VELA Meywalk 4 Large
ND WEIGHT				
Length	77 cm	85 cm	86 cm	100 cm
Width, outer measures	51 cm	62 cm	62 cm	68 cm
Width, inner measures	42 cm	49 cm	50 cm	57 cm
Seat width:				
Front	7 cm	9 cm	10 cm	10 cm
Middle	5,2 cm	7,5 cm	7,2 cm	7,2 cm
Back	11,5 cm	13,5 cm	15,5 cm	15,5 cm
Seat depth	24 cm	27 cm	32,5 cm	32,5 cm
Seat height/stride length*	25-45 cm	44-67 cm	60-84 cm	80-104 cm
Height of body ring*	52-65 cm	74-99 cm	101-133 cm	120-152 cm


Meywalk® MK4



Description	Item No.		
-------------	----------	--	--

Accessories for Meywalk MK 4®

Large Basket for MK4, including mounting brackets	11810		Not available for trial	<input type="checkbox"/>
Small Basket for MK4, including mounting brackets	11811		Not available for trial	<input type="checkbox"/>
Universal Brick unit for trunk support (to increase by 4")	100505		Not available for trial	<input type="checkbox"/>
Cover pad for trunk support, silicone	100506		Not available for trial	<input type="checkbox"/>
NB: Please order both 100505 & 100506 to extend Trunk Support by 4"				

Belt, wide for Universal Trunk support	100511		Not available for trial	<input type="checkbox"/>
Trunk Pad, Small Model, Complete (7" x 3.5" opening)	011710		Not available for trial	<input type="checkbox"/>
Trunk Pad, Medium Model, Complete (8" x 4.5" opening)	011715		Not available for trial	<input type="checkbox"/>
Trunk Pad, Large Model, Complete (10" x 5.5" opening)	011720		Not available for trial	<input type="checkbox"/>
Trunk Pad, Adjustment Set	011725		Not available for trial	<input type="checkbox"/>
Leg Divider, Silver, Micro, Mini & Small	013130		Not available for trial	<input type="checkbox"/>
Leg Divider, Cover, Micro, Mini & Small	013150		Not available for trial	<input type="checkbox"/>
Leg Divider, Silver, Medium & Large	011150		Not available for trial	<input type="checkbox"/>
Leg Divider, Cover, Medium & Large	011740		Not available for trial	<input type="checkbox"/>
Leg Guides, Micro and Mini - SET	013148		Not available for trial	<input type="checkbox"/>
Leg Guides, Small - SET	013168	Not available for trial	<input type="checkbox"/>	
Leg Guides, Medium - SET	011160	Not available for trial	<input type="checkbox"/>	
Leg Guides, Large - SET	011161	Not available for trial	<input type="checkbox"/>	
Leg Guides, Medium/Large - black plastic - SET	104303	Not available for trial	<input type="checkbox"/>	
Leg Guides, Small - black plastic - SET	013703	Not available for trial	<input type="checkbox"/>	
Leg Guides, Mini - black plastic - SET	012703	Not available for trial	<input type="checkbox"/>	
Leg Guides, Micro - black plastic - SET	014703	Not available for trial	<input type="checkbox"/>	
Hip Pads, Complete - SET	013075	Not available for trial	<input type="checkbox"/>	
Anti-tip Support, Micro, Mini, Small & Medium - SET	011147			<input type="checkbox"/>
Anti-tip Support, Large - SET	011148			<input type="checkbox"/>
Height Reducing Fittings, Small & Medium - SET	011163		Not available for trial	<input type="checkbox"/>
Height Reducing Fittings, Large - SET	011164		Not available for trial	<input type="checkbox"/>
Extension of brake arm (define length)	DJCUSTOM		Not available for trial	<input type="checkbox"/>
Comfort Saddle for MK4	011334		Not available for trial	<input type="checkbox"/>
Tray, Mini & Small	011735		Not available for trial	<input type="checkbox"/>
Tray, Medium & Large (23"W x 12.5" D)	011730		Not available for trial	<input type="checkbox"/>
Mini Saddle seat (with bolts and washers (mini seat for small) **	106056		Not available for trial	<input type="checkbox"/>
Saddle bracket for mini (order with 106056) **	012370		Not available for trial	<input type="checkbox"/>
Push and drag bar for MK 4 small and mini	011760		Not available for trial	<input type="checkbox"/>
Push and drag bar for MK 4 micro	011763		Not available for trial	<input type="checkbox"/>
Sling Seat (M/L) with brackets, complete, medium and Large - SET	011800		Not available for trial	<input type="checkbox"/>
Sling seat (S) with brackets, complete, Micro, Mini and Small	011805		Not available for trial	<input type="checkbox"/>