



**ALLIED  
MEDICAL**  
LIMITED

# MyWam Mewa Stroller

DEMW-MEWA-GRY





# Mewa

The **MEWA special needs stroller** is designed for the youngest users who are unable to walk or sit in a correct position as a result of injuries, diseases or motor dysfunctions. The stroller is suitable for children at a height of **90 – 125 cm** and a maximum weight of **40 kg**.

Thanks to the backrest tilt and recline feature as well as adjustable torso and head pads you may easily adjust **MEWA** to the required position and ensure the user with optimal protection and stability.

The shock-absorbing suspension on all four, large, gel-filled wheels together with the tilt-foot facilitate travel on uneven terrain. Modern and stylish design is an additional advantage of this stroller.

### INDICATIONS

paralysis and paresis of lower limbs of various origins, paralysis and paresis of the trunk, weakening of postural muscles responsible for maintaining the correct body position, cerebral palsy, myelomeningocele, muscular dystrophies.

### CONTRAINDICATIONS

In case of any contraindications to sitting position &/or passive verticalization stroller adjustments should be consulted with physiotherapist or physician-in-charge.



MEDICAL DEVICE - CLASS I

## Colour.



**UPHOLSTERY**  
synthetic



**FRAME**  
spray painted / glossy



# features. safety and comfort

## TILT IN SPACE SEAT

simple mechanism



5,5° - 21,5°

4 positions

## RECLINING BACKREST

wide range of adjustments



8° - 86°

4 positions

## ADJUSTABLE AND ANGLE TILT FOOTREST

simple and quick adjustment



64° - 180°

5 positions

## SEAT MINIMIZING KIT

reducing inserts



The SEAT MINIMIZING KIT in MEWA is a set of 4 comfortable, independent inserts that allows you to reduce both width and depth of the seat.

## STABILIZATION SYSTEM

adjustable pads



The head and torso support consists of two pairs of independent pads that are height, width and angle adjustable. Its wide range of adjustments and simplicity gives you a fast and effective adaptation to the user's individual needs.

## LYING POSITION

wide range of backrest recline



With a fully adjustable seat and footrest you may quickly alter the stroller into a cozy bed perfect for a short rest or a nap.



## REVERSIBLE SEAT

arrange your travel easily



The design of the stroller allows you to easily rotate the seat opposite to the driving direction. Thanks to this feature you can keep an eye on your child all the time.

## FOLDING LIGHTWEIGHT FRAME

ease of storage,  
transportation and dismantling



The construction of main elements allows you to dismantle MEWA into 3 separate units (frame, seat and wheels).

## DRIVING SUPPORT

comfort for the bumpiest  
journeys



The stroller is equipped with a shock-absorbing suspension on all four wheels and an adjustable rear suspension in the rear wheels, swivel front castor wheels with option to lock, central brake and a tilt-foot.

# functional accessories. in practice



# accessories. wide range of equipment



seat minimizing kit



seat minimizing kit



seat minimizing kit



head / torso adjustable pads



abduction belts



adjustable wedge



safety barrier



therapy tray



tourist bag



shopping basket



rain cover



mosquito net

## TECHNICAL PARAMETERS OF THE STROLLER

seat width: 25\*/32 cm

effective seat depth: 24 – 28 cm

backrest height: 62 cm

backrest width: 20\*/30 cm

distance between footplate and seat: 17 – 32 cm, 5 positions

\*values are given with pads

total weight: 20 – 28 kg

stroller length: 94 – 126 cm

stroller height: 124 – 129 cm

stroller width: 65 cm

**IMPORTANT NOTE!** The choice of a medical device appropriate for a specific disease entity and condition of the patient must be made with the help of a specialist doctor or physiotherapist. Before use, read the operating manual and the label or consult your doctor, as this medical device may not be suitable for you.