

Quantum Outback 4x4 Powerchair Script



Customer Details

Client						
Assessor						
Organisation						
User Weight						
Funder	Whaikaha		ACC		Other	

For assistance with scripting please contact your **Technical Product Specialist** or call **0800 31 61 81** to be put directly through to them.

Standard Specification

Expandable QL3 Electronics for Lights and Indicators	
Smooth Ride Suspension (SRS)	Base Mounted LED Lights
4 High Torque IP67 Motors	Pneumatic All-Terrain (Knobby)

Seat Size Width x Depth

Width		Depth	
12	19	12	19
13	20	13	20
14	21	14	21
15	22	15	22
16		16	
17		17	
18		18	



Rigidiser Bar

Yes	No
-----	----

*Not compatible with power recline

Charger

Yes	No
-----	----

Powered Seating

Non-Tilt	Lift & Tilt
Tilt	Recline

Seat to Floor Height

Lift and Tilt	Tilt Only
17.25"	18.5"
18.25"	19.5"

Colour

Tan	Sage
-----	------

Controller

Controller Type	Standard
	Attendant Only
	Q-Logic 3e
	Other (please specify in notes)

Controller Side

Left hand side
Right hand side
Midline
N/A

Controller Mount

Fixed in-line
Swing away
Retract4
N/A

Legrest Type

Centre Mount
Powered Centre Mount
70 Degree Swing Away
60-90 Degree Swing Away
Manual Elevating Swing Away
Manual Elevating Centre Mount
Powered Swing Away
Paediatric

Stump Support

Left
Right
Both
None

Armrest Type

Two post flip back
Flip-up armrests (not compatible with recline)
Quick height adjust drop-in
Heavy duty
Paediatric
Recline armrests (only compatible with recline)
None

Footplate Type & Size

Standard

Angle Adjustable

6 x 8
8 x 11

Centre Mount

5 x 8
5 x 11

Paediatric

Not required

Armrest Length & Type

Length

Desk
Full

Type

Standard
Gel
Channel arms

Hip Belt Type & Size

Size

Small	Medium	Large
-------	--------	-------

Type

Standard
2 Point
4 Point
Other (please specify)

Seat Pan to Footplate Length

Left

14
15
16
17
18
19
20
Not required
Other (please specify)

Right

14
15
16
17
18
19
20
Not required
Other (please specify)

Heel Loops

Yes

No

Headrest

Not required

Stealth

10" Non-flip down mount
10" Flip down mount
14" Non-flip down mount
14" Flip down mount
Other (please specify in notes)

Seating - Cushion

Not Required
Varilite Stratus
Varilite Evolution
Ride Java
Ride Forward
Ki Mobility Axiom P
Ki Mobility Axiom SP Visco
Ki Mobility Axiom SP Fluid
Ottobock Cloud
Other (please specify)

Seating - Back Support

	Not required
	Varilite Icon Mid
	Varilite Icon Tall
	Varilite Icon Deep
	Ride Designs Java Back Regular
	Ride Designs Java Back Tall
	Ride Designs Java Back Tall Deep
	Ride Designs Java Back Grande
	Ride Designs Java Back Grande Deep
	Axiom Back
	ADI Aluminium Back Active Contour (inc height in notes)
	ADI Aluminium Back Deep Contour (inc height in notes)
	ADI Q-Back Back Active Contour (inc height in notes)
	ADI Q-Back Back Deep Contour (inc height in notes)
	ADI Flat Top Back Active Contour (inc height in notes)
	ADI Flat Top Back Deep Contour (inc height in notes)
	ADI Carbon Fibre Back Active Contour
	Other (please specify)

Seating - Laterals (Please Specify)

	Not required		Thoracic		Pelvic
--	--------------	--	----------	--	--------

Tray

	Swing Away		Custom
	Drop In		None

Accessories

--

Set-up Instructions & Notes

--