

PIPER Comfort

Hybrid semi adjustable standard rehab-buggy



- extended backrest recline
- extended backrest height
- large adjustable footrest
- very compact foldable buggy
- one-piece foldable handlebar
- load capacity up to 50 kg
- suitable for already very small children (with seat size reduction set)

Size

UNI



PIPER Comfort is the new state of the art comfort rehab-buggy This buggy is determined for short stay usage of children with special needs. Fully foldable undercarriage is produced from lightweight aluminum profiles. Buggy is equipped with seat base with adjustable calf rest, extended reclining backrest, solid wheels of 175/225 mm, big detachable and adjustable footrest, 5-point immo belt with shoulder pads. Backrest is equipped with flexible multi-anchoring system for optimal accessory shoulder alignment. Those features make the products very well suitable for intown activities and public transportation, but as well as for light outdoor activities. The variability of the buggy and the ability to attach a wide range of accessories, ensure fulfilling of majority requirements of user's wide spectrum.

The correct usage of the product should be consulted with a seating specialist. Products corresponds to the newest EC standards as well as hygienic, technical and safety requirements.

Standard equipment:

- seat depth non-adjustable
- extended backrest height non-adjustable
- backrest recline adjustable
- footrest height adjustable, detachable
- seat width adjustable by FLEXI-pads set (*) (optional eq.)
- 5-point immo.harness with crotch belt and shoulder padding
- 175/225 mm solid wheels, front with swivel lock
- detachable wheels, forks
- one-piece foldable handle
- basket, transportation lock
- frame colours and cover patterns according to the current offer



folded buggy

Optional equipment:



Flexipads set



Headrest adjustable



Abduction Block



Feet fixation



Canopy Classic



Raincover




Footmuff - Summer



Footmuff - Winter



Seat inlay (seat size reduction)



Footrest raiser

.... and more according to actual list

Growing functions



with the Seat inlay and footrest raiser set is possible to use for small children from two years of age



Positioning



seat positioning



footrest positioning



DESIGN OPTIONS



G06



G07



G12



G13

with Seat inlay / footrest raiser

PIPER Comfort

Backrest height [mm]	680	680
Seat depth [mm]	240	290
Seat width [mm]	230/270	350
Seat to footrest [mm]	220-225	250-375
Backrest angle	90° - 135°	90° - 135°
Seat angle	26°	26°
Total length [mm]	1080	1080
Total width [mm]	590	590
Total height/with handle [mm]	100	1000
Brake	foot oper.	foot oper.
Weight	14,5 kg	14,5 kg
Maximum load	50,0 kg	50,0 kg
Folded buggy length [mm]	1020	1020
Folded buggy width [mm]	310	310
Folded buggy height [mm]	300	300