



THE ADJUSTABLE RAMP SYSTEM



helpis@alliedmedical.co.nz | www.alliedmedical.co.nz | 0800 31 61 81

FOR BARRIER-FREE ACCESS

The adaptable ramp system PONDUS impresses with its lightweight and compact design. The heat-treated, high strength aluminum construction with the innovative perforated surface is extremely durable and allows safe use throughout the year, even in wet conditions.

FEATURES AT A GLANCE

- The ramps are fitted with an ergonomic carrying handle and a transport lock.
On the underside, an anti-slip coating prevents the ramps from slipping.
The PONDUS rail ramps are available in 3 versions: foldable, telescopic and as a combination of both variants so foldable and telescopic.
- The ramps are equipped with plastic plain bearings, which guarantee smooth functioning of the telescopic mechanism.
The ramps are maintenance-free. The anodised surface provides a very high resistance to scratches and gives an attractive appearance.



Details



Transportation lock



Ergonomic carrying handle on both sides



Plastic sliding bearings to ensure smooth operation



With a non-slip and perforated surface

COLOURS

Aluminum & Blue



Available in five versions:


S – Set (2 pcs.) 1-part

F2 – Set (2 pcs.) foldable: 2-part

T2 – Set (2 pcs.) telescopic: 2-part

T3 – Set (2 pcs.) telescopic: 3-part

FT3 – Set (2 pcs.) foldable & telescopic: 3-part



SIZES

PONDUS S: SET (2 PCS.) 1-PART:

LENGTH	ITEM NO.Set	AXLE WIDTH	WEIGHT Set	LOAD CAPACITY _{max}	STEP HEIGHT _{max}
55 cm	PS000550	21 cm	2,9 kg	400 kg	11 cm
120 cm	PS001170	21 cm	5,9 kg	400 kg	23 cm

PONDUS F2: SET (2 PCS.) FOLDABLE: 2-PART

Approved medical aid number: 22.50.01 .0063

LENGTH	ITEM NO.Set	AXLE WIDTH	WEIGHT Set	LOAD CAPACITY _{max}	STEP HEIGHT _{max}
110 cm	PF201100	21 cm	7,1 kg	400 kg	22 cm
150 cm	PF201500	21 cm	9,2 kg	400 kg	30 cm
200 cm	PF202000	21 cm	11,9 kg	300 kg	40 cm

PONDUS T2: SET (2 PCS.) TELESCOPIC: 2-PART

Approved medical aid number: 22.50.01 .0062

LENGTH	ITEM NO.Set	AXLE WIDTH	WEIGHT Set	LOAD CAPACITY _{max}	STEP HEIGHT _{max}
120 cm	PT201190	19,5 cm	7,4 kg	350 kg	24 cm
200 cm	PT202030	19,5 cm	11,4 kg	300 kg	40 cm

PONDUS T3: SET (2 PCS.) TELESCOPIC: 3-PART

Approved medical aid number: 22.50.01 .0062

LENGTH	ITEM NO.Set	AXLE WIDTH	WEIGHT Set	LOAD CAPACITY _{max}	STEP HEIGHT _{max}
200 cm	PT301990	18 cm	12,2 kg	300 kg	40 cm
290 cm	PT302890	18 cm	16,5 kg	250 kg	58 cm

PONDUS FT3: SET (2 PCS.) FOLDABLE & TELESCOPIC: 3-PART

LENGTH	ITEM NO.Set	AXLE WIDTH	WEIGHT Set	LOAD CAPACITY _{max}	STEP HEIGHT _{max}
280 cm	PFT32770	19,5 cm	16,6 kg	250 kg	60 cm