

Powerchair Script



Customer Details

Client			
Assessor			
Organisation			
ACC	<input type="checkbox"/>	DSS	<input type="checkbox"/>
User Weight			

Powerchair Script

For assistance with scripting please contact your **Technical Product Specialist** or call **0800 31 61 81** to be put directly through to them.

Base Type (Select PWC)



Edge HD - Mid wheel, 204kg max user weight.



J623 2.0 - Mid wheel, 136kg max user weight.



Q1450 - Front wheel, 272kg max user weight.



R44 - Rear wheel, 136kg max user weight.

Powered Seating

<input type="checkbox"/>	Non-Tilt
<input type="checkbox"/>	Tilt

<input type="checkbox"/>	Lift & Tilt
<input type="checkbox"/>	Recline

Tyres

<input type="checkbox"/>	Pneumatic
<input type="checkbox"/>	Solid - Standard for Z, HD, 1450
<input type="checkbox"/>	Knobby Tyre

Seat Size Width x Depth

Width

<input type="checkbox"/>	12"
<input type="checkbox"/>	13"
<input type="checkbox"/>	14"
<input type="checkbox"/>	15"
<input type="checkbox"/>	16"
<input type="checkbox"/>	17"
<input type="checkbox"/>	18"
<input type="checkbox"/>	19"
<input type="checkbox"/>	20"
<input type="checkbox"/>	21"
<input type="checkbox"/>	22"
<input type="checkbox"/>	23" <small>not available on J623/R44</small>
<input type="checkbox"/>	24" <small>not available on J623/R44</small>
<input type="checkbox"/>	Other - Please specify

Depth

<input type="checkbox"/>	12"
<input type="checkbox"/>	13"
<input type="checkbox"/>	14"
<input type="checkbox"/>	15"
<input type="checkbox"/>	16"
<input type="checkbox"/>	17"
<input type="checkbox"/>	18"
<input type="checkbox"/>	19"
<input type="checkbox"/>	20"
<input type="checkbox"/>	21"
<input type="checkbox"/>	22"
<input type="checkbox"/>	23" <small>not available on J623/R44</small>
<input type="checkbox"/>	24" <small>not available on J623/R44</small>
<input type="checkbox"/>	Other - Please specify

Rigidiser Bar

<input type="checkbox"/>	Yes	<input type="checkbox"/>	No
--------------------------	-----	--------------------------	----

Charger

<input type="checkbox"/>	Yes	<input type="checkbox"/>	No
--------------------------	-----	--------------------------	----

Controller

Controller Type	
Standard - QL3 EX	VR2 - J623 2.0 only
Standard - QL3e	Other (please specify in notes)
Attendant Only	

Controller Side		Controller Mount	
Left hand side		Fixed in-line	
Right hand side		Swing away	
Midline		Retract4	
N/A		N/A	

Seat to Floor Height

17.5"	19.5"
18.5"	Other

Legrest Type

Centre Mount
Powered Centre Mount
Articulating Footplate Centre Mount - Not available on J623 2.0
Manual Elevating Centre Mount

70 Degree Swing Away	Stump Support
60-90 Degree Swing Away	
Manual Elevating, Swing Away	
Powered Elevating, Swing Away	
Paediatric	
	Left
	Right
	Both
	None

Footplate Type & Size

Standard	
Angle Adjustable	Centre Mount
6" x 8"	5" x 8"
8" x 11"	5" x 11"
Paediatric	Not required

If Centremount selected, select footplate option below.

Footplate Options for Centre Mount

No F/Plate selected
Tapered F/Plate Width Ext
Tapered F/Plate Outside Edge
Tapered F/Plate Width Ext & Outside Edge

Heel Loops

Yes	No
-----	----

Armrest Type

2 Post Flip Back
Flip-up Armrests - Not compatible with recline
Quick Height Adj Drop-In
Heavy Duty
Paediatric
Recline Armrests - Only compatible with recline
Recline Armrests with Down post - Only compatible with recline
None

Armpad Style and Length - Channel

Channel Sml	OB436A2=3-7
Channel Med	OB436A2=2-7
Channel Lge	OB436A2=1-7
Channel One Piece	OB436A3=2-7

Armpad Style and Length - Gel

Gel 2"x10"	GO210G
Gel 2"x12"	GO212G
Gel 2"x14"	GO214G
Gel 3.5"x10"	GO310G
Gel 3.5"x12"	GO312G
Gel 3.5"x14"	GO314G
Gel 3.5"x16"	GO316G
Gel Waterfall 4"x14"	GO414GWT

Armpad Style and Length - Standard/Recline

Standard Desk 10"	PRIDSETLRAM1115
Standard Full 14"	PRIDSETLRAM1116
Recline Desk 10"	PRIDEDGN1801407
Recline Full 14"	PRIDEDGN1801411
Recline Full 17"	PRIDEDGN1801530

Seat Pan to Footplate Length

Left	Right
14"	14"
15"	15"
16"	16"
17"	17"
18"	18"
19"	19"
20"	20"
Not required	Not required
Other (please specify:)	Other (please specify:)

Seating - Cushion

<input type="checkbox"/>	Not Required
<input type="checkbox"/>	Vicair Vector O2
<input type="checkbox"/>	Vicair Adjuster O2
<input type="checkbox"/>	Vicair Active O2
<input type="checkbox"/>	Varilite Stratus
<input type="checkbox"/>	Varilite Evolution
<input type="checkbox"/>	Ride Java
<input type="checkbox"/>	Ride Forward
<input type="checkbox"/>	Ki Mobility Axiom P
<input type="checkbox"/>	Ki Mobility Axiom SP Visco
<input type="checkbox"/>	Ki Mobility Axiom SP Fluid
<input type="checkbox"/>	Ottobock Cloud
<input type="checkbox"/>	Other (please specify:)

Headrest

Stealth	
<input type="checkbox"/>	Not required
<input type="checkbox"/>	10" Non-flip down mount
<input type="checkbox"/>	10" Flip down mount
<input type="checkbox"/>	14" Non-flip down mount
<input type="checkbox"/>	14" Flip down mount
<input type="checkbox"/>	Other (please specify:)

Thigh Guides

<input type="checkbox"/>	Unilink Thigh Guide Tall Height Adj	SPUTT6A1260000
<input type="checkbox"/>	Unilink Pelvic Pad 4x12	SPUTP54D056100

Hip Belt Type & Size

Size

<input type="checkbox"/>	Small	<input type="checkbox"/>	Medium	<input type="checkbox"/>	Large
--------------------------	-------	--------------------------	--------	--------------------------	-------

Type

<input type="checkbox"/>	Standard
<input type="checkbox"/>	2 Point
<input type="checkbox"/>	4 Point
<input type="checkbox"/>	Other (please specify:)

Accessories

Seating - Back Support

<input type="checkbox"/>	Not required
<input type="checkbox"/>	Varilite Icon Mid
<input type="checkbox"/>	Varilite Icon Tall
<input type="checkbox"/>	Varilite Icon Deep
<input type="checkbox"/>	Ride Designs Java Back Regular
<input type="checkbox"/>	Ride Designs Java Back Tall
<input type="checkbox"/>	Ride Designs Java Back Grande
<input type="checkbox"/>	Ride Designs Java Back Grande Deep
<input type="checkbox"/>	Axiom Back
<input type="checkbox"/>	ADI Aluminium Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Aluminium Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Q-Back Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Active Contour (inc height in notes)
<input type="checkbox"/>	ADI Flat Top Back Deep Contour (inc height in notes)
<input type="checkbox"/>	ADI Carbon Fibre Back Active Contour
<input type="checkbox"/>	Other (please specify:)

Seating - Laterals (Please Specify)

<input type="checkbox"/>	Not required	<input type="checkbox"/>	Thoracic	<input type="checkbox"/>	Pelvic
--------------------------	--------------	--------------------------	----------	--------------------------	--------

Tray

<input type="checkbox"/>	Swing Away	<input type="checkbox"/>	Custom
<input type="checkbox"/>	Drop In	<input type="checkbox"/>	None

Set-up Instructions & Notes