



Change of the mounting screws on the rear wheel with the aluminum push rim



Previous Appearance



Updated Appearance



- ✔ Karma has changed the mounting screws from L35mm hex socket binding head screw to L40mm phillips pan head screw for easier assembly. Note that this change will NOT affect the total width of the wheelchair.
- ✔ The updated version is backward compatible with the existing one.
- ✔ The updated version was implemented by the batch of serial no. TM22-056430.
- ✔ The updated version was implemented on Flexx, Flexx HD, Flexx Tall, Flexx Junior, and KM-8000 series.

Order Information

Article number	156000000262	RearTire Assy with Aluminium Rim 20"x1 3/8 GrayPneumatic Tire Q AluminiumRound Tube Handrim
	420400000075	Alu. Handrim Assy with Rivet 20" 6 Holes
	420400000076	Alu. Handrim w/rivet silver 22" 6 Holes
	420400000078	Alu. Handrim w/rivet champagne 22" 6 Holes
	420400000080	Handrim Screw Kit Assy 25L
	420400000171	RearTire Assy with Aluminium Rim 22"x1 3/8 GrayPU Tire Q AluminiumRound Tube Handrim
	420400000173	RearTire Assy with Aluminium Rim 20"x1 3/8 GrayPU Tire Q AluminiumRound Tube Handrim

ECN-00003651