



QUANTUM® OUTBACK® 4x4
shown in *Tan* with 10" (25.4cm)
safe seat elevation and 14"
(38.1cm) pneumatic tires



The Quantum® Outback® 4x4 Powerchair is your ultimate adventure companion, ready to tackle almost any terrain with confidence. Designed for superior performance, this rugged powerchair features powerful 4-pole sealed motors and all-wheel drive for unbeatable traction and stability. Its robust construction ensures a smooth and comfortable ride on uneven surfaces, featuring advanced (SRS) suspension and off-road pneumatic 14" knobby tyres. Whether you're navigating gravel paths, bush tracks, or country properties, the Outback® 4x4 powerchair makes it a breeze. Experience unmatched freedom and independence with a powerchair that combines strength, durability, and cutting-edge technology to conquer any challenge!

Q
QUANTUM®
#1 FOR REHAB POWER

Specifications

DRIVE WHEELS	Standard: 14" knobby pneumatic (38.1cm)
DRIVETRAIN	4 high torque motors
SUSPENSION	Automotive-grade smooth ride suspension
MAXIMUM SPEED ¹	Up to 10 km/h (5.6 km/h when elevated)
MAXIMUM OBSTACLE CLIMBING ABILITY	10.2 cm
MAXIMUM RATED SLOPE	8° (14.1%)
GROUND CLEARANCE	10.7 cm (at motors)
TURNING RADIUS ²	125 cm
OVERALL LENGTH ²	106.7 cm
BASE WIDTH	73.2 cm
SEATING SIZES	
TRU-Balance® 3 Power Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
TRU-Balance® 3 Power Lift & Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Power Tilt ^{5,9}	19-20" (48.3 cm - 50.8 cm)
TRU-Balance® 3 Power Lift & Tilt ^{5,9}	20-21" (50.8 cm - 53.4 cm)
BATTERY SIZE	Group 27
BATTERY WEIGHT ⁷	26.6 kg each
RANGE PER CHARGE ³	Up to 25.75 km at 136 kg Up to 27.4 km at 99.8 kg
AVAILABLE ELECTRONICS	120A Q-Logic 3
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-pole motors
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT	76.5 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded.

QUANTUM®
#1 FOR REHAB POWER




Standard Features

- 4-pole sealed motors
- Maximum speed of up to 10km/h 5.6 km/h when elevated)
- LED lighting package
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions

Options

- 10" (25.4 cm) of power adjustable safe seat elevation at speeds of up to 5.6 km/h
- Quantum® backup camera
- Glove box
- Backpack holder
- Multi-use accessory bag
- Hydration system
- Accessory bar

Colors



0800 31 61 81 | helpis@alliedmedical.co.nz | alliedmedical.co.nz