


Quantum TRU-Balance[®] 3 Power Positioning System Range of Motion



There is a power inhibit on your TRU-Balance[®] 3 Power Positioning System which will not allow you to drive your power chair while in the tilt position. A full inhibit will activate at 20°-50°. See below.

Inhibit Matrix

	Recline Lockout	Tilt Lockout	Elevate Lockout	1/4 Speed	Full Drive Lockout	iLevel [®]
Elevated				X		X
Tilted more than 20°			X		X	
Reclined more than 20°			X		X	
Tilt/reclined more than 125°	X	X	X		X	