



**ALLIED
MEDICAL**
LIMITED

USER MANUAL

SANCTUARY COMMODE

NRSS08237



Sanctuary™ Commode

Product Description

Designed to resemble a premium lounge chair, this commode combines function with an attractive form. This comfortable, high-quality chair turns into a commode by lifting off the vanity cushion and the lid of the commode pan. The wooden lower frame works to conceal the commode pan.

Supplied with the legs and commode pan detached, assembly simply needs the legs and lower frame fastening to the chair with nuts and washers, then the pan fitting into position. The commode includes a simple tool for assembly.

Product Technical Details

- Solid wood lower frame and legs. Metal upper frame surrounded by fire-retardant upholstery.
- Foam padding to chair and removable vanity cushion.
- Removable 3.5 litre capacity commode pan with lid and carry handle.
- Product Weight: 12.2 kg (2 st)
- Maximum User Weight: 160 kg (25 st).
- For use in dry areas only - not showers or wet-rooms.



Replaces

- F11924 Traditional Commode

Spares

- S38219 Potty Kit Spare
- S38220 Spare PU Seat

Standard Features

- NRS guarantee - refer to the Terms and Conditions on our website
- Designed and developed in the UK

Features and Benefits

- Attractive contemporary design to blend with most decor
 - Wooden lower frame and removable vanity cushion to conceal the commode pan
 - Quick assembly - use the tool supplied to fit the legs and frame to prepare for use
 - Easy-clean commode pan with handle and lid for easy transport
 - Washable fire-retardant upholstery, solid wood legs and a steel upper frame for safety, durability and ease of cleaning
-