

# Serial Number Location Guide



Product	Location
<b>J623 2.0</b>	Label inside the LHS of battery casing or on the inside of the rear cover (on the bottom of the frame down below the powerbase module). White label with 14 characters, starting with JB.
<b>J623</b>	Label inside the LHS of Battery casing, white label with 14 characters, starting with JB. Also on the chairs manual.
<b>Edge 2.0</b>	Label inside the LHS of Battery casing, white label with 14 characters, starting with a JD. Also on the chairs manual.
<b>Edge 3</b>	Label inside rear compartment, under neath the power module.
<b>R44</b>	Front of base, inside of frame, by power module.
<b>EdgeZ</b>	Label on outside battery frame, LHS , above the motor, starts with JC
<b>Edge HD</b>	Label in outside battery frame, LHS , above the motor, starts with JD
<b>Q610</b>	Remove rear cover on base, label is just on the frame inside Edge Stretto Inside the rear door, starts with JF
<b>Edge Stretto</b>	Inside the rear door, starts with JF
<b>Q1450</b>	Rear of frame, LHS, starts with a JC
<b>4Front / 4Front 2 / 4Front 2 HD</b>	Under the power module, at the rear of the base. Please remove the rear shroud, to inspect the sub frame for the details
<b>KP25.2</b>	LHS lower frame rail
<b>KP25</b>	Upper RHS frame, under seat rail
<b>Karma 8020</b>	RHS, lower frame
<b>Karma Flexx</b>	RHS, lower frame
<b>Karma Flexx Jnr</b>	Cross brace under the seat upholstery
<b>Ki Mobility Rogue / Rogue XP</b>	Cross brace at rear
<b>Ki Mobility Focus</b>	LHS frame, by rear tyre
<b>Ki Mobility Cat 5</b>	Lower RHS frame, by rear tyre
<b>Ki Mobility Liberty</b>	Lower RHS frame, in front of the rear tyre

