

END-TO-END

COMPARISON CHART

Scan the QR Code to view
our range in your home



ELEVATE YOUR LIFE



Model	Abingdon	Alperton (Petite)	Brixton	Ealing	Eton	Kennington	Ludlow	Marcos	Mercer	Mercer (Petite)	Pembroke (Petite)	Seagrove	Somerset	Studio	Studio (Petite)	Winslow	Winslow (Petite)
Features	Quad motor lift chair with headrest & lumbar function	Petite quad motor lift lay-flat recliner with headrest & lumbar function	Triple motor lift lay-flat recliner with independent head & legrest	Quad motor lift lay-flat recliner with headrest & lumbar function	‘Box-in-box’ triple motor lift chair with headrest function	Duo seat, quad motor lift lay-flat recliner with headrest & lumbar function	Quad motor lift chair with headrest & lumbar function	Quad motor lift lay-flat recliner with headrest, leg extension & lumbar function	Quad motor lift chair with headrest & lumbar function	Petite quad motor lift chair with headrest & lumbar function	Petite triple motor lift lay-flat recliner with independent head & legrest	Quad motor lift lay-flat recliner with headrest & lumbar function	Quad motor lift chair recliner with headrest & lumbar function	Quad motor lift chair with headrest & lumbar function	Petite quad motor lift chair with headrest & lumbar function	Five motor lift chair with headrest & lumbar function	Five motor lift chair with headrest & lumbar function
Weight Capacity	158.7 kg	158.7 kg	158.7 kg	158.7 kg	180 kg	158.7 kg	158.7 kg	158.7 kg	158.7 kg	136 kg	136 kg	158.7 kg	158.7 kg	158.7 kg	136 kg	158.7 kg	136 kg
Seat Width	58.8 cm	43.1 cm	53.3 cm	53.3 cm	53.3 cm	57.1 cm	53.3 cm	57.1 cm	50.8 cm	45.7 cm	45.7 cm	55.8 cm	66 cm	57.1 cm	50.8 cm	53.3 cm	45.7 cm
Seat Depth	53.3 cm	46.9 cm	48.2 cm	48.2 cm	52 cm	54.6 cm	53.3 cm	55.8 cm	54.6 cm	45.7 cm	45.7 cm	53.3 cm	55.8 cm	54.6 cm	45.7 cm	48.2 cm	45.7 cm
Seat-to-Floor	47 cm	45.7 cm	49.5 cm	49.5 cm	48.2 cm	49.5 cm	46.9 cm	49.5 cm	46.9 cm	45.7 cm	45.7 cm	49.5 cm	51 cm	46.9 cm	45.7 cm	49.5 cm	45.7 cm
Chair Width	99 cm	74.9 cm	78.7 cm	78.7 cm	87.6 cm	193 cm	83.8 cm	88.9 cm	86.4 cm	81.2 cm	80 cm	99 cm	104.1 cm	96.5 cm	91.4 cm	78.7 cm	81.2 cm
Chair Depth	106.6 cm	93.9 cm	101.6 cm	101.6 cm	95.3 cm	106.6 cm	104.2 cm	105.4 cm	104 cm	96.5 cm	93.9 cm	107.9 cm	111.7 cm	106.6 cm	96.5 cm	101.6 cm	96.5 cm
Chair Height	106.6 cm	99 cm	106.6 cm	106.6 cm	106.7 cm	106.6 cm	105.5 cm	107.9 cm	106.7 cm	101.6 cm	101.6 cm	106.6 cm	83.8 cm	106.6 cm	101.6 cm	106.6 cm	101.6 cm
Operation	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	USB Hand Wand	App & USB Hand Wand	App & USB Hand Wand
Fabric Type & Colour Options	Hannah Coffee	Merino Fresh Latte	Santa Barbara Latte	Bonanza Walnut Stark Onyx	Stonewash Dune Amber	Stonewash Gun Metal	Stonewash Dove Fudge	Mosaic Fudge Charcoal	Badlands Saddle Eclipse	Badlands Saddle Eclipse	Mosaic Fudge Truffle	Picasso Walnut Midnight	Rainier Ocean	Santa Barbara Fog Latte Bark	Santa Barbara Fog Latte Bark	Cappricio Smoke	Cappricio Smoke