

FOR THE ACTIVE WHEELCHAIR USER

Trochanter support

firmer density to provide lateral stability | redistributes pressure to help create an ischial well

Ischial well - Two compartments
for superior pressure reduction in the sacral and IT-coccyx area & extra positioning possibilities

SmartCell – Foam barrier
limits sliding forward

Ergonomically-shaped foam
for easy transfers

**PRESSURE
INJURY
TREATMENT
PLAN**

**100%
BREATHABLE**

**SOLID
SHAPE**

**MACHINE
WASHABLE**

READY FOR USE

**145 STANDARD
SIZES**

- ▲ **Perfect for users who transfer frequently**
- ▲ **Excellent lateral and frontal stability**
- ▲ **Features a nifty storage pouch for personal items**



The Vicair Active 02 is a pre-contoured cushion which consists of 4 individual compartments filled with Smartcells at the back of the cushion and an anatomically contoured foam inlay at the front. Striking a unique balance between preserving skin at high risk areas at the ischials, coccyx and sacrum by redistribution of pressures, stabilizing the pelvis and facilitating easy transfers. For specific needs the filling of each of the 4 SmartCell filled compartments can be easily adjusted. Allowing you to change the density of each compartment benefitting users who: (1) Need correction or accommodation of an assymetric pelvis. (2) Want (partial) off-loading of the ischial, sacral and coccyx area. (3) Benefit from additional lateral support.

AVAILABLE COVERS

Not applicable. The Vicair Active 02 only has one cover serving as inner and outer cover.

PRODUCT PERFORMANCE

Skin protection / Anti decubitus

Pressure redistribution	▲▲▲▲△
Shear force reduction	▲▲▲▲△
Microclimate control	▲▲▲▲△
Impact & vibration damping	▲▲▲▲△

Positioning

Lateral stability	▲▲▲▲▲
Frontal stability	▲▲▲▲▲
Sliding reduction	▲▲▲▲▲
Allows manual adjustments	▲▲▲▲△
Increasing endurance	▲▲▲▲▲

Maintenance & Safety

Durability	▲▲▲▲▲
Reliability	▲▲▲▲▲
Functionality when misused	▲▲▲△△
Easy to maintain	▲▲▲▲△
Level of hygiene	▲▲▲▲△

Functionality & Application

Comfort	▲▲▲▲▲
Ultra lightweight	▲▲▲▲△
Allows easy transfers	▲▲▲▲▲
Available standard sizes	▲▲▲▲▲
Maximum user weight	▲▲▲▲△

TECHNICAL SPECIFICATIONS

Thickness	6 and 9 cm
Average weight	1250 grams
Max. user weight	200 kg
Machine washable	60°C (without foam)

CUSHION SIZES

■ **Standard** (Sizes >53 cm thickness 9cm only)

■ **On request**

		width												
		14	15	16	17	18	19	20	21	22	23	24		
Inch	CM	35	38	40	43	45	48	50	53	55	58	60		
14	35													
15	38													
16	40													
17	43													
18	45													
19	48													
20	50													
21	53													
22	55													
23	58													
24	60													

070069-OCT-EN-REV01/2025