

Especially developed for bariatrics

- ▶ User weight up to 500kg
- ▶ Shaped like the Vicair Adjuster O2
- ▶ Wider size range than other Vicair products

Follows the body's contours
for equal pressure distribution | support IT-coccyx and sacral area | managing pelvic asymmetry

Pre-ischial bar limits sliding forward | converts shear into normal forces



Thigh support

PRESSURE INJURY TREATMENT PLAN

SOLID SHAPE

READY FOR USE

100% BREATHABLE

MACHINE WASHABLE

The Vicair XXtra O2 consists of 5 individual compartments, designed to contour under loading and contains especially designed – large – SmartCells to be able to support bariatric individuals. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials, coccyx, sacrum and trochanters, positioning the pelvis and supporting the femurs. The wheelchair cushion is able to support a body weight up to 500KG. For specific needs the filling of each of the 5 compartments can be easily adjusted, allowing you to change the density of each compartment.

Product performance

Skin protection / Anti decubitus	
Pressure redistribution	◀◀◀◀◀
Shear force reduction	◀◀◀◀◀
Microclimate control	◀◀◀◀◀
Impact & vibration damping	◀◀◀◀◀

Positioning	
Lateral stability	◀◀◀◀◀◀
Frontal stability	◀◀◀◀◀◀
Sliding reduction	◀◀◀◀◀◀
Allows manual adjustments	◀◀◀◀◀◀
Increasing endurance	◀◀◀◀◀◀

Maintenance & Safety	
Durability	◀◀◀◀◀◀
Reliability	◀◀◀◀◀◀
Functionality when misused	◀◀◀◀◀◀
Easy to maintain	◀◀◀◀◀◀
Level of hygiene	◀◀◀◀◀◀

Functionality & application	
Comfort	◀◀◀◀◀◀
Ultra lightweight	◀◀◀◀◀◀
Allows easy transfer	◀◀◀◀◀◀
Available standard sized	◀◀◀◀◀◀
Maximum user weight	◀◀◀◀◀◀

Technical specifications

Thickness	10 cm
Average weight	900 grams
Max. user weight	500 kg
Machine washable	60°C max.

Cushion sizes

Standard
 On request

Inch	width				
	22	24	26	28	30
18	45				
20	50				
22	55				
24	60				
26	65				
28	70				
30	75				

Available covers

Vicair O2 Top Cover, Comfair Cover and Incotec Cover.