



WHILL Model F



The Perfect Travel Powerchair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.*



When not in use, the Model F can be conveniently stored in a closet or out-of-the-way corner of the room.



*Model F is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

Lightweight

Weighing only 24 kgs, the Model F is one of the lightest power chairs on the market today.**



Top Speed



Range



Obstacle Clearance



780 mm turning radius for tight turns

Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



**Weight is 26.7 kg. if battery is not removed. Battery weighs 2.7 kg.

Note: Driving range of 20 km will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

Color Options




- White ● Black ● Light Green
- Red ● Light Blue



Controller can be positioned on either the left or right side.

Accessories



Specifications

Driving Range	Up to 20 km	Maximum Speed***	6 km/h
Foldable	Yes	Max Incline	10°
Turning Radius	780 mm	Obstacle Clearance	35 mm
Seat Sizes Available	18"	Width Between Armrests	Adjustable: 2 positions (17", 19")
Device Width	555 mm / 605 mm	Device Length	935 mm (465 mm when folded)
Device Weight	26,7 kg	Charging Time	~ 5 hours
App Enabled	Yes, with remote control	Air Travel Approved	Yes
Max Weight Capacity	Up to 115 kg	Swappable Arm Covers	5 Colors
Battery	Lithium-Ion (25.3V 10Ah)	Water Resistance	IPX4

***Acceleration and/or deceleration may change when driving on steep slopes.

- This device is intended to provide mobility to individuals who are comfortable with operating a power chair.
 - WHILL Model F is classified as Wheelchair Class A.

