

Low Deep (10" height / 4" contour)																						
Width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
A	4.25	10.8	4.5	11.4	5	12.7	5.25	13.3	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0	5.5	14.0
B	5	12.7	5	12.7	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9	6.25	15.9
C	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4	10	25.4

Medium Deep (13" height / 4" contour)																						
Width	12"		13"		14"		15"		16"		17"		18"		19"		20"					
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm				
A	6.5	16.5	6.5	16.5	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1	6.75	17.1
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8
C	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0	13	33.0

Standard Deep (16" height / 4" contour)															
Width	14"		15"		16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	
C	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	


Standard Extra Deep (16" height / 6" contour)											
Width	16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	7.25	18.4	
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2	
C	16	40.6	16	40.6	16	40.6	16	40.6	16	40.6	

Tall Deep (18" height / 4" contour)											
Width	16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	
C	18	45.7	18	45.7	18	45.7	18	45.7	18	45.7	

Tall Extra Deep (18" height / 6" contour)											
Width	16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	9.5	24.1	
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2	
C	18	45.7	18	45.7	18	45.7	18	45.7	18	45.7	

Extra Tall Deep (20" height / 4" contour)											
Width	16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	11	27.9	11	27.9	11	27.9	11	27.9	11	27.9	
B	7	17.8	7	17.8	7	17.8	7	17.8	7	17.8	
C	20	50.8	20	50.8	20	50.8	20	50.8	20	50.8	

Extra Tall Extra Deep (20" height / 6" contour)											
Width	16"		17"		18"		19"		20"		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
A	11	27.9	11	27.9	11	27.9	11	27.9	11	27.9	
B	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2	
C	20	50.8	20	50.8	20	50.8	20	50.8	20	50.8	



Please note: Some sizes are available for purchase only