

ADI

ADI[®]

Designed with you in mind

Stealth Products
Solutions For Life[®]



Lightweight Positioning Backs

The lightest back supports on the market

ALUMINIUM SERIES

5 Heights
3 Contours



Fits Chair Widths

10in (25.40cm) through 20in (50.80cm)

Support Heights**

Low (10in - 25.40cm)

Medium (13in - 33.02cm)

Standard (16in - 40.64cm)

Tall (18in - 45.72cm)

Extra Tall (20in - 50.80cm)

Contour Depth

Active (2.5in - 6.35cm)

Deep (4.5in - 11.43cm)

Extra Deep (6in - 15.24cm)

HCPCS Coding (USA only)

Active: E2613

Deep / Extra Deep: E2615

Material

Lightweight aircraft quality aluminum

Weight

1.5lb (680 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

ALUMINIUM FLAT TOP SERIES

5 Heights
3 Contours



Material

Carbon Fiber Material

Weight

.820lb (370 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

CARBON FIBRE SERIES

3 Heights / 2 Contours



Fits Chair Widths

12in (30.48 cm) to 19in (48.26 cm)

Support Heights**

Low (10in - 25.40cm)

Medium (13in - 33.02cm)

Standard (16in - 40.64cm)

Contour Depth

Active (2.5in - 6.35cm)

Deep (4.5in - 11.43cm)

Material

Carbon Fiber Material

Weight

.820lb (370 grams)

approx. weight of a 16in -40.64cm- standard height back

Note:

Custom options available

** Support Heights are measured on a base shell, cushion and cover will add approximately 2.5in (6.35cm) of overall height.

MOUNTING

REMOVABLE



Quick Release

All

1.3 lbs (589gm)

Removable mounting with adjustable positioning



Basic Release

AL/AFT Series Only

1.8 lbs (816gm)

ADJUSTABLE



2-Point PRO™

All

1.25 lbs (566gm)

Stationary mounting with adjustable positioning



4-Point PRO™

ALL

1.32 lbs (596gm)

FIXED



Fixed ELITE™

CF Series Only

.60 lbs (272gm)

Stationary mounting with fixed positioning

No Adjustments



Width Adjustable ELITE

ALL

.63 lbs (285gm)

Width: +/- 1in (2.54cm)

Angle: 20° | Depth: 1.25in (31.75cm) | Width: +/- 1in (2.54cm)

For folding frame chair users that require angle, depth and width adjustment.

Introducing the NEW ADI Q-Back

Product Description

The idea behind the design of the Q-Back is to allow for seamless integration of the ADI back support line to the Quantum TB3 Power Positioning System. This combination will grant the consumer an array of back support height, width, depth and profile options for their multifunction-al Quantum TB3 System, resulting in a highly sophisticated seating package.


Features

- The new Q-Back includes a high quality Lycra cover to complement the anti-shear features of the TB3 system, as well as offering water resistance and breathability.
- The hole pattern on the Q-Back is designed for ease of installation to the TB3 back plate as well as offering height adjustability in .5" increments equaling a total of 3.5" of height adjustment.
- Additional slots have been included on the Q-Back to allow for integration of Stealth laterals and positioning straps.
- Pre-drilled holes for mounting head-rests.
- New XX-Tall (22") height option.

Sizes

- Heights – Standard (16"), Tall (18"), X-Tall (20"), XX-Tall (22")
- Widths – 16", 18", 20"
- Contours – Active (2.5"), Deep (4"), Extra Deep (6")
- Profiles – AL Series (Scapular Cut-Out), AFT Series (Flat Top)
- Selected sizes are Indent Only
- Custom options available

Please note that all backs (except the Carbon Fibre series) come standard with Quick Release Mounting Hardware. Carbon Fibre backs come with the Fixed Elite Mounting Hardware. Alternative hardware can be purchased separately.



Easy flush mounting



Pre-drilled holes for headrest



Flush mounting



3.5" in height adjustment



High quality Lycra cover



Designed with you in mind



0800 31 61 81
www.alliedmedical.co.nz
helpis@alliedmedical.co.nz

**ALLIED
MEDICAL**
L I M I T E D