



**ALLIED
MEDICAL
LIMITED**

Allied Medical AFO Quick Reference Guide

	Supralite AFO	Blue Rocker	Kiddiegait	ToeOFF	NaviGAIT	Dictus Band	Step-Smart
							
Product Indications	Lightweight AFO designed for those patients suffering mild to moderate neurological deficit resulting in drop foot.	Supports foot with impaired ability to actively dorsiflex and when ability to perform push off is reduced. Ideal for drop foot from neurological disorders. Can also be used for partial foot amputations.	An orthosis for toddlers that will allow their little feet to move in a more fluid and natural biomechanical gait pattern. Designed to support footdrop, gait deviations, toe-walking with no midfoot collapse, and low tone crouch gait.	Designed for users who have mild to severe drop foot as a result of stroke, Multiple Sclerosis, Polio, Spinal Cord injury or other neuromuscular conditions.	Ideal for people with mild drop foot who need to gain better dorsiflexion. It is intended to support the foot when the ability to actively dorsiflex is reduced due to the impact of a stroke, MS or Post Polio.	Dynamic support for drop foot or similar for which support of dorsiflexion is desirable.	A drop foot brace and AFO ankle brace solution for indications related to Charcot Marie Tooth disease, Diabetic drop foot, peripheral neuropathy, CVA, Stroke and Peroneal Nerve injury.
Construction	Thermoplastic	Carbon Fibre	Carbon Fibre	Carbon Fibre	Padded Thermoplastic	Padded Vinyl	Polypropylene
Colour	White	Blue	Black	Black	Black	Black, Natural	Black
Sizes	Male or Female	Small, Medium, Large or X-Large	Small, Medium, Large or X-Large	Small, Medium, Large or X-Large	Small / Medium or Large / X-Large	Junior, Medium or Large	Small or Large
Orientation	Left or Right	Left or Right	Left or Right	Left or Right	Universal	Universal	Left or Right
User Weight Limit	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Accessories	N/A	Soft Liner Kit Strap Kit	Soft Liner Kit Strap Kit	Soft Liner Kit Strap Kit	Foot Sock (for use without shoes or with shoes without laces), Extra Strap with Coil	Indoor Dictus Band (for use without shoes) Mono Hook, additional elastic bands	Additional Calf Straps, I Straps and additional Compression Components
Product Code	AMLSUPAFOxx	CAMPSTRxx	AL28901x01x	CAMPSTOxx	AL28410001X	ERDSxx	IPSSxx



allard | NZ

allard | NZ

allard | NZ

allard | NZ

DICTUS
The Dictus Band

