



**ALLIED  
MEDICAL**  
LIMITED

# USER MANUAL

Perching Stool with Arms and Backrest

AML30305



## Parts of the product:



1	Backrest
2	Seat
3	Leg
4	Tip
5	Armrest

## Warranty:

1 year for main frame only, not including any consuming parts( wheels, brakes, and so on...)

## Intended purposes

Perching stools provide a stable and supportive seat for individuals who have difficulty standing for long periods but do not require full sitting support. They are particularly useful for individuals who need support while performing tasks at higher surfaces, such as kitchen counters or workbenches. Overall, perching stools offer a balance between sitting and standing, providing support and stability for individuals who require assistance with tasks performed at elevated surfaces.

## Intended users

1. Elderly people who may have mobility issues or difficulty standing for long periods of time but still want to remain active and engaged in daily activities.
2. Individuals with mobility limitations due to conditions such as arthritis, back pain, or balance disorders, which require support while standing for tasks.
3. Patients undergoing rehabilitation from injuries, surgeries, or medical conditions

that affect their ability to stand or sit comfortably for extended periods.

4. People with disabilities who may benefit from the stability and support provided by perching stools while performing daily tasks at higher surfaces.
5. Perching stools can also be used by caregivers who assist individuals with mobility challenges or disabilities in performing tasks that require standing, such as personal care or household chores.

### Warning

Please tighten all screws securely before sitting.

### Specification Sheet

<b>Product name</b>	KITCHEN STOOL
<b>Item No.</b>	HE12680041
<b>Weight capacity(KG)</b>	136KG
<b>Product N.W.(KG)</b>	13.1KG
<b>Product dimension LxWxH(mm)</b>	(480-520)x498x(810-965)mm
<b>Seat dimension LxWxH(mm)</b>	342x294x25mm
<b>Frame material</b>	Steel
<b>Backrest material</b>	Foam
<b>Seat material</b>	Foam
<b>Frame processing</b>	Liquid coating white
<b>Package</b>	1PC/CTN
<b>Product G.W.(KG)</b>	7.36KG
<b>Carton dimension LxWxW(mm)</b>	570*525*570mm
<b>MOQ</b>	500PCS
<b>1x40ft HQ Loading Q'ty</b>	1312PCS
<b>FOB port</b>	YANTIAN, CHINA