



**ALLIED  
MEDICAL**  
LIMITED

Client: \_\_\_\_\_

Assessor: \_\_\_\_\_

Organisation: \_\_\_\_\_

ACC:



MOH:

Other:



Fully Configured Kit	Part Number
Dynamic QR Arm Support, L, Small Hand Block Flex	AS201L
Dynamic QR Arm Support, R, Small Hand Block Flex	AS201R
Dynamic QR Arm Support, L, Large Hand Block Flex	AS202L
Dynamic QR Arm Support, R, Large Hand Block Flex	AS202R
Attachment Hardware (REQUIRED)	Part Number
Armrest Mounting Clamp, Ø19mm (Ø3/4")	AS119
Armrest Mounting Clamp, Ø22mm & 25mm (Ø7/8" & 1")	AS118
Armrest Mounting Plate	AS120
Assessment Kit	Part Number
Dynamic Arm Support Assessment Kit Includes Dynamic QR Arm Support Base Kit, 2 Large Hand Blocks, 2 Small Hand Blocks, all Armrest Mounting Clamps, one Armrest Mounting Plate and one of each fastener kits to mount to all above listed wheelchair manufacturers.	AS001

### Hand Block Sizing Chart

	 <p><b>FINGER LENGTH</b> Measure length of middle finger from distal portion of the first phalangeal joint to tip of finger.</p>	 <p><b>HAND WIDTH</b> Measure hand widths across widest part of first phalangeal joints.</p>
Small	8.5cm - 10.5cm (3 3/8" - 4 1/8")	<= 9.0cm (3 1/2")
Large	10.5cm - 12.0cm (4 1/8" - 4 3/4")	<= 11.0cm (4 1/4")