



**ALLIED
MEDICAL
LIMITED**

Cervical Collar Solutions



NEW!

	SOFT FOAM CERVICAL COLLAR	PATRIOT EXTRICATION COLLAR	PHILADELPHIA TRACHEOTOMY RIGID COLLAR	PHILADELPHIA ADJUSTABLE COLLAR	ATLAS COLLAR	MIAMI J COLLAR	MIAMI J COLLAR THORACIC EXTENSION	MIAMI J SELECT COLLAR
BRIEF DESCRIPTION	Neck Extension injuries e.g.: whiplash. RA and OA pf cervical spine	Immobilisation during emergency transportation. Routine patient care	The Philadelphia tracheotomy collar provides strong, open, anterior support for greater stability of the neck.	The Philadelphia Adjustable collar is constructed of the same high quality materials as the original Philadelphia Tracheotomy Collar, that has been an industry standard for decades, but combines convenient adjustability with immobilisation and comfort.	Immobilisation after trauma / surgery. Degenerative conditions. Spinal Stenosis. Spondylolisthesis.	C-Spine precaution for trauma patients. Stable non-displaced fractures. Immobilisation for pre and post C-Spine surgery. Degenerative conditions. Spinal Stenosis. Spondylolisthesis. Motor Neurone disease.	C-Spine precaution for trauma patients. Stable non-displaced fractures. Immobilisation for pre and post C-Spine surgery. Degenerative conditions. Spinal Stenosis. Spondylolisthesis. Motor Neurone disease.	The Miami J Select cervical collar offers the widest range of height adjustability in the market with a single, universal design that replicates the comfort, immobilisation and quality you expect from a Miami J product.
CONSTRUCTION	Foam covered with unbleached stockinet	Durable high density polyethylene lined with a soft closed cell foam	Laminated foam plastic	Laminated foam plastic	Lightweight polyethylene with removable lined foam pads	Lightweight polyethylene with Sorbatex padding	Lightweight polyethylene with Sorbatex padding	Lightweight polyethylene with anti-microbial Sorbatex padding
COLOUR	Cream	White / Navy	Beige	Beige	White / Grey	White / Navy Blue	White / Navy Blue	White / Grey / Blue
NECK CIRCUMFERENCE SIZING	Small: 28 - 38 cm Medium: 33 - 45.7 cm Large: 38 - 50.8 cm Large Long: 40 - 66 cm	Paediatric: 20.3 - 45.7 cm Adult: 27.9 - 58.4 cm	Small: 25 - 32.5cm Medium: 32.5 - 40cm Large: 40 - 47.5cm	Universal Sizing	X Short: 20.3 - 27.9 cm Short: 27.9 - 43 cm Regular: 33 - 50.8 cm Tall: 33 - 50.8 cm Extra Tall: 33 - 63.5 cm	Super Short, Stout, Bariatric, XS, Short, Regular, Tall	Universal to cover Stout, Tall, Regular, Short	Universal - with 7 adjustable fittings
FEATURES & BENEFITS	<ul style="list-style-type: none"> Mild control, high comfort Simple Velcro closure Machine washable 	<ul style="list-style-type: none"> One piece extrication style collar Ventilation openings to reduce heat and moisture Multi-height adjustment: adult 4 positions, paediatric 3 positions Nasal cannula holders Open rear ventilation panel for fluid drainage, airflow and cervical spine palpitation 	<ul style="list-style-type: none"> Comes in two pieces for easy fitting Available in three sizes and four heights (measured from based of chin to pocket between clavicle bones: 2.25", 3.25" and 4.25") This product is recommended for single patient use 	<ul style="list-style-type: none"> Anterior height adjustability accommodates neck height from 2.25" to 4.25" Tear-off tabs accommodate neck circumference from 10 to 19 inches Restricts cervical spine flexion, extension and rotation This product is recommended for single patient use, reducing required stock holding 	<ul style="list-style-type: none"> Large trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies Easy to customise 2-piece design with shoulder and mandible contours to increase stability and comfort 	<ul style="list-style-type: none"> MRI, CT and X-ray lucent Bioengineered to minimise pressure points in key known areas: chin, occiput, trapezius and clavicle Ant/Post adjustment buttons allow for customisation of chin and occipital areas, maximising comfort Removable Sternal pad for spine lying and swallowing 	<ul style="list-style-type: none"> Same features as the Miami J, but with added advantage of adjustable height to reduce inventory required and over moulded edges that inhibit the development of pressure points 	<ul style="list-style-type: none"> Tamper resistant height-adjustment mechanism Sternal-thoracic extension with micro-adjustment to enhance coverage and compliance while maintaining immobilization Large tracheal opening – improved function for ICU patients Optional reproducible fit tabs and blue patient touch-points to send patients home with confidence – Blue is You
PRODUCT CODE	CEX	PHP-PP - Paediatric PHP-PA - Adult	OSPHP-Txx	OSPHP-AA	OSPHP-A330	OSMJR-XXX	OSMTR900	OSMJSR-101
BRAND								

Information correct as of June 2019. From time to time, product specifications are changed by manufacturers without prior notice